

# Interventional Urology

Endourology - Men's Health - Women's Health - Surgery

---

Product Catalogue 2025



Altis®

Aris®

Bioflex®

Biosoft®

BoNee®

Cystodrain®

Dormia®

Elefant®

Folysil®

Genesis®

Imajin®

Kolibri®

K-Box®

N.stone®

Poseidon®

Pyelostent®

Restorelle®

ReTrace®

Rossello®

Soprano®

Stenostent®

Titan®

Virtue®

Vortek®

X-Flow®

### Coloplast History

Coloplast has a long history of providing quality, innovative products to support endourologists and their patients. Our story began in 1893 in Sarlat, France where Porges Laboratories specialized in the manufacturing of catheters and medical devices for urology. Today, the Coloplast endourology brand is recognized globally for delivering differentiated and clinically-focused disposable stone management products that help make life easier for people with intimate healthcare needs.

Our business includes Ostomy Care, Continence Care, Wound & Skin Care, Interventional Urology and Voice & Respiratory Care. We work closely with you to help patients who have intimate healthcare needs. We listen to your needs and respond with products that help make their lives easier.

Present throughout Europe, we are always by your side.

**Our customer service and sales representatives are available whenever you need them**  
to help you find the right solutions to your specific requirements.

Feel free to visit our website to discover all our solutions:

[www.coloplast.com](http://www.coloplast.com)



\* Paris Area :  
Le Plessis-Robinson (European Headquarters)  
Le Plessis-Pâté (Worldwide Warehouse)  
Fontenay-sous-Bois (French Sales & Marketing Subsidiary)





## Table of Contents

|                                      |           |
|--------------------------------------|-----------|
| <b>Coloplast Professional</b>        | <b>6</b>  |
| <hr/>                                |           |
| <b>Urology</b>                       | <b>9</b>  |
| 1. Endourology                       | 10        |
| 2. Percutaneous access               | 47        |
| 3. Connectors and spigot             | 56        |
| 4. Lower urinary tract               | 58        |
| <hr/>                                |           |
| <b>Surgery</b>                       | <b>73</b> |
| 1. Laparoscopic surgery              | 75        |
| 2. Surgical drainage                 | 76        |
| <hr/>                                |           |
| <b>Men's Health</b>                  | <b>81</b> |
| 1. Erectile Dysfunction              | 82        |
| 2. Male Urinary Incontinence         | 86        |
| <hr/>                                |           |
| <b>Women's Health</b>                | <b>89</b> |
| 1. Stress Urinary Incontinence (SUI) | 90        |
| 2. Pelvic Organ Prolapse (POP)       | 92        |

# Coloplast Professional

## Sharing Our Expertise, Empowering Your Success\_

Our educational programs for doctors offer a structured learning path opportunity to enhance knowledge, skills, and continuous development. These programs range from continuing medical education training courses, e-learning, and hands-on workshops to masterclasses. By joining our European urologist network, you will not only enhance your technical expertise but also play a major role in transforming patient care.

Explore the valuable benefits and opportunities these programs can offer to enhance your practice:

### Knowledge Enhancement

- **Up-to-date Information:**
  - Access the latest research, EAU Guidelines, advancements in endourology, and treatments for Men's and Women's Health.
- **Specialized Expertise:** Deepen your knowledge in:
  - Endourology, laser treatments for kidney stones, and BPH/soft tissue enucleation.
  - Erectile restoration/dysfunction, stress urinary incontinence, and pelvic organ prolapse.
- **Procedure-Specific Learning:**
  - Ureteroscopy: Removal of kidney stones and treatment of ureteral issues.
  - Percutaneous Nephrolithotomy (PCNL): Treatment for large kidney stones.
  - Endoscopic Tumor Resection: Management of bladder and upper tract urothelial cancer.
  - Minimally Invasive BPH Treatments: Techniques such as ThuLEP (Thulium Laser Enucleation of the Prostate).

### Skill Development

- **Clinical Skills:** Improve diagnostic, surgical, and therapeutic techniques.
- **Soft Skills:** Enhance communication, leadership, and teamwork abilities.
- **Technological Proficiency:** Master tools such as robotic surgery, telemedicine, or electronic health records.

### Improved Patient Outcomes

- **Evidence-Based Practice:** Apply the latest findings to provide high-quality care.
- **Holistic Approach:** Better understand patient needs through advanced diagnostics and treatments.
- **Reduced Errors:** Minimize mistakes with updated safety protocols.
- **Faster Recovery Times:** Leverage minimally invasive methods to reduce hospital stays, postoperative pain, and complications.
- **Precision Treatment:** Use endourology procedures for highly targeted interventions, minimizing risks to surrounding tissues.
- **Patient Satisfaction:** Achieve better outcomes, less scarring, and quicker recovery, enhancing patient trust and satisfaction.

## Networking Opportunities

- **Peer Collaboration:** Interact with other professionals to share knowledge and best practices.
- **Mentorship:** Learn from leading experts and mentors in the field.
- **Global Connections:** Join international conferences, research collaborations, and global medical forums.

## Research and Innovation Opportunities

- **Cutting-Edge Advancements:** Contribute to innovations in laser technology, imaging, and robotic-assisted surgeries.
- **Pioneering Research:** Help develop improved diagnostic methods, surgical techniques, and urological devices.
- **Teaching and Mentorship:** Train the next generation of urologists and endourologists.
- **Publication and Recognition:** Publish studies and share your expertise at international congresses.



# Urology - Kidney & Bladder Health



## Urology

|                                           |           |
|-------------------------------------------|-----------|
| <b>1 • Endourology</b>                    | <b>10</b> |
| 1.1. Thulium Fiber Laser                  | 10        |
| 1.2. Visualization                        | 12        |
| 1.3. Stone extractors                     | 16        |
| 1.4. Access sheaths                       | 20        |
| 1.5. K-Box                                | 21        |
| 1.6. Double loop ureteral stents          | 22        |
| 1.7. Single loop ureteral stents          | 37        |
| 1.8. Guidewires                           | 39        |
| 1.9. Ureteral catheters                   | 43        |
| 1.10. Ureteral dilators                   | 46        |
| <b>2 • Percutaneous access</b>            | <b>47</b> |
| 2.1. Percutaneous nephrostomy set         | 47        |
| 2.2. Percutaneous nephrostomy accessories | 49        |
| 2.3. Direct puncture                      | 52        |
| 2.4. Suprapubic drainage                  | 53        |
| 2.5. Urinary Diversion                    | 55        |
| <b>3 • Connectors and spigot</b>          | <b>56</b> |
| 3.1. Connectors and Spigot                | 56        |
| <b>4 • Lower urinary tract</b>            | <b>58</b> |
| 4.1. Prostatic catheters                  | 58        |
| 4.2. Bladder evacuator                    | 64        |
| 4.3. Bladder injection                    | 65        |
| 4.4. Urethral dilation                    | 66        |
| 4.5. Foley catheters                      | 68        |
| 4.6. Catheters without balloon            | 70        |



# Coloplast TFL Drive

Power made Safer\*

## Laser System

**Coloplast TFL Drive is a 60W laser system intended for multiple indications:** Lithotripsy, BPH and Soft Tissue. It's the only TFL laser to have **safety monitoring thanks** to preset safety ranges and defined power limits. Intuitive Interface are designed by surgeon for surgeons and aiming to adapt power to the size of the fiber for specific clinical condition.

The Coloplast TFL Drive interface allows you to adapt energy, frequency and power values to the size of fiber used.

The **TFL technology** is the latest cutting-edge laser technology, bringing highly effective intra-operative performance. Finer and Faster dusting capabilities combined with less stone retropulsion of TFL allow decrease treatment time.

## Coloplast TFL Drive in Key Points

### All-in-one Solution

- Stone
- BPH
- Soft Tissue

### Intuitive Interface

- Defined aiming at optimizing the effectiveness

### Ergonomic design

- Easy transportability
- Limited footprint

### Efficiency

- Finer and faster dusting capabilities
- Less retropulsion
- High ablation speed

### More Focus

- Quiet System
- Noise level <70dBA

### Profesional education

- Support you to learn this technology with Europe's leading experts

## Technical Data

### Laser information

- Type TFL (Thulium Fiber Laser)
- Pulse energy 0.02 – 6 J
- Repetition rate 1 – 2500Hz
- Max. Power 60W (CW and pulse mode)
- Pulse Duration 50µs – CW

### Electrical Data

- Voltage 100-240 vac
- Frequency 50/60Hz
- Power Consumption 1000 VA

- Wavelength 1940nm ± 20nm
- Noise Level max. <70dBA
- Aiming Beam Green 532nm (adjustable) <5mW, class 3R

### Dimensions

- Size 470 (W) x 810 (D) x 1156 (H) mm  
470 (W) x 810 (D) x 940 (H) mm (display closed)
- Weight 100kg
- Cooling system Air cooling system



\* User interface designed with warning limits for safer use.

## 1 • ENDOUROLOGY


### 1.1. Thulium Fiber Laser

## Ordering Information

| Order Ref. | Description                                                                                                                                                                                                                                                                                                                                                                                |
|------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| TFLD01     | Coloplast TFL Drive Laser including: Coloplast TFL Drive Laser System - Pedal Double footswitch - Optical blast shield - Safety goggles (x2) - Power cord - User manual - Accessory case containing: Fiber stripper for optical fibers 300 to 1000µm, Fiber stripper for optical fibers 100 to 400µm, interlock connector, optical blast shield, and key for use with Coloplast TFL Drive. |

## Fibers


TFL is suitable with use with the smallest fiber 150 µm, allowing to free up space inside the ureterscope working channel to increase irrigation and thus visibility.



| Fiber size | Order Ref.         |             |              | Ordre Unit |
|------------|--------------------|-------------|--------------|------------|
|            | Single use         | Reusable 5x | Reusable 10x |            |
| 150µm      | FSI150             |             |              | 5 per box  |
| 200µm      | FSI200             | FSV200      | FSX200       | 5 per box  |
| 272µm      | FSI272             | FSV272      | FSX272       | 5 per box  |
| 272µm      | FBI272<br>Ball tip |             |              | 5 per box  |
| 365µm      | FSI365             | FSV365      | FSX365       | 5 per box  |
| 550µm      | FSI550             | FSV550      | FSX550       | 5 per box  |
| 800µm      | FSI800             | FSV800      | FSX800       | 5 per box  |
| 1000µm     | FSI000             | FSV000      | FSX000       | 5 per box  |


## Spare Parts

| Order Ref.  | Description                         |
|-------------|-------------------------------------|
| 000PCPGL01  | Laser protection goggles            |
| 000PCPCL01  | Ceramic fiber cleaver with slipcase |
| 000PCPIC01  | Interlock connector                 |
| 000PFS0104  | Fiber Stripper 100 to 400µm         |
| 000PFS0310  | Fiber Stripper 300 to 1000µm        |
| 000PCPDFS1  | Pedal Double footswitch             |
| 000PCPOBS1  | Optical Blast Shield                |
| 000PCPCASE  | Accessories case                    |
| 000PKEY001  | Key                                 |
| 000PFLGCASE | TFL Flight case                     |



# CoralView - Single-Use Digital Flexible Ureteroscope\*

- The CoralView is used by physicians to access, visualize and perform procedures in the urinary tract.
- Each unit is supplied with a Leak resistant and rotatable T-valve.
- 120° field of view, 3-50mm depth of view.
- Accurate color reproduction.
- Autohold technology self-locking lever.
- **Box of 5** • Sterile.




| Order Ref. | Description                       | Distal end outer diameter | Nominal insertion tube outer diameter | Maximum insertion portion width | Deflection mode |
|------------|-----------------------------------|---------------------------|---------------------------------------|---------------------------------|-----------------|
| 000P1U10B  | CoralView Ureteroscope – Standard | 7.2 Ch/Fr                 | 8.25 Ch/Fr                            | 2.9 mm                          | Standard        |
| 000P1U10BR | CoralView Ureteroscope – Reverse  | 7.2 Ch/Fr                 | 8.25 Ch/Fr                            | 2.9 mm                          | Reverse         |
| 000P1U20B  | CoralView Ureteroscope – Standard | 6.7 Ch/Fr                 | 7.5 Ch/Fr                             | 2.7 mm                          | Standard        |
| 000P1U20BR | CoralView Ureteroscope – Reverse  | 6.7 Ch/Fr                 | 7.5 Ch/Fr                             | 2.7 mm                          | Reverse         |

\* For France and Spain only.

## Deflection

**Reverse deflection** Moving the lever downward will make the tip bend down




## 1 • ENDOUROLOGY

### 1.2 Visualization

# BubbleView - Single-Use

## Digital Flexible Cystoscope\*


- The BubbleView is used by physicians to access, visualize and perform procedure in the bladder.
- Depth of view: 3mm - 100mm.
- Field of view: 140°.
- Digital scope, HD 920k pixels.
- **Box of 5 • Sterile.**



| Order Ref. | Description                      | Nominal Insertion tube outer diameter | Maximum insertion portion width | Working channel width | Deflection mode |
|------------|----------------------------------|---------------------------------------|---------------------------------|-----------------------|-----------------|
| 000P1C50C  | BubbleView Cystoscope – Standard | 14.4 Ch/Fr                            | 17.4 Ch/Fr                      | 6.6 Ch/Fr             | Standard        |
| 000P1C50CR | BubbleView Cystoscope – Reverse  | 14.4 Ch/Fr                            | 17.4 Ch/Fr                      | 6.6 Ch/Fr             | Reverse         |

\* For France and Spain only.

**Standard deflection** *Moving the lever downward will make the tip bend up*




# ViewHub S13-T

## Video Processor\*

Compatible with both CoralView ureteroscope and BubbleView cystoscope, adapting to various scenarios in outpatient and operating rooms.

### Features


- 1920 x 1080 pixels.
- Snapshot.
- Image zooming.
- Video recording.
- Image freezing.
- Brightness control.
- White balance.



| Order Ref. | Description                   |
|------------|-------------------------------|
| 000P1S13T  | ViewHub S13 - Video processor |

\* For France and Spain only.

### Accessories



| Order Ref.       | Description                          |
|------------------|--------------------------------------|
| 000P1BRACKET     | Bracket to put screen on IV pole     |
| 000P1TROLLEY     | Trolley without holder               |
| 000P1TROLLEYPLUS | Trolley with holder (removable part) |





# Dormia® DualPerc

## A unique concept\_

Dormia® DualPerc offers dual functionality thanks to its quick and easy switch between flexible and rigid configurations.

When used with both configurations, the device is designed to remove stones through a unique percutaneous access, reducing multiple renal punctures and/or the need of additional devices.



### Rigid configuration

Working length: 35 cm




#### Nitinol basket

- Super elasticity and maximum shape memory
- 4 wires cage for capture and retrieval of stones of different sizes and volumes
- Tipless design for direct contact against tissues with minimal trauma
- Specific design for adapted opening in small cavities

#### Stainless-steel cannula

- Rounded distal tip for smoother contact with the tissue and instruments, designed to minimize trauma
- Radiopaque tube for accurate visualization under x-ray



#### Flexible mandrel

- Diameter of the sheath and mandrel permits full deflection of the flexible cystoscope during stone extraction
- Working length designed for optimal compatibility with flexible cystoscope


### Flexible configuration

Full length: 64 cm  
3 Ch/Fr flexible sheath

# 1 • ENDOUROLOGY

## 1.3. Stone extractors

| Order Ref. | Description                                | Rigid Shaft Length | Rigid Shaft Diameter | Flexible Shaft Length | Flexible Shaft Diameter | Basket Diameter | Basket Length | Ordre Unit |
|------------|--------------------------------------------|--------------------|----------------------|-----------------------|-------------------------|-----------------|---------------|------------|
| EXP448     | Dormia® DualPerc Stone Extractor 4.8 Ch/Fr | 35cm               | 4.8 Ch/Fr            | 64cm                  | 3 Ch/Fr                 | 12.5mm          | 23mm          | 5 per box  |
| EXP470     | Dormia® DualPerc Stone Extractor 7 Ch/Fr   | 35cm               | 7 Ch/Fr              | 64cm                  | 3 Ch/Fr                 | 12.5mm          | 23mm          | 5 per box  |




## **Dormia® No-Tip** Nitinol • 4 twisted wires

- Individually packaged • Sterile • Detachable handle.

| Order Ref. | Ch/Fr | Length | Basket size |
|------------|-------|--------|-------------|
| EXN434     | 3     | 90cm   | 14mm        |

Designed for semi-rigid and rigid ureteroscopy.





### Nitinol • 4 round wires.

- Individually packaged • Sterile • Detachable handle.

| Order Ref. | Ch/Fr | Length | Basket size |
|------------|-------|--------|-------------|
| EXN734     | 2.2   | 120cm  | 11mm        |
| EXN934     | 1.5   | 120cm  | 9mm         |

Designed for flexible, semi-rigid and rigid ureteroscopy.  
Supplied with a waterproof connector (Tuohy-Borst).



## **Dormia® N.stone®** Nitinol • 4 helical wires

- Individually packaged • Sterile • Detachable handle.

| Order Ref. | Ch/Fr | Length | Basket size |
|------------|-------|--------|-------------|
| EXT224     | 4     | 90cm   | 15mm        |
| EXT424     | 3     | 90cm   | 15mm        |
| EXT624     | 2.5   | 90cm   | 12.5mm      |



## 1 • ENDOUROLOGY

### 1.3. Stone extractors

# **Dormia® Front** **Nitinol • Frontal design**

• Individually packaged • Sterile.

| Order Ref. | Ch/Fr | Length | Basket size |
|------------|-------|--------|-------------|
| DOG158     | 1,5   | 120cm  | 8mm         |
| DOG151     | 1,5   | 120cm  | 11mm        |
| DOG208     | 2     | 120cm  | 8mm         |
| DOG201     | 2     | 120cm  | 11mm        |


**2 basket sizes opening, 8mm and 11mm to treat multiple stone sizes and volumes**

Reinforced tip with shrink tubing made of polyester to assure durability.

Large opening to the front as in the grasper



Secure wrapping like a basket



# ReTrace®

## Ureteral access sheaths


- Radiopaque ring.
- Hydro-coated reinforced sheath.
- Compatible with a 0.038" stiff guidewire.
- Individually packaged • Sterile.



**Single wire concept: the working guidewire turns into a safety wire.**

**Flexible tip concept: designed to respect the anatomy and prevent ureteral injuries.**

| Order Ref. | Sheath inner diameter (Fr) | Sheath outer diameter (Fr) | Length (cm) |
|------------|----------------------------|----------------------------|-------------|
| ASXL10     | 10                         | 12                         | 28          |
| ACXL10     | 10                         | 12                         | 35          |
| AXXL10     | 10                         | 12                         | 45          |
| ALXL10     | 10                         | 12                         | 55          |
| ASXL12     | 12                         | 14                         | 28          |
| ACXL12     | 12                         | 14                         | 35          |
| AXXL12     | 12                         | 14                         | 45          |
| ALXL12     | 12                         | 14                         | 55          |



- Connector: ergonomic shape.



## 1 • ENDOUROLOGY

### 1.5. K-Box

# K-Box®

## The training model

Unique 3D model for flexible ureteroscopy simulation and training.


It allows to learn:

- To manipulate the scope.
- To manipulate the scope with single use devices as baskets to catch and relocate stones.
- To perform step by step techniques.
- To use laser in the different settings.



Order Ref.

000PKBOX01





## Silicone double loop ureteral stents

The remarkable biocompatibility and flexibility of silicone makes it the material of choice for long-term implantation: improved patient comfort\* and less encrustation material compared to Polyurethane.

Duration of implantation: up to 12 months.



Kits with  
simple pusher

- **Kit composition:** each stent is supplied in a kit with a pusher (simple or steerable) and with or without a guidewire depending on the model.  
Ø Guide: 0.035".  
Steerable pusher: with or without radiopaque ring.

- **Box of 1** (with a guidewire) • Sterile.
- **Box of 5** (without a guidewire) • Sterile.



Kits with  
steerable pusher

| Ch/Fr | Length | O/O                  | O/O                  |                        |                  |                     |                   | O/C                  |                        | C/C                  |
|-------|--------|----------------------|----------------------|------------------------|------------------|---------------------|-------------------|----------------------|------------------------|----------------------|
|       |        | Fixed core guidewire | Fixed core guidewire | Movable core guidewire | Terumo guidewire | Poseidon® guidewire | Without guidewire | Fixed core guidewire | Movable core guidewire | Fixed core guidewire |
| 6     | 16cm   | AJ4A61               | AJ4261               |                        |                  | AJ5061              |                   | AJ4861               | AJ4761                 | AJ4161               |
|       | 20cm   |                      | AJ4262               | AJ4562                 | AJ42D2           | AJ5062              | AJ4C62            |                      |                        |                      |
|       | 22cm   | AJ4A67               | AJ4267               |                        | AJ42D7           |                     | AJ4C67            |                      |                        |                      |
|       | 24cm   | AJ4A63               | AJ4263               | AJ4563                 | AJ42D3           | AJ5063              | AJ4C63            | AJ4863               | AJ4763                 | AJ4163               |
|       | 26cm   | AJ4A64               | AJ4264               | AJ4564                 | AJ42D4           | AJ5064              | AJ4C64            | AJ4864               | AJ4764                 | AJ4164               |
|       | 28cm   | AJ4A65               | AJ4265               | AJ4565                 | AJ42D5           | AJ5065              | AJ4C65            | AJ4865               | AJ4765                 | AJ4165               |
|       | 30cm   |                      | AJ4266               |                        | AJ42D6           | AJ5066              | AJ4C66            |                      |                        |                      |
| 7     | 16cm   | AJ4A71               | AJ4271               |                        |                  | AJ5071              |                   | AJ4871               |                        | AJ4171               |
|       | 20cm   | AJ4A72               | AJ4272               |                        | AJ42E2           | AJ5072              | AJ4C72            |                      |                        |                      |
|       | 22cm   | AJ4A77               | AJ4277               |                        | AJ42E7           |                     | AJ4C77            |                      |                        |                      |
|       | 24cm   | AJ4A73               | AJ4273               | AJ4573                 | AJ42E3           | AJ5073              | AJ4C73            | AJ4873               | AJ4773                 | AJ4173               |
|       | 26cm   | AJ4A74               | AJ4274               | AJ4574                 | AJ42E4           | AJ5074              | AJ4C74            | AJ4874               | AJ4774                 | AJ4174               |
|       | 28cm   | AJ4A75               | AJ4275               | AJ4575                 | AJ42E5           | AJ5075              | AJ4C75            | AJ4875               | AJ4775                 | AJ4175               |
|       | 30cm   | AJ4A76               | AJ4276               | AJ4576                 | AJ42E6           | AJ5076              | AJ4C76            | AJ4876               | AJ4776                 | AJ4176               |
| 8     | 16cm   |                      | AJ4281               |                        |                  | AJ5081              |                   |                      |                        |                      |
|       | 20cm   |                      | AJ4282               |                        |                  | AJ5082              |                   |                      |                        |                      |
|       | 22cm   | AJ4A87               | AJ4287               |                        | AJ42F7           |                     | AJ4C87            |                      |                        |                      |
|       | 24cm   | AJ4A83               | AJ4283               |                        | AJ42F3           | AJ5083              | AJ4C83            |                      |                        |                      |
|       | 26cm   | AJ4A84               | AJ4284               | AJ4584                 | AJ42F4           | AJ5084              | AJ4C84            | AJ4884               |                        |                      |
|       | 28cm   | AJ4A85               | AJ4285               | AJ4585                 | AJ42F5           | AJ5085              | AJ4C85            | AJ4885               |                        |                      |
|       | 30cm   |                      | AJ4286               | AJ4586                 | AJ42F6           | AJ5086              | AJ4C86            |                      |                        |                      |

\* Compared to Percuflex™ plus (Boston Scientific) stent.

# 1 • ENDOUROLOGY

## 1.6. Double loop ureteral stents

# Imajin® hydro-coated

## Silicone double loop ureteral stents

The remarkable biocompatibility of silicone makes it the material of choice for long-term implantation: improved patient comfort\* and less encrustation material compared to Polyurethane.

Hydro-coated stents to facilitate advancement.

Duration of implantation: up to 12 months.

• **Kit composition:** each stent is supplied in a kit with a steerable pusher and with or without a guidewire depending on the model.

• **Ø Guide:** 0.035".

• **Steerable pusher:** with or without radiopaque ring.

• **Box of 1** (with a guidewire) • Sterile.

• **Box of 5** (without a guidewire) • Sterile.



|       |        | O/O                  |                  |                     |                   |
|-------|--------|----------------------|------------------|---------------------|-------------------|
| Ch/Fr | Length | Fixed core guidewire | Terumo guidewire | Poseidon® guidewire | Without guidewire |
| 6     | 16cm   | BCHA61               | BCHG61           | BCHS61              | BCHJ61            |
|       | 20cm   | BCHA62               | BCHG62           | BCHS62              | BCHJ62            |
|       | 22cm   | BCHA67               | BCHG67           | BCHS67              | BCHJ67            |
|       | 24cm   | BCHA63               | BCHG63           | BCHS63              | BCHJ63            |
|       | 26cm   | BCHA64               | BCHG64           | BCHS64              | BCHJ64            |
|       | 28cm   | BCHA65               | BCHG65           | BCHS65              | BCHJ65            |
|       | 30cm   | BCHA66               | BCHG66           | BCHS66              | BCHJ66            |
| 7     | 16cm   | BCHA71               | BCHG71           | BCHS71              | BCHJ71            |
|       | 20cm   | BCHA72               | BCHG72           | BCHS72              | BCHJ72            |
|       | 22cm   | BCHA77               | BCHG77           | BCHS77              | BCHJ77            |
|       | 24cm   | BCHA73               | BCHG73           | BCHS73              | BCHJ73            |
|       | 26cm   | BCHA74               | BCHG74           | BCHS74              | BCHJ74            |
|       | 28cm   | BCHA75               | BCHG75           | BCHS75              | BCHJ75            |
|       | 30cm   | BCHA76               | BCHG76           | BCHS76              | BCHJ76            |
| 8     | 16cm   | BCHA81               | BCHG81           | BCHS81              | BCHJ81            |
|       | 20cm   | BCHA82               | BCHG82           | BCHS82              | BCHJ82            |
|       | 24cm   | BCHA83               | BCHG83           | BCHS83              | BCHJ83            |
|       | 26cm   | BCHA84               | BCHG84           | BCHS84              | BCHJ84            |
|       | 28cm   | BCHA85               | BCHG85           | BCHS85              | BCHJ85            |
|       | 30cm   | BCHA86               | BCHG86           | BCHS86              | BCHJ86            |

\* Compared to Percuflex™ plus (Boston Scientific) stent.

Duration  
of implantation:  
up to 12 months.

**Vortek®**

## Double loop ureteral stents

The double layer structure makes it easy to insert over the guidewire, the stiffness of each size has been adapted in order to keep the mechanical properties.




- **Kit composition:** each stent is supplied in a kit with a pusher (simple or steerable), with or without a guidewire and with or without a Flush catheter depending on the model.

Ø Guide: 0.035".

**Steerable pusher:** with or without radiopaque ring.

- **Box of 1**, except for Order Ref. ACB6xx (box of 5) • Sterile.



Kits with  
simple pusher

| Ch/Fr                                   | Length | O/O                  |                   |                               | O/C                  |
|-----------------------------------------|--------|----------------------|-------------------|-------------------------------|----------------------|
|                                         |        | Fixed core guidewire | Without guidewire | Terumo guidewire + Flush ACP6 | Fixed core guidewire |
| Stent 4.8 Ch/Fr + Pusher 4.8 Ch/Fr 75cm | 12cm   | ACB550               |                   |                               | ACBC50               |
|                                         | 16cm   | ACB551               |                   |                               | ACBC51               |
|                                         | 20cm   | ACB552               |                   |                               | ACBC52               |
|                                         | 22cm   | ACB557               |                   |                               | ACBC57               |
|                                         | 24cm   | ACB553               |                   |                               | ACBC53               |
|                                         | 26cm   | ACB554               |                   |                               | ACBC54               |
|                                         | 28cm   | ACB555               |                   |                               | ACBC55               |
|                                         | 30cm   | ACB556               |                   |                               | ACBC56               |
| Stent 6 Ch/Fr + Pusher 6 Ch/Fr 40cm     | 24cm   | ACB563               | ACB663            |                               | ACBC63               |
|                                         | 26cm   | ACB564               | ACB664            |                               | ACBC64               |
|                                         | 28cm   | ACB565               | ACB665            |                               | ACBC65               |
|                                         | 30cm   | ACB566               | ACB666            |                               | ACBC66               |
| Stent 7 Ch/Fr + Pusher 7 Ch/Fr 40cm     | 24cm   | ACB573               |                   |                               | ACBC73               |
|                                         | 26cm   | ACB574               | ACB674            | ACBG74                        | ACBC74               |
|                                         | 28cm   | ACB575               | ACB675            | ACBG75                        | ACBC75               |
|                                         | 30cm   | ACB576               | ACB676            |                               | ACBC76               |
| Stent 8 Ch/Fr + Pusher 7 Ch/Fr 40 cm    | 26cm   | ACB584               |                   |                               |                      |
|                                         | 28cm   | ACB585               |                   |                               |                      |
|                                         | 30cm   | ACB586               |                   |                               |                      |

# 1 • ENDOUROLOGY

## 1.6. Double loop ureteral stents

|                                                                |        | <br>Kits with<br>simple pusher |                                            |                              |                                                 | <br>Kits with<br>steerable pusher |                              |                      |                                         |
|----------------------------------------------------------------|--------|-----------------------------------------------------------------------------------------------------------------|--------------------------------------------|------------------------------|-------------------------------------------------|----------------------------------------------------------------------------------------------------------------------|------------------------------|----------------------|-----------------------------------------|
|                                                                |        | O/O                                                                                                             |                                            |                              |                                                 |                                                                                                                      |                              |                      |                                         |
| Ch/Fr                                                          | Length | Fixed<br>core<br>guidewire                                                                                      | Fixed core<br>guidewire<br>+ Flush<br>ACP2 | Movable<br>core<br>guidewire | Movable<br>core<br>guidewire<br>+ Flush<br>ACP6 | Terumo<br>guidewire                                                                                                  | Poseidon<br>URO<br>guidewire | Without<br>guidewire | Without<br>guidewire<br>+ Flush<br>ACP2 |
| Stent<br>4.8 Ch/Fr<br><br>+<br><br>Pusher<br>4.8 Ch/Fr<br>75cm | 12cm   |                                                                                                                 |                                            |                              |                                                 | ACB1B0                                                                                                               | ACBS40                       |                      |                                         |
|                                                                | 16cm   |                                                                                                                 |                                            |                              |                                                 | ACB1B1                                                                                                               | ACBS41                       |                      |                                         |
|                                                                | 20cm   |                                                                                                                 |                                            |                              |                                                 | ACB1B2                                                                                                               | ACBS42                       |                      |                                         |
|                                                                | 22cm   |                                                                                                                 |                                            |                              |                                                 | ACB1B7                                                                                                               | ACBS47                       |                      |                                         |
|                                                                | 24cm   |                                                                                                                 | ACB253                                     |                              |                                                 | ACB1B3                                                                                                               | ACBS43                       |                      |                                         |
|                                                                | 26cm   |                                                                                                                 | ACB254                                     |                              |                                                 | ACB1B4                                                                                                               | ACBS44                       |                      | ACB454                                  |
|                                                                | 28cm   |                                                                                                                 | ACB255                                     |                              |                                                 | ACB1B5                                                                                                               | ACBS45                       |                      | ACB455                                  |
|                                                                | 30cm   |                                                                                                                 | ACB256                                     |                              |                                                 | ACB1B6                                                                                                               | ACBS46                       |                      |                                         |
| Stent<br>4.8 Ch/Fr<br>+<br>Pusher<br>6 Ch/Fr<br>40cm           | 12cm   | ACB1C0                                                                                                          |                                            |                              |                                                 | ACB1G0                                                                                                               | ACBS50                       |                      |                                         |
|                                                                | 16cm   | ACB1C1                                                                                                          |                                            |                              |                                                 | ACB1G1                                                                                                               | ACBS51                       |                      |                                         |
|                                                                | 20cm   | ACB1C2                                                                                                          |                                            |                              |                                                 | ACB1G2                                                                                                               | ACBS52                       |                      |                                         |
|                                                                | 22cm   | ACB1C7                                                                                                          |                                            |                              |                                                 | ACB1G7                                                                                                               | ACBS57                       |                      |                                         |
|                                                                | 24cm   | ACB1C3                                                                                                          | ACB2C3                                     |                              |                                                 | ACB1G3                                                                                                               | ACBS53                       |                      | ACB4C3                                  |
|                                                                | 26cm   | ACB1C4                                                                                                          | ACB2C4                                     |                              |                                                 | ACB1G4                                                                                                               | ACBS54                       |                      | ACB4C4                                  |
|                                                                | 28cm   | ACB1C5                                                                                                          | ACB2C5                                     |                              |                                                 | ACB1G5                                                                                                               | ACBS55                       |                      | ACB4C5                                  |
|                                                                | 30cm   | ACB1C6                                                                                                          |                                            |                              |                                                 | ACB1G6                                                                                                               | ACBS56                       |                      | ACB4C6                                  |
| Stent<br>6 Ch/Fr<br>+<br>Pusher<br>6 Ch/Fr<br>40cm             | 20cm   | ACB162                                                                                                          |                                            |                              |                                                 | ACB1D2                                                                                                               | ACBS62                       |                      |                                         |
|                                                                | 22cm   | ACB167                                                                                                          |                                            |                              |                                                 | ACB1D7                                                                                                               | ACBS67                       | ACBM67               |                                         |
|                                                                | 24cm   | ACB163                                                                                                          | ACB263                                     | ACB363                       |                                                 | ACB1D3                                                                                                               | ACBS63                       | ACBM63               | ACB463                                  |
|                                                                | 26cm   | ACB164                                                                                                          | ACB264                                     | ACB364                       |                                                 | ACB1D4                                                                                                               | ACBS64                       | ACBM64               | ACB464                                  |
|                                                                | 28cm   | ACB165                                                                                                          | ACB265                                     | ACB365                       |                                                 | ACB1D5                                                                                                               | ACBS65                       | ACBM65               | ACB465                                  |
|                                                                | 30cm   | ACB166                                                                                                          | ACB266                                     | ACB366                       |                                                 | ACB1D6                                                                                                               | ACBS66                       | ACBM66               |                                         |
| Stent<br>7 Ch/Fr<br>+<br>Pusher<br>7 Ch/Fr<br>40cm             | 20cm   | ACB172                                                                                                          |                                            |                              |                                                 | ACB1E2                                                                                                               | ACBS72                       |                      |                                         |
|                                                                | 22cm   | ACB177                                                                                                          |                                            |                              |                                                 | ACB1E7                                                                                                               | ACBS77                       |                      |                                         |
|                                                                | 24cm   | ACB173                                                                                                          | ACB273                                     | ACB373                       |                                                 | ACB1E3                                                                                                               | ACBS73                       | ACBM73               | ACB473                                  |
|                                                                | 26cm   | ACB174                                                                                                          | ACB274                                     | ACB374                       | ACB3C4                                          | ACB1E4                                                                                                               | ACBS74                       | ACBM74               | ACB474                                  |
|                                                                | 28cm   | ACB175                                                                                                          | ACB275                                     | ACB375                       | ACB3C5                                          | ACB1E5                                                                                                               | ACBS75                       | ACBM75               | ACB475                                  |
|                                                                | 30cm   | ACB176                                                                                                          | ACB276                                     | ACB376                       |                                                 | ACB1E6                                                                                                               | ACBS76                       | ACBM76               | ACB476                                  |
| Stent<br>8 Ch/Fr<br>+<br>Pusher<br>7 Ch/Fr<br>40cm             | 24cm   | ACB183                                                                                                          | ACB283                                     |                              |                                                 | ACB1F3                                                                                                               | ACBS83                       |                      |                                         |
|                                                                | 26cm   | ACB184                                                                                                          | ACB284                                     |                              |                                                 | ACB1F4                                                                                                               | ACBS84                       |                      |                                         |
|                                                                | 28cm   | ACB185                                                                                                          | ACB285                                     |                              |                                                 | ACB1F5                                                                                                               | ACBS85                       |                      |                                         |
|                                                                | 30cm   | ACB186                                                                                                          |                                            |                              |                                                 | ACB1F6                                                                                                               | ACBS86                       |                      |                                         |

Duration  
of implantation:  
up to 12 months.

**Vortek®**

## Double loop ureteral stents

Soft outer layer for the patient comfort.

The double layer structure makes it easy to insert over the guidewire, the stiffness of each size has been adapted in order to keep the mechanical properties.



|                                                        |      | O/C                  |
|--------------------------------------------------------|------|----------------------|
|                                                        |      | Fixed core guidewire |
| Stent<br>4.8 Ch/Fr<br>+<br>Pusher<br>4.8 Ch/Fr<br>75cm | 12cm | ACB950               |
|                                                        | 16cm | ACB951               |
|                                                        | 20cm | ACB952               |
|                                                        | 22cm | ACB957               |
|                                                        | 24cm | ACB953               |
|                                                        | 26cm | ACB954               |
|                                                        | 28cm | ACB955               |
|                                                        | 30cm | ACB956               |
| Stent 4.8 Ch/Fr<br>+<br>Pusher 6 Ch/Fr<br>40cm         | 24cm | ACB9C3               |
|                                                        | 26cm | ACB9C4               |
|                                                        | 28cm | ACB9C5               |
| Stent 6 Ch/Fr<br>+<br>Pusher 6 Ch/Fr<br>40cm           | 24cm | ACB963               |
|                                                        | 26cm | ACB964               |
|                                                        | 28cm | ACB965               |
|                                                        | 30cm | ACB966               |
| Stent 7 Ch/Fr<br>+<br>Pusher 7 Ch/Fr<br>40cm           | 24cm | ACB973               |
|                                                        | 26cm | ACB974               |
|                                                        | 28cm | ACB975               |
|                                                        | 30cm | ACB976               |
| Stent 8 Ch/Fr<br>+ Pusher 7 Ch/Fr<br>40cm              | 26cm | ACB984               |
|                                                        | 28cm | ACB985               |

Kits with  
simple pusher

Kits with  
steerable pusher



# 1 • ENDOUROLOGY

## 1.6. Double loop ureteral stents


# Vortek® hydro-coated

## Double loop ureteral stents

Soft outer layer for the patient comfort. The double layer structure makes it easy to insert over the guidewire.

The hydrophilic coating facilitates the advancement of the stent.

Duration of implantation: up to 12 months.




- **Kit composition:** each stent is supplied in a kit with a pusher (simple or steerable), with or without a guidewire and with or without a Flush catheter depending on the model.

Ø Guide: 0.035".

**Steerable pusher:** with or without radiopaque ring.

- **Box of 1** (with a guidewire) • Sterile.
- **Box of 5** (without a guidewire) • Sterile.



|                                         |        | Kits with simple pusher |                  |                               | Kits with steerable pusher |                                   |                  |                               |                        |                   |
|-----------------------------------------|--------|-------------------------|------------------|-------------------------------|----------------------------|-----------------------------------|------------------|-------------------------------|------------------------|-------------------|
|                                         |        | O/O                     |                  |                               | O/O                        |                                   |                  |                               |                        |                   |
| Ch/Fr                                   | Length | Fixed core guidewire    | Terumo guidewire | Terumo guidewire + Flush ACP6 | Fixed core guidewire       | Fixed core guidewire + Flush ACP2 | Terumo guidewire | Terumo guidewire + Flush ACP2 | Poseidon URO guidewire | Without guidewire |
| Stent 4.8 Ch/Fr + Pusher 4.8 Ch/Fr 75cm | 12cm   |                         |                  |                               | BCFA50                     |                                   | BCFG50           |                               |                        | BCFJ50            |
|                                         | 16cm   |                         |                  |                               | BCFA51                     |                                   | BCFG51           |                               |                        | BCFJ51            |
|                                         | 20cm   |                         |                  |                               | BCFA52                     |                                   | BCFG52           |                               |                        | BCFJ52            |
|                                         | 22cm   |                         |                  |                               | BCFA57                     |                                   | BCFG57           |                               | BCFS57                 | BCFJ57            |
|                                         | 24cm   |                         |                  |                               | BCFA53                     |                                   | BCFG53           |                               | BCFS53                 | BCFJ53            |
|                                         | 26cm   |                         |                  |                               | BCFA54                     |                                   | BCFG54           |                               | BCFS54                 | BCFJ54            |
|                                         | 28cm   |                         |                  |                               | BCFA55                     |                                   | BCFG55           |                               | BCFS55                 | BCFJ55            |
|                                         | 30cm   |                         |                  |                               | BCFA55                     |                                   | BCFG56           |                               |                        | BCFJ56            |
| Stent 4.8 Ch/Fr + Pusher 6 Ch/Fr 40cm   | 12cm   |                         |                  |                               | BCFA40                     |                                   | BCFG40           |                               |                        | BCFJ40            |
|                                         | 16cm   |                         |                  |                               | BCFA41                     |                                   | BCFG41           |                               |                        | BCFJ41            |
|                                         | 20cm   |                         |                  |                               | BCFA42                     |                                   | BCFG42           |                               |                        | BCFJ42            |
|                                         | 22cm   |                         |                  |                               | BCFA47                     |                                   | BCFG47           |                               | BCFS47                 | BCFJ47            |
|                                         | 24cm   |                         |                  |                               | BCFA43                     |                                   | BCFG43           |                               | BCFS43                 | BCFJ43            |
|                                         | 26cm   |                         |                  |                               | BCFA44                     |                                   | BCFG44           |                               | BCFS44                 | BCFJ44            |
|                                         | 28cm   |                         |                  |                               | BCFA45                     |                                   | BCFG45           |                               | BCFS45                 | BCFJ45            |
|                                         | 30cm   |                         |                  |                               | BCFA46                     |                                   | BCFG46           |                               |                        | BCFJ46            |
| Stent 6 Ch/Fr + Pusher 6 Ch/Fr 40cm     | 20cm   |                         |                  |                               | BCFA62                     |                                   | BCFG62           |                               |                        | BCFJ62            |
|                                         | 22cm   |                         |                  |                               | BCFA67                     | BCFM67                            | BCFG67           |                               | BCFS67                 | BCFJ67            |
|                                         | 24cm   |                         |                  |                               | BCFA63                     | BCFM63                            | BCFG63           |                               | BCFS63                 | BCFJ63            |
|                                         | 26cm   | BNFA64                  | BNFG64           | BNFQ64                        | BCFA64                     | BCFM64                            | BCFG64           | BCFP64                        | BCFS64                 | BCFJ64            |
|                                         | 28cm   | BNFA65                  | BNFG65           | BNFQ65                        | BCFA65                     | BCFM65                            | BCFG65           | BCFP65                        | BCFS65                 | BCFJ65            |
|                                         | 30cm   |                         | BNFG66           | BNFQ66                        | BCFA66                     |                                   | BCFG66           |                               | BCFS66                 | BCFJ66            |
| Stent 7 Ch/Fr + Pusher 7 Ch/Fr 40cm     | 20cm   |                         |                  |                               | BCFA72                     |                                   | BCFG72           |                               |                        | BCFJ72            |
|                                         | 22cm   |                         |                  |                               | BCFA77                     | BCFM77                            | BCFG77           |                               | BCFS77                 | BCFJ77            |
|                                         | 24cm   |                         |                  |                               | BCFA73                     | BCFM73                            | BCFG73           |                               | BCFS73                 | BCFJ73            |
|                                         | 26cm   | BNFA74                  | BNFG74           | BNFQ74                        | BCFA74                     | BCFM74                            | BCFG74           | BCFP74                        | BCFS74                 | BCFJ74            |
|                                         | 28cm   | BNFA75                  | BNFG75           | BNFQ75                        | BCFA75                     | BCFM75                            | BCFG75           | BCFP75                        | BCFS75                 | BCFJ75            |
|                                         | 30cm   |                         | BNFG76           | BNFQ76                        | BCFA76                     |                                   | BCFG76           |                               | BCFS76                 | BCFJ76            |

Duration  
of implantation:  
up to 12 months.

# Biosoft® duo

## Double loop ureteral stents

Produced by coextrusion of a flexible material and a rigid material, Biosoft duo was designed for easy insertion, while retaining great flexibility for the patient comfort.



- **Kit composition:** each stent is supplied in a kit with a pusher (simple or steerable), with or without a guidewire and with or without a Flush catheter depending on the model.  
Ø Guide: 0.035".  
**Steerable pusher:** with or without radiopaque ring.
- **With a guidewire:** box of 1 • **Without a guidewire:** Box of 5.
- Sterile.

Kits with  
simple pusher

Kits with  
steerable pusher

|                                         |             | O/O                  | O/O                  |                                     |                        |                  |                                 |                                 |                         |                   | O/C                  |                        |
|-----------------------------------------|-------------|----------------------|----------------------|-------------------------------------|------------------------|------------------|---------------------------------|---------------------------------|-------------------------|-------------------|----------------------|------------------------|
| Ch/Fr                                   | Length (cm) | Fixed core guidewire | Fixed core guidewire | Fixed core guidewire + Flush ACP2xx | Movable core guidewire | Terumo guidewire | Terumo guidewire + Flush ACP2xx | Terumo guidewire + Flush ACP6xx | Poseidon URO guidewire* | Without guidewire | Fixed core guidewire | Movable core guidewire |
| Stent 4.8 Ch/Fr + Pusher 4.8 Ch/Fr 75cm | 12cm        | BNAA50               |                      | BCAM50                              | BCAD50                 | BCAGB0           | BCAG50                          | BCAN50                          | BCAS50                  | BCAJ50            | BCAB50               | BCAE50                 |
|                                         | 16cm        | BNAA51               |                      | BCAM51                              | BCAD51                 | BCAGB1           | BCAG51                          | BCAN51                          | BCAS51                  | BCAJ51            | BCAB51               | BCAE51                 |
|                                         | 20cm        | BNAA52               |                      | BCAM52                              | BCAD52                 | BCAGB2           | BCAG52                          | BCAN52                          | BCAS52                  | BCAJ52            | BCAB52               | BCAE52                 |
|                                         | 22cm        | BNAA57               |                      | BCAM57                              | BCAD57                 | BCAGB7           | BCAG57                          | BCAN57                          | BCAS57                  | BCAJ57            | BCAB57               | BCAE57                 |
|                                         | 24cm        | BNAA53               |                      | BCAM53                              | BCAD53                 | BCAGB3           | BCAG53                          | BCAN53                          | BCAS53                  | BCAJ53            | BCAB53               | BCAE53                 |
|                                         | 26cm        | BNAA54               |                      | BCAM54                              | BCAD54                 | BCAGB4           | BCAG54                          | BCAN54                          | BCAS54                  | BCAJ54            | BCAB54               | BCAE54                 |
|                                         | 28cm        | BNAA55               |                      | BCAM55                              | BCAD55                 | BCAGB5           | BCAG55                          | BCAN55                          | BCAS55                  | BCAJ55            | BCAB55               | BCAE55                 |
|                                         | 30cm        | BNAA56               |                      | BCAM56                              | BCAD56                 | BCAGB6           | BCAG56                          | BCAN56                          | BCAS56                  | BCAJ56            | BCAB56               | BCAE56                 |
| Stent 4.8 Ch/Fr + Pusher 6 Ch/Fr 40cm   | 12cm        |                      | BCAA40               | BCAM40                              | BCAD40                 | BCAGC0           | BCAG40                          | BCAN40                          | BCAS40                  | BCAJ40            | BCAB40               | BCAE40                 |
|                                         | 16cm        |                      | BCAA41               | BCAM41                              | BCAD41                 | BCAGC1           | BCAG41                          | BCAN41                          | BCAS41                  | BCAJ41            | BCAB41               | BCAE41                 |
|                                         | 20cm        |                      | BCAA42               | BCAM42                              | BCAD42                 | BCAGC2           | BCAG42                          | BCAN42                          | BCAS42                  | BCAJ42            | BCAB42               | BCAE42                 |
|                                         | 22cm        |                      | BCAA47               | BCAM47                              | BCAD47                 | BCAGC7           | BCAG47                          | BCAN47                          | BCAS47                  | BCAJ47            | BCAB47               | BCAE47                 |
|                                         | 24cm        |                      | BCAA43               | BCAM43                              | BCAD43                 | BCAGC3           | BCAG43                          | BCAN43                          | BCAS43                  | BCAJ43            | BCAB43               | BCAE43                 |
|                                         | 26cm        |                      | BCAA44               | BCAM44                              | BCAD44                 | BCAGC4           | BCAG44                          | BCAN44                          | BCAS44                  | BCAJ44            | BCAB44               | BCAE44                 |
|                                         | 28cm        |                      | BCAA45               | BCAM45                              | BCAD45                 | BCAGC5           | BCAG45                          | BCAN45                          | BCAS45                  | BCAJ45            | BCAB45               | BCAE45                 |
|                                         | 30cm        |                      | BCAA46               | BCAM46                              | BCAD46                 | BCAGC6           | BCAG46                          | BCAN46                          | BCAS46                  | BCAJ46            | BCAB46               | BCAE46                 |
| Stent 6 Ch/Fr + Pusher 6 Ch/Fr 40cm     | 20cm        |                      | BCAA62               |                                     |                        | BCAGD2           | BCAG62                          |                                 | BCAS62                  | BCAJ62            | BCAB62               |                        |
|                                         | 22cm        |                      | BCAA67               |                                     |                        | BCAGD7           | BCAG67                          |                                 | BCAS67                  | BCAJ67            |                      |                        |
|                                         | 24cm        | BNAA63               | BCAA63               | BCAM63                              |                        | BCAGD3           | BCAG63                          |                                 | BCAS63                  | BCAJ63            | BCAB63               | BCAE63                 |
|                                         | 26cm        | BNAA64               | BCAA64               | BCAM64                              | BCAD64                 | BCAGD4           | BCAG64                          | BCAN64                          | BCAS64                  | BCAJ64            | BCAB64               | BCAE64                 |
|                                         | 28cm        | BNAA65               | BCAA65               | BCAM65                              | BCAD65                 | BCAGD5           | BCAG65                          | BCAN65                          | BCAS65                  | BCAJ65            | BCAB65               | BCAE65                 |
|                                         | 30cm        |                      | BCAA66               |                                     | BCAD66                 | BCAGD6           | BCAG66                          |                                 | BCAS66                  | BCAJ66            | BCAB66               |                        |
| Stent 7 Ch/Fr + Pusher 7 Ch/Fr 40cm     | 20cm        |                      | BCAA72               |                                     |                        | BCAGE2           | BCAG72                          |                                 | BCAS72                  | BCAJ72            |                      |                        |
|                                         | 22cm        |                      | BCAA77               |                                     |                        | BCAGE7           | BCAG77                          |                                 | BCAS77                  | BCAJ77            | BCAB77               |                        |
|                                         | 24cm        | BNAA73               | BCAA73               | BCAM73                              |                        | BCAGE3           | BCAG73                          |                                 | BCAS73                  | BCAJ73            | BCAB73               |                        |
|                                         | 26cm        | BNAA74               | BCAA74               | BCAM74                              | BCAD74                 | BCAGE4           | BCAG74                          | BCAN74                          | BCAS74                  | BCAJ74            | BCAB74               | BCAE74                 |
|                                         | 28cm        | BNAA75               | BCAA75               | BCAM75                              | BCAD75                 | BCAGE5           | BCAG75                          | BCAN75                          | BCAS75                  | BCAJ75            | BCAB75               | BCAE75                 |
|                                         | 30cm        |                      | BCAA76               |                                     | BCAD76                 | BCAGE6           | BCAG76                          |                                 | BCAS76                  | BCAJ76            | BCAB76               |                        |
| Stent 8 Ch/Fr + Pusher 8 Ch/Fr 40 cm    | 26cm        | BNAA84               | BCAA84               |                                     |                        | BCAGF4           | BCAG84                          |                                 | BCAS84                  | BCAJ84            | BCAB84               |                        |
|                                         | 28cm        | BNAA85               | BCAA85               |                                     | BCAD85                 | BCAGF5           | BCAG85                          |                                 | BCAS85                  | BCAJ85            | BCAB85               |                        |
|                                         | 30cm        |                      | BCAA86               |                                     | BCAD86                 | BCAGF6           | BCAG86                          |                                 | BCAS86                  | BCAJ86            | BCAB86               |                        |
| Stent 9 Ch/Fr + Pusher 8 Ch/Fr 40 cm    | 26cm        |                      | BCAA94               |                                     |                        |                  |                                 |                                 |                         |                   |                      |                        |
|                                         | 28cm        |                      | BCAA95               |                                     | BCAD95                 |                  |                                 |                                 |                         |                   |                      |                        |
|                                         | 30cm        |                      | BCAA96               |                                     | BCAD96                 |                  |                                 |                                 |                         |                   |                      |                        |

\* Only Italy Poseidon.

# 1 • ENDOUROLOGY

## 1.6. Double loop ureteral stents



Duration  
of implantation:  
up to 12 months.

# Biosoft® duo

## Multi-length hydro-coated double loop ureteral stents

Produced by coextrusion of a flexible material and a rigid material, Biosoft duo was designed for easy insertion, while retaining great flexibility for the patient comfort.

- **Kit composition:** each stent is supplied in a kit with a pusher (simple or steerable) and a guidewire.  
 Ø Guide: 0.035".  
 Steerable pusher: with radiopaque ring.
- **Box of 1 • Sterile.**




Kits with  
simple pusher


| O/O         |              |               |
|-------------|--------------|---------------|
| Ch/Fr Stent | Ch/Fr Pusher | Length Pusher |
| 4.8         | 4.8          | 75cm          |
| 6           | 6            | 40cm          |
| 7           | 7            | 40cm          |

Kits with  
steerable pusher

| O/O         |              |               |
|-------------|--------------|---------------|
| Ch/Fr Stent | Ch/Fr Pusher | Length Pusher |
| 4.8         | 6            | 40cm          |
| 6           | 6            | 40cm          |
| 7           | 7            | 40cm          |

The global stent is between 22 to 28 cm when the multi-length is straightened.





## 1 • ENDOUROLOGY

### 1.6. Double loop ureteral stents

Duration  
of implantation:  
up to 30 days.

# PU-R Double loop ureteral stents

## Double loop ureteral stents in rigid polyurethane

- **Kit composition:** each stent is supplied in a kit with a pusher (simple or steerable) and with or without a guidewire depending on the model.

Ø Guide: 0.035".

**Steerable pusher:** with or without radiopaque ring.

- **Box of 1** (with a guidewire) • Sterile.
- **Box of 5** (without a guidewire) • Sterile.

|       |        | Kits with<br>simple pusher |                         | Kits with<br>steerable pusher |                      |                         |
|-------|--------|----------------------------|-------------------------|-------------------------------|----------------------|-------------------------|
|       |        | O/O                        | O/O                     |                               |                      | O/C                     |
| Ch/Fr | Length | Fixed core<br>guidewire    | Fixed core<br>guidewire | Movable core<br>guidewire     | Without<br>guidewire | Fixed core<br>guidewire |
| 6     | 24cm   | AC4B63                     | AC4663                  | AC4D63                        |                      | AC4563                  |
|       | 26cm   | AC4B64                     | AC4664                  | AC4D64                        | AC4E64               | AC4564                  |
|       | 28cm   | AC4B65                     | AC4665                  | AC4D65                        | AC4E65               | AC4565                  |
|       | 30cm   | AC4B66                     | AC4666                  | AC4D66                        |                      | AC4566                  |
| 7     | 24cm   | AC4B73                     | AC4673                  | AC4D73                        | AC4E73               |                         |
|       | 26cm   | AC4B74                     | AC4674                  | AC4D74                        | AC4E74               | AC4574                  |
|       | 28cm   | AC4B75                     | AC4675                  | AC4D75                        | AC4E75               | AC4575                  |
|       | 30cm   | AC4B76                     | AC4676                  |                               | AC4E76               | AC4576                  |
| 8     | 28cm   | AC4B85                     |                         |                               |                      |                         |
|       | 30cm   | AC4B86                     |                         |                               |                      |                         |

# PU-S Double loop ureteral stents

## Double loop ureteral stents in soft polyurethane

Duration  
of implantation:  
up to 90 days.

- **Kit composition:** each stent is supplied in a kit with a pusher (simple or steerable) and a guidewire.

Ø Guide: 0.035".

**Steerable pusher:** with or without radiopaque ring.

- **Box of 5** • Sterile.

|       |        | Kits with<br>simple pusher |                         | Kits with<br>steerable pusher |
|-------|--------|----------------------------|-------------------------|-------------------------------|
|       |        | O/O                        |                         | O/C                           |
| Ch/Fr | Length | Fixed core<br>guidewire    | Fixed core<br>guidewire | Fixed core<br>guidewire       |
| 6     | 24cm   |                            | AC4263                  |                               |
|       | 26cm   | ACBH64                     | AC4264                  | AC4364                        |
|       | 28cm   | ACBH65                     | AC4265                  |                               |
| 7     | 24cm   |                            | AC4273                  |                               |
|       | 26cm   | ACBH74                     | AC4274                  | AC4374                        |
|       | 28cm   | ACBH75                     | AC4275                  |                               |



Duration  
of implantation:  
up to 30 days.

# Polyamide & Polyurethane double loop ureteral stent sets



- **Set composition:** each stent is supplied in the same peel pouch with a simple pusher and with or without a guidewire depending on the model.

Ø Guide: 0.035".

- **Box of 5** • Sterile.

Sets with  
simple pusher

| Ch/Fr | Length | O/O                  |                   | O/C                  |                   |
|-------|--------|----------------------|-------------------|----------------------|-------------------|
|       |        | Fixed core guidewire | Without guidewire | Fixed core guidewire | Without guidewire |
| 4.8   | 24cm   | BNBA53               | BNBJ53            | BNBB53               | BNBK53            |
|       | 26cm   | BNBA54               | BNBJ54            | BNBB54               | BNBK54            |
|       | 28cm   | BNBA55               | BNBJ55            | BNBB55               | BNBK55            |
|       | 30cm   | BNBA56               | BNBJ56            |                      |                   |
| 6     | 24cm   | BNBA63               | BNBJ63            | BNBB63               | BNBK63            |
|       | 26cm   | BNBA64               | BNBJ64            | BNBB64               | BNBK64            |
|       | 28cm   | BNBA65               | BNBJ65            | BNBB65               | BNBK65            |
|       | 30cm   | BNBA66               | BNBJ66            |                      |                   |
| 7     | 24cm   | BNBA73               | BNBJ73            | BNBB73               | BNBK73            |
|       | 26cm   | BNBA74               | BNBJ74            | BNBB74               | BNBK74            |
|       | 28cm   | BNBA75               | BNBJ75            | BNBB75               | BNBK75            |
|       | 30cm   | BNBA76               | BNBJ76            |                      |                   |
| 8     | 26cm   | BNBA84               | BNBJ84            | BNBB84               | BNBK84            |
|       | 28cm   | BNBA85               | BNBJ85            | BNBB85               | BNBK85            |
|       | 30cm   | BNBA86               | BNBJ86            | BNBB86               | BNBK86            |



1 • ENDOUROLOGY

1.6. Double loop ureteral stents

Duration  
of implantation:  
up to 12 months.


# Double loop ureteral stents

## For ureteral stenosis

**STENOSTENT** • Silicone for long-term implantation.  
12 Ch/Fr reinforced body for maximum resistance to stenosis, 8 Ch/Fr loops for patient comfort.

- **Kit composition:** each stent is supplied in a kit with a steerable pusher (with radiopaque ring) and a guidewire.
- **Ø Guide:** 0.035" (fixed core guidewire).
- **Box of 1** • Sterile.



 Kits with  
steerable pusher

|            | Q/O    |        |        |        |        |
|------------|--------|--------|--------|--------|--------|
| Length     | 16cm   | 24cm   | 26cm   | 28cm   | 30cm   |
| Order Ref. | AJ4W81 | AJ4W83 | AJ4W84 | AJ4W85 | AJ4W86 |

Duration  
of implantation:  
up to 12 months.

# Double loop ureteral stents

## For pyeloplasty\_

Pyelostent silicone for long-term implantation.

12 Ch/Fr reinforced part in renal pelvis for a better healing, 8 Ch/Fr loops for patient comfort.

- **Kit composition:** each stent is supplied in a kit with a steerable pusher (with radiopaque ring) and a guidewire.  
Ø Guide: 0.035" (fixed core guidewire).
- **Box of 1** • Sterile.



Kits with  
steerable pusher

|            | O/O    |        |        |
|------------|--------|--------|--------|
| Length     | 26cm   | 28cm   | 30cm   |
| Order Ref. | AJ4Y84 | AJ4Y85 | AJ4Y86 |

# 1 • ENDOUROLOGY

## 1.6. Double loop ureteral stents

Duration  
of implantation:  
up to 12 months.

# Double loop ureteral stents

## Without eyes on the straight section

**Vortek:** its two-layer structure makes it easy to insert.

**Silicone:** the biocompatibility of silicone makes it the material of choice for long-term implantation.

- **Kit composition:** each stent is supplied in a kit with a pusher (simple or steerable) and a guidewire.

Ø Guide: 0.035" (fixed core guidewire).

**Steerable pusher:** with or without radiopaque ring.

- **Box of 1** • Sterile.



Sets with  
simple pusher

|       |        | O/O Vortek          |                  |
|-------|--------|---------------------|------------------|
| Ch/Fr | Length | Steerable<br>pusher | Simple<br>pusher |
| 6     | 16cm   |                     |                  |
|       | 20cm   |                     |                  |
|       | 22cm   |                     |                  |
|       | 24cm   | ACB763              |                  |
|       | 26cm   | ACB764              | ACB864           |
|       | 28cm   | ACB765              | ACB865           |
|       | 30cm   | ACB766              | ACB866           |
| 7     | 16cm   |                     |                  |
|       | 20cm   |                     |                  |
|       | 22cm   |                     |                  |
|       | 24cm   | ACB773              |                  |
|       | 26cm   | ACB774              | ACB874           |
|       | 28cm   | ACB775              | ACB875           |
|       | 30cm   | ACB776              | ACB876           |
| 8     | 16cm   |                     |                  |
|       | 20cm   |                     |                  |
|       | 22cm   |                     |                  |
|       | 24cm   | ACB783              |                  |
|       | 26cm   | ACB784              |                  |
|       | 28cm   | ACB785              |                  |
|       | 30cm   | ACB786              |                  |

|  |  | O/O Silicone        |                  |
|--|--|---------------------|------------------|
|  |  | Steerable<br>pusher | Simple<br>pusher |
|  |  | AJ4U61              | AJ4B61           |
|  |  | AJ4U62              | AJ4B62           |
|  |  | AJ4U67              | AJ4B67           |
|  |  | AJ4U63              | AJ4B63           |
|  |  | AJ4U64              | AJ4B64           |
|  |  | AJ4U65              | AJ4B65           |
|  |  | AJ4U66              | AJ4B66           |
|  |  | AJ4U71              | AJ4B71           |
|  |  | AJ4U72              | AJ4B72           |
|  |  | AJ4U77              | AJ4B77           |
|  |  | AJ4U73              | AJ4B73           |
|  |  | AJ4U74              | AJ4B74           |
|  |  | AJ4U75              | AJ4B75           |
|  |  | AJ4U76              | AJ4B76           |
|  |  | AJ4U81              | AJ4B81           |
|  |  | AJ4U82              | AJ4B82           |
|  |  | AJ4U87              | AJ4B87           |
|  |  | AJ4U83              | AJ4B83           |
|  |  | AJ4U84              | AJ4B84           |
|  |  | AJ4U85              | AJ4B85           |
|  |  | AJ4U86              | AJ4B86           |


Duration  
of implantation:  
up to 12 months.

## Vortek® Tumor stents

### Without eyes on the straight section

The inner layer is reinforced compared to the standard Vortek ureteral stent. It allows passing through stricture and the stiffness of material offers an excellent resistance to compression.

- **Kit composition:** each stent is supplied in a kit with a steerable pusher (with radiopaque ring) and a guidewire.  
Ø Guide: 0.035" (hydrophilic nitinol guidewire).
- **Box of 1** • Sterile.



Kits with  
steerable pusher

| Ch/Fr | Length | O/O              |                     |
|-------|--------|------------------|---------------------|
|       |        | Terumo guidewire | Poseidon® guidewire |
| 7     | 26cm   | BCCG74           | BCCS74              |
|       | 28cm   | BCCG75           | BCCS75              |
|       | 30cm   | BCCG76           | BCCS76              |
| 8     | 26cm   | BCCG84           | BCCS84              |
|       | 28cm   | BCCG85           | BCCS85              |
|       | 30cm   | BCCG86           | BCCS86              |


1 • ENDOUROLOGY

1.7. Single loop ureteral stents

# Vortek® Single loop ureteral stents

## Device for drainage and urine collection

- **Set composition:**
  - Single loop ureteral stent 90cm.
  - Seldinger fixed core 0.035" guidewire.
  - Clamp.
  - A connector with a Luer end for connection to a Luer-bag.
- **Box of 1 • Sterile.**



|            | Q/C    |        |        |
|------------|--------|--------|--------|
| Ch/Fr      | 6      | 7      | 8      |
| Order Ref. | AC4406 | AC4407 | AC4408 |

Eyes on loop and body.

|            | Q/O    |        |        |
|------------|--------|--------|--------|
| Ch/Fr      | 6      | 7      | 8      |
| Order Ref. | ACA206 | ACA207 | ACA208 |

Eyes on loop and body.

|            | Q/O    |        |        |
|------------|--------|--------|--------|
| Ch/Fr      | 6      | 7      | 8      |
| Order Ref. | ACA106 | ACA107 | ACA108 |

Eyes on loop only.



1 • ENDOUROLOGY  
1.8. Guidewires

# Soprano®

## Exclusive Hybrid Design

Soprano® is a guidewire with an exclusive hybrid design that combines the access of a nitinol wire with the handling of a PTFE-coated wire

### Soprano®

The hydrophilic flexible tip reduces surface friction, and its round distal tip creates a lubricious surface for ease of insertion and navigation.



### Soprano® Dual Flex

The Soprano Dual Flex features a flexible proximal end designed to ease passage and to minimize the risk of scope damage.



| Order Ref. | Description        | Diameter        | Length | Tip shape | Tip style                                                                                 | Quantity |
|------------|--------------------|-----------------|--------|-----------|-------------------------------------------------------------------------------------------|----------|
| AEHA35     | Soprano®           | 0.035" (0.89mm) | 150cm  | Straight  | 5cm hydrophilic straight distal tip with the first 3cm flexible                           | Box of 5 |
| AEHB35     | Soprano® Dual Flex | 0.035" (0.89mm) | 150 cm | Straight  | 5cm hydrophilic straight distal tip with the first 3cm flexible 7.5cm flexible (proximal) | Box of 5 |



# Hydrophilic Nitinol guidewires

## Radifocus™ Terumo guidewires\*

Its configuration and the combination of Nitinol and hydrophilic properties allow the Terumo guidewires to reach difficult access sites.

- Flexible tip length: 3cm.
- Box of 5 • Sterile.


| Diameter | Length | Rigidity | Tip      | Order Ref. |
|----------|--------|----------|----------|------------|
| 0.035"   | 150cm  | Standard | Straight | AE0N35     |
|          |        |          | Angled   | AE0T35     |
|          |        | Stiff    | Straight | AE0P35     |
|          |        |          | Angled   | AE0X35     |



## Poseidon guidewires\*\*

A hydrophilic coating is applied to the guide wire to lower the friction and improve the navigation and trackability.

- Flexible tip length: 3cm.
- Box of 5 • Sterile.



| Diameter | Length | Rigidity | Tip      | Order Ref.  |
|----------|--------|----------|----------|-------------|
| 0.035"   | 150cm  | Standard | Straight | 000P1HYDS35 |
|          |        |          | Angled   | 000P1HYDA35 |
|          |        | Stiff    | Straight | 000P1HYSS35 |
|          |        |          | Angled   | 000P1HYSA35 |



\* Not authorized in Italy.  
\*\* For Italy only.

# PTFE-coated guidewires

## Seldinger guidewires

- Box of 5
- Sterile.

### Fixed core guidewires

Straight tip



|          | 150cm    | 150cm  | 90cm     |
|----------|----------|--------|----------|
| Diameter | Standard | Stiff  | Standard |
| 0.032"   | AE0032   |        | AE0E32   |
| 0.035"   | AE0035   | AE0C35 |          |
| 0.038"   | AE0038   |        | AE0E38   |

5cm flexible tip except AE0C35 11cm.

### Fixed core guidewires

"J" tip



|          | 80cm   |
|----------|--------|
| Diameter | Stiff  |
| 0.038"   | AE1138 |

5cm flexible tip.

### Movable core guidewires

Straight tip



|          | 150cm    |        |
|----------|----------|--------|
| Diameter | Standard | Stiff  |
| 0.032"   | AE0A32   |        |
| 0.035"   | AE0A35   | AE0D35 |
| 0.038"   | AE0A38   |        |

5cm flexible tip.

# Stainless steel guidewires

## Schüller and Lunderquist guidewires

• Box of 5 • Sterile.



### Schüller

| 100cm    |          |        |
|----------|----------|--------|
| Diameter | Straight | "J"    |
| 0.035"   | AE0435   | AEJ435 |

7cm flexible tip.



### Lunderquist

| 80cm     |          |        |
|----------|----------|--------|
| Diameter | Straight | "J"    |
| 0.035"   | AE0535   | AEJ535 |

6cm flexible tip for AE0535.  
7cm flexible tip for AEJ535.

## 1 • ENDOUROLOGY

### 1.9. Ureteral catheters

# Drainage catheters

## Wide range of tips

- Each catheter is supplied with a stylet (metal or polyester) and a colour code connector.
- Material: Neoplex.
- Length: 70cm for AC55xx, AC52xx, ACN6xx and 69cm for AC54xx, AC51xx, AC53xx and AC50xx.
- Catheters are individually packed.
- **Box of 10** • Sterile.

### Straight open whistle tip

| Eyes: 2 - Stylet: metal |        |        |        |        |        |        |        |
|-------------------------|--------|--------|--------|--------|--------|--------|--------|
| Ch/Fr                   | 3      | 4      | 5      | 6      | 7      | 8      | 9      |
| Order Ref.              | AC5503 | AC5504 | AC5505 | AC5506 | AC5507 | AC5508 | AC5509 |

| Eyes: 6 - Stylet: metal |        |        |        |        |        |        |        |
|-------------------------|--------|--------|--------|--------|--------|--------|--------|
| Ch/Fr                   | 3      | 4      | 5      | 6      | 7      | 8      | 9      |
| Order Ref.              | ACN603 | ACN604 | ACN605 | ACN606 | ACN607 | ACN608 | ACN609 |

| Eyes: 2 - Stylet: polyester |        |        |        |        |        |        |        |
|-----------------------------|--------|--------|--------|--------|--------|--------|--------|
| Ch/Fr                       | 3      | 4      | 5      | 6      | 7      | 8      | 9      |
| Order Ref.                  | AC5203 | AC5204 | AC5205 | AC5206 | AC5207 | AC5208 | AC5209 |



### Straight closed olive tip

| Eyes: 2 - Stylet: metal |        |        |        |        |        |        |
|-------------------------|--------|--------|--------|--------|--------|--------|
| Ch/Fr                   | 3      | 4      | 5      | 6      | 7      | 8      |
| Order Ref.              | AC5403 | AC5404 | AC5405 | AC5406 | AC5407 | AC5408 |

| Eyes: 2 - Stylet: polyester |        |        |        |        |        |        |
|-----------------------------|--------|--------|--------|--------|--------|--------|
| Ch/Fr                       | 3      | 4      | 5      | 6      | 7      | 8      |
| Order Ref.                  | AC5103 | AC5104 | AC5105 | AC5106 | AC5107 | AC5108 |



### Straight closed cylindrical tip

| Eyes: 2 - Stylet: metal |        |        |        |        |        |        |
|-------------------------|--------|--------|--------|--------|--------|--------|
| Ch/Fr                   | 3      | 4      | 5      | 6      | 7      | 8      |
| Order Ref.              | AC5303 | AC5304 | AC5305 | AC5306 | AC5307 | AC5308 |

| Eyes: 2 - Stylet: polyester |        |        |        |        |        |        |
|-----------------------------|--------|--------|--------|--------|--------|--------|
| Ch/Fr                       | 3      | 4      | 5      | 6      | 7      | 8      |
| Order Ref.                  | AC5003 | AC5004 | AC5005 | AC5006 | AC5007 | AC5008 |



## Interventional catheters

These catheters are suitable for use with a guidewire during endoscopic procedures

Made of PEBA.

- Each catheter is supplied with a metal stylet and a Luer syringe connector.
- Length: 74cm for ACP20x, APC30x and 73cm for ACP40x, ACP60x.
- Catheters are individually packed.
- **Box of 5** • Sterile.

### Flush straight open tip ureteric catheters without side eye

| Ch/Fr                 | 4.8    | 6      | 7      | 8      |
|-----------------------|--------|--------|--------|--------|
| Guidewire recommended | 0.038" | 0.038" | 0.038" | 0.038" |
| Order Ref.            | ACP205 | ACP206 | ACP207 | ACP208 |



### Straight open tip ureteric catheters with 2 side eyes

| Ch/Fr                 | 3      | 4      | 4.8    | 6      | 7      |
|-----------------------|--------|--------|--------|--------|--------|
| Guidewire recommended | 0.021" | 0.035" | 0.038" | 0.038" | 0.038" |
| Order Ref.            | ACP303 | ACP304 | ACP305 | ACP306 | ACP307 |




### Open coude tip ureteric catheters with 2 eyes

| Ch/Fr                 | 3      | 4      | 4.8    | 6      | 7      |
|-----------------------|--------|--------|--------|--------|--------|
| Guidewire recommended | 0.021" | 0.035" | 0.038" | 0.038" | 0.038" |
| Order Ref.            | ACP403 | ACP404 | ACP405 | ACP406 | ACP407 |



### Open coude tip ureteric catheters without side eye

| Ch/Fr                 | 3      | 4      | 4.8    | 6      | 7      |
|-----------------------|--------|--------|--------|--------|--------|
| Guidewire recommended | 0.021" | 0.035" | 0.038" | 0.038" | 0.038" |
| Order Ref.            | ACP603 | ACP604 | ACP605 | ACP606 | ACP607 |



## 1 • ENDOUROLOGY


### 1.9. Ureteral catheters

# Interventional catheters

## Floppy tip hydro-coated ureteric catheter

The flexible open tip is designed for smooth passage through the ureter.  
The hydrophilic coating over 15cm makes catheter insertion and advancement easier.  
The grey part at the distal end improves the visualisation of the catheter tip.

- Each catheter is supplied with a polyester stylet and a Luer syringe connector.
- Material: Neoplex.
- Length: 70cm.
- Catheters are individually packed.
- **Box of 5** • Sterile.



### Straight open ureteric catheters without side eye

| Ch/Fr                 | 5      | 6      | 7      |
|-----------------------|--------|--------|--------|
| Guidewire recommended | 0.035" | 0.035" | 0.035" |
| Order Ref.            | ACR205 | ACR206 | ACR207 |

## Chevassu catheters for Retrograde UreteroPyelography (RUP)

The shaped bulb of Chevassu catheters allows effective injection of contrast medium.

- Each catheter is supplied with a metal stylet and a Luer syringe connector.
- Material: Neoplex.
- 73cm for AC5C07 and 70cm for AC59xx.
- Catheters are individually packed.
- **Box of 5** • Sterile.

### Chevassu ureteric catheters compatible with a 0.035" guidewire

|            |        |
|------------|--------|
| Bulb Ch/Fr | 7      |
| Body Ch/Fr | 5      |
| Order Ref. | AC5C07 |



### Chevassu ureteric catheters


| Bulb Ch/Fr | 4      | 5      | 6      | 7      | 8      | 9      | 10     | 11     | 12     |
|------------|--------|--------|--------|--------|--------|--------|--------|--------|--------|
| Body Ch/Fr | 3      | 3      | 4      | 5      | 5      | 6      | 7      | 7      | 7      |
| Order Ref. | AC5904 | AC5905 | AC5906 | AC5907 | AC5908 | AC5909 | AC5910 | AC5911 | AC5912 |



# Ureteral dilators

Used for dilatation of the ureter during ureteroscopic procedures

- Compatible with a 0.038" guidewire.
- Diameter (Ch/Fr) indicated at each tip.
- **Box of 10** • Sterile.



| Order Ref. | Ch/Fr     | Length |
|------------|-----------|--------|
| RBD010     | 8 and 10  | 48cm   |
| RBD014     | 12 and 14 | 48cm   |




2 • PERCUTANEOUS ACCESS

2.1. Percutaneous nephrostomy set

Kolibri® set with J catheter

Complete range of drainage set with J catheter

- **Kit composition:**
  - 2 Chiba needles: 18 and 22G.
  - 1 Seldinger or Schüller J-tip guidewire.
  - 1 set of dilators according to the catheter size with or without working sheath.
  - 1 Vortek J catheter with sleeve (29cm) except for 14Ch/Fr (32cm).
  - 1 Luer-lock stopcock and PVC connector tubing.
- **Box of 1 • Sterile.**




|                            |                                               | Ch/Fr   |         |                   |        |        |        |        |
|----------------------------|-----------------------------------------------|---------|---------|-------------------|--------|--------|--------|--------|
| J catheter                 |                                               | 6       | 7       | 8                 | 9      | 10     | 12     | 14     |
| Seldinger guidewire 0.038" | Simple dilator                                | RJA006* |         | RJA008*<br>RJA108 |        | RJA110 | RJA112 | RJA114 |
|                            | Simple dilators + dilator with working sheath |         |         | RJA208            |        | RJA210 | RJA212 |        |
| Schüller guidewire 0.035"  | Simple dilator                                | RJB006* | RJB007* | RJB108            | RJB109 | RJB110 |        |        |
|                            | Simple dilators + dilator with working sheath |         |         | RJB208            |        | RJB210 |        |        |

\* Catheter without stylet.

# Kolibri® set with Malecot catheter

## Complete range of drainage set with Malecot catheter

- **Kit composition:**
  - 2 Chiba needles: 18 and 22G.
  - 1 Seldinger J-tip guidewire 0.038".
  - 1 set of dilators according to the catheter size with or without working sheath.
  - 1 Vortek Malecot catheter with sleeve (36cm with the connector, except for 8Ch/Fr 34cm).
  - 1 Luer-lock stopcock and PVC connector tubing.
- **Individually packaged** • Sterile.




|                     |                                               | Ch/Fr  |        |        |        |
|---------------------|-----------------------------------------------|--------|--------|--------|--------|
| Malecot catheter    |                                               | 8      | 10     | 12     | 14     |
| Seldinger guidewire | Simple dilator                                | RMA108 | RMA110 | RMA112 | RMA114 |
|                     | Simple dilators + dilator with working sheath | RMA208 | RMA210 | RMA212 | RMA214 |


2 • PERCUTANEOUS ACCESS  
 2.2. Percutaneous nephrostomy accessories

Guidewires  
 Catheter insertion and changing\_


• Box of 5 • Sterile.



| Type              | Tip      | Diameter | Length | Order Ref. |
|-------------------|----------|----------|--------|------------|
| Seldinger (stiff) | "J"      | 0.038"   | 80cm   | AE1138     |
|                   | Straight | 0.035"   | 150cm  | AE0C35     |



| Type     | Tip      | Diameter | Length | Order Ref. |
|----------|----------|----------|--------|------------|
| Schüller | "J"      | 0.035"   | 100cm  | AEJ435     |
|          | Straight | 0.035"   | 100cm  | AE0435     |




| Type        | Tip      | Diameter | Length | Order Ref. |
|-------------|----------|----------|--------|------------|
| Lunderquist | "J"      | 0.035"   | 80cm   | AEJ535     |
|             | Straight | 0.035"   | 80cm   | AE0535     |

## Needles

### Puncture of the renal cavities

- Box of 5 • Sterile.

| Order Ref. | Tip   | Diameter |
|------------|-------|----------|
| RAC018     | Chiba | 18G      |
| RAC022     | Chiba | 22G      |



Length: 20cm.


## Dilators

### Tract preparation simple or splittable working sheath to protect the parenchyma and make withdrawal easier

- Box of 5 • Sterile.



|            | Simple dilator |        |        |          |                  |
|------------|----------------|--------|--------|----------|------------------|
| Ch/Fr      | 8              | 10     | 12     | 6, 8, 10 | 6, 8, 10, 12, 14 |
| Order Ref. | RBA008         | RBA010 | RBA012 | RBB010   | RBB014           |



|            | Simple dilator and dilator with working splittable sheath |                     |         |                     |
|------------|-----------------------------------------------------------|---------------------|---------|---------------------|
|            | Simple                                                    | With working sheath | Simple  | With working sheath |
| Ch/Fr      | 6                                                         | 8                   | 6 and 8 | 10                  |
| Order Ref. | RBC008                                                    |                     | RBC010  |                     |

## 2 • PERCUTANEOUS ACCESS

### 2.2. Percutaneous nephrostomy accessories

# J and Malecot catheters

For percutaneous drainage of the renal cavities

- Box of 5 • Sterile.

## J catheters

J catheters are supplied with sleeve + 1 Luer-lock stopcock and PVC connector tubing.

|            | Without stylet |        |        | With stylet |        |        |        |        |
|------------|----------------|--------|--------|-------------|--------|--------|--------|--------|
| Ch/Fr      | 6              | 7      | 8      | 8           | 9      | 10     | 12     | 14     |
| Order Ref. | RCJ006         | RCJ007 | RCJ008 | RCJ108      | RCJ109 | RCJ110 | RCJ112 | RCJ114 |

Length: without connector: 27cm except for RCJ114: 30cm.

## Malecot catheters

Malecot catheter is supplied with sleeve + 1 Luer-lock stopcock and PVC connector tubing. Except RCM016/18/20

|            | With stylet |        |        |        |        |        |        |
|------------|-------------|--------|--------|--------|--------|--------|--------|
| Ch/Fr      | 8           | 10     | 12     | 14     | 16     | 18     | 20     |
| Order Ref. | RCM008      | RCM010 | RCM012 | RCM014 | RCM016 | RCM018 | RCM020 |

Length: without connector 36cm except 34cm for RCM008 and RCM016/18/20.

# Balloon catheters in silicone

For increased comfort during drainage

- All silicone for biocompatibility and patient comfort.
- Radiopaque line helps to see product on X-Ray.
- Box of 10 • Sterile.

| Ch/Fr                                            | 10     | 12     | 14     | 16     | 18     | 20     |
|--------------------------------------------------|--------|--------|--------|--------|--------|--------|
| Catheter with 2 eyes above the balloon           | AJ6610 | AJ6612 | AJ6614 | AJ6616 | AJ6618 | AJ6620 |
| Catheter with 2 eyes above and under the balloon | AJ6710 | AJ6712 | AJ6714 | AJ6716 | AJ6718 | AJ6720 |

Length: 40cm.

# Connection

- Box of 10 • Sterile.

|                                                     |
|-----------------------------------------------------|
| Removable Luer-lock and universal urine bag adapter |
| REA000                                              |



## Set with J catheter

Percutaneous set. Each component of the kit is supplied in the same peel pouch\_

- **Set composition:**
  - Chiba needle 18G.
  - Seldinger guidewire 0.038".
  - 1 set of dilators depending on the catheter size.
  - Vortek J catheter with stylet (except RJE106).
  - Length:29cm (except RJE114: 32.5cm).
- **Box of 5 • Sterile.**

| Ch/Fr      | 6      | 8      | 10     | 12     | 14     |
|------------|--------|--------|--------|--------|--------|
| Order Ref. | RJE106 | RJE108 | RJE110 | RJE112 | RJE114 |

## 2 • PERCUTANEOUS ACCESS

### 2.4. Suprapubic drainage

# Cystodrain® integral Set

## Suprapubic drainage sets with Silicone Integrated balloon catheter

| Punction set                       | Light punction set                 | Dilation set                       | Replacement set                    |
|------------------------------------|------------------------------------|------------------------------------|------------------------------------|
| Silicone Integral balloon catheter | Silicone Integral balloon catheter | Silicone Integral balloon catheter | Silicone Integral balloon catheter |
| Splittable trocar (12cm)           | Splittable trocar (12cm)           | Dilator                            | Adhesive tape                      |
| Locking sleeve                     |                                    | PTFE Guidewire (90 cm)             | Spigot                             |
| Adhesive tape                      |                                    | Adhesive tape                      | Optional PTFE guidewire (90 cm)    |
| Spigot                             |                                    | Spigot                             |                                    |
| Scalpel                            |                                    |                                    |                                    |
| Box of 5 • Sterile                 | Box of 5 • Sterile                 | Box of 5 • Sterile                 | Box of 1 or 5 • Sterile            |

Silicone Integral balloon catheter: Channels: 2 – Eyes: 2 – Balloon size: 3 (Ch.10) / 5 ml – PTFE Guidewire: 0.038" (0.97mm).



### Punction sets references

|                               |        |        |        |
|-------------------------------|--------|--------|--------|
| Trocar Ø int/ext Ch,Fr        | 12/14  | 13/16  | 16/17  |
| Catheter Ø Ch/Fr              | 10     | 12     | 14     |
| Punction set Order Ref.       | CD1010 | CD1012 | CD1014 |
| Light punction set Order Ref. | CD1210 | CD1212 | CD1214 |

### Dilatation set references

|                |        |        |
|----------------|--------|--------|
| Catheter Ch/Fr | 12     | 14     |
| Set Order Ref. | CD1112 | CD1114 |

### Simple replacement sets

|                                                |        |        |        |        |        |        |
|------------------------------------------------|--------|--------|--------|--------|--------|--------|
| Catheter Ch/Fr                                 | 10     | 12     | 14     | 16     | 18     | 20     |
| Replacement set Order Ref. (with guidewire)    | CD1A10 | CD1A12 | CD1A14 | CD1A16 |        |        |
| Replacement set Order Ref. (without guidewire) | CD1B10 | CD1B12 | CD1B14 | CD1B16 | CD1B18 | CD1B20 |
|                                                | CD1C10 | CD1C12 | CD1C14 | CD1C16 | CD1C18 | CD1C20 |

Box of 1

Box of 5

## Supraflow


The benefits of suprapubic puncture  
with the simplicity of a balloon catheter\_

- Complete set.
- Silicone Foley catheter.
- A peelable sheath trocar.
- Scalpel.
- **Box of 5** • Sterile.

### Complete set references

|                     |         |          |
|---------------------|---------|----------|
| Trocar Ø int/ext mm | 5.5/6.8 | 6.75/8.0 |
| Catheter Ø Ch/Fr    | 12      | 16       |
| Ballon size (ml)    | 10      | 15       |
| Set Order Ref.      | AJ9212  | AJ9216   |

Silicone Foley catheter additional info: Length: 40cm – Channels: 2 – Eyes: 2.





## 2 • PERCUTANEOUS ACCESS

### 2.5. Urinary Diversion

# Urinary diversion

- Individually packaged, except for AC41xx and AC65xx (box of 5) • Sterile.



| Length: 155cm • Material: PVC |        |        |        |        |        |        |        |        |
|-------------------------------|--------|--------|--------|--------|--------|--------|--------|--------|
| Ch/Fr                         | 5      | 6      | 7      | 8      | 9      | 10     | 11     | 12     |
| Order Ref.                    | AC4105 | AC4106 | AC4107 | AC4108 | AC4109 | AC4110 | AC4111 | AC4112 |

Catheter for uretero-sigmoidostomy, 6 eyes, straight proximal end, radiopaque line.



| Length: 120cm • Material: PVC |        |        |        |        |
|-------------------------------|--------|--------|--------|--------|
| Ch/Fr                         | 8      | 10     | 12     | 14     |
| Order Ref.                    | AC4708 | AC4710 | AC4712 | AC4714 |

Catheter for entero-cystoplasty, with eyes, funnel proximal end, radiopaque line.



| Length: 45cm • Material: Silicone |        |        |        |        |        |        |        |        |
|-----------------------------------|--------|--------|--------|--------|--------|--------|--------|--------|
| Ch/Fr                             | 8      | 10     | 12     | 14     | 16     | 18     | 20     | 22     |
| Order Ref.                        | AC6708 | AC6710 | AC6712 | AC6714 | AC6716 | AC6718 | AC6720 | AC6722 |

Ureterostomy catheter, 5 eyes, bevel tip, funnel proximal end, radiopaque line. Round flange.



| Length: 45cm • Material: Silicone |        |        |        |        |        |        |        |
|-----------------------------------|--------|--------|--------|--------|--------|--------|--------|
| Ch/Fr                             | 8      | 10     | 12     | 14     | 16     | 20     | 22     |
| Order Ref.                        | AC6808 | AC6810 | AC6812 | AC6814 | AC6816 | AC6820 | AC6822 |

Ureterostomy catheter, 5 eyes, bevel tip, funnel proximal end, radiopaque line. Oval flange.



| Length: 45cm • Material: Silicone |        |        |        |
|-----------------------------------|--------|--------|--------|
| Ch/Fr                             | 8      | 10     | 12     |
| Order Ref.                        | AC7008 | AC7010 | AC7012 |

Ureterostomy catheter, 5 eyes, bevel tip, straight proximal end, radiopaque line. Round flange.



| Total Length: 52cm • Material: Silicone |        |        |        |
|-----------------------------------------|--------|--------|--------|
| Ch/Fr                                   | 8      | 10     | 12     |
| Order Ref.                              | BJ3008 | BJ3010 | BJ3012 |

Ureterostomy and entero-cystoplasty catheter, single loop, open and tapered distal tip, radiopaque. Oval flange.

# Connectors and spigot

- Box of 10 • Sterile.
- CP18200000: Box of 50 • Sterile.

## Connectors for tube



| Material: Polypropylene |        |        |          |
|-------------------------|--------|--------|----------|
| External diameter (mm)  | 4-6-4  | 7-11-7 | 10-14-10 |
| Internal diameter (mm)  | 3-3    | 5-5    | 8-8      |
| Order Ref.              | PB0101 | PB0102 | PB0103   |



| Material: Polypropylene |            |             |              |
|-------------------------|------------|-------------|--------------|
| External diameter (mm)  | 4-6 / 7-11 | 4-6 / 14-10 | 7-11 / 14-10 |
| Internal diameter (mm)  | 3-5        | 3-8         | 5-8          |
| Order Ref.              | PB0201     | PB0202      | PB0203       |



| Material: Polypropylene |         |         |          |           |
|-------------------------|---------|---------|----------|-----------|
| External diameter (mm)  | 3 x 4-5 | 3 x 6-8 | 3 x 9-11 | 3 x 11-13 |
| Internal diameter (mm)  | 2-2-2   | 4-4-4   | 7-7-7    | 8-8-8     |
| Order Ref.              | PB1201  | PB1203  | PB1204   | PB1205    |



## Tuohy Borst Adapter

|                                                 |        |
|-------------------------------------------------|--------|
| Material:<br>Polycarbonate • Acrylic • Silicone |        |
| Order Ref.                                      | AK2000 |

## Spigot



|               |            |
|---------------|------------|
| Material: ABS |            |
| Order Ref.    | CP18200000 |

## 3 • CONNECTORS AND SPIGOT

### 3.1. Connectors and Spigot

#### Connectors for catheters



| Material: Thermoplastic Elastomer Styrene (TPE) |        |        |
|-------------------------------------------------|--------|--------|
| Ch/Fr                                           | 3-10   | 9-12   |
| Order Ref.                                      | AK4010 | AK4012 |



| Material: PVC |        |
|---------------|--------|
| Order Ref.    | AK5300 |



| Material : Thermoplastic Elastomer Styrene (Sheath)<br>Acrylic copolymer (Luer male) • Thermoplastic (Luer female) |        |
|--------------------------------------------------------------------------------------------------------------------|--------|
| Order Ref.                                                                                                         | AK4100 |



| Material: Thermoplastic Elastomer Styrene (Sheath)<br>Acrylic copolymer (Luer male) • Thermoplastic and Silicone (Luer connector) |        |        |
|-----------------------------------------------------------------------------------------------------------------------------------|--------|--------|
| Ch/Fr                                                                                                                             | 3-8    | 9-12   |
| Order Ref.                                                                                                                        | AK4200 | AK4201 |



| Material: Thermoplastic Elastomer Styrene (Sheath) • Acrylic copolymer (Luer male) |        |
|------------------------------------------------------------------------------------|--------|
| Order Ref.                                                                         | AK4300 |



| Material: Polyamide • Silicone |        |        |
|--------------------------------|--------|--------|
| Ch/Fr                          | 3-8    | 9-12   |
| Order Ref.                     | AK3400 | AK3401 |

# Dufour tip



- Length: 41cm.
- Box of 5 • Sterile.



Silicone X-Flow  
Hydro X-Flow



Silicone X-Flow  
Short Tip

|        |              |                                  | Ch/Fr  |        |        |        |
|--------|--------------|----------------------------------|--------|--------|--------|--------|
| Ways   | Balloon size | Material                         | 18     | 20     | 22     | 24     |
| 3 Ways | 30ml         | X-Flow Silicone Dufour short tip | AB6S18 | AB6S20 | AB6S22 | AB6S24 |
|        | 50ml         | X-Flow Silicone Dufour short tip | AB6R18 | AB6R20 | AB6R22 | AB6R24 |
|        |              | X-Flow Silicone Dufour tip       | AB6318 | AB6320 | AB6322 | AB6324 |
|        |              | Hydro X-Flow Silicone Dufour tip | XB6318 | XB6320 | XB6322 | XB6324 |
|        | 80ml         | Hydro X-Flow Silicone Dufour tip | XB6L18 | XB6L20 | XB6L22 | XB6L24 |
| 2 Ways | 50ml         | X-Flow Silicone Dufour tip       | AB6418 | AB6420 | AB6422 | AB6424 |

4 • LOWER URINARY TRACT  
 4.1. Prostatic catheters

Delinotte tip  
 Coudé tip\_

- Length: 41cm.
- Box of 5 • Sterile.



Silicone X-Flow  
 Hydro X-Flow

|        |              |                                     | Ch/Fr  |        |        |        |
|--------|--------------|-------------------------------------|--------|--------|--------|--------|
| Ways   | Balloon size | Material                            | 18     | 20     | 22     | 24     |
| 3 Ways | 50ml         | X-Flow Silicone Delinotte tip       | AB6218 | AB6220 | AB6222 | AB6224 |
|        |              | Hydro X-Flow Silicone Delinotte tip | XB6218 | XB6220 | XB6222 | XB6224 |
|        | 80ml         | Hydro X-Flow Silicone Delinotte tip | XB6N18 | XB6N20 | XB6N22 | XB6N24 |

# Couvelaire tip

## Whistle tip\_

- Length: 41cm.
- Box of 5 • Sterile.



Silicone X-Flow  
Hydro X-Flow

|        |              |                                      | Ch/Fr  |        |        |        |
|--------|--------------|--------------------------------------|--------|--------|--------|--------|
| Ways   | Balloon size | Material                             | 18     | 20     | 22     | 24     |
| 3 Ways | 50ml         | X-Flow Silicone Couvelaire tip       | AB6118 | AB6120 | AB6122 | AB6124 |
|        |              | Hydro X-Flow Silicone Couvelaire tip | XB6118 | XB6120 | XB6122 | XB6124 |
|        | 80ml         | Hydro X-Flow Silicone Couvelaire tip | XB6M18 | XB6M20 | XB6M22 | XB6M24 |
| 2 Ways | 50ml         | X-Flow Silicone Couvelaire tip       | AB6518 | AB6520 | AB6522 | AB6524 |

4 • LOWER URINARY TRACT

4.1. Prostatic catheters

# Straight tip



- Length: 41cm.
- **Box of 5** • Sterile.



Silicone  
X-Flow



Silicone short tip  
X-Flow



Over guidewire

| Ways   | Balloon size | Material                                  | Ch/Fr  |        |        |        |
|--------|--------------|-------------------------------------------|--------|--------|--------|--------|
|        |              |                                           | 18     | 20     | 22     | 24     |
| 3 Ways | 50ml         | X-Flow Silicone straight tip              | AB6018 | AB6020 | AB6022 | AB6024 |
|        |              | X-Flow Silicone straight short tip        | AB6H18 | AB6H20 | AB6H22 | AB6H24 |
|        |              | X-Flow Silicone over the guide short tip* | AB6G18 | AB6G20 | AB6G22 | AB6G24 |

\* Compatible with guidewire.

## Freudenberg Stylet

- Intended to stiffen or elbow a prostatic catheter in case of urethral insertion difficulty.
- **Individually packaged** • Sterile • Reusable.



| Freudenberg stylet |        |
|--------------------|--------|
| Order Ref.         | AL0000 |



4 • LOWER URINARY TRACT

4.1. Prostatic catheters

Catheters

With 30ml balloon



- Sterile.



Silicone  
X-Flow



Silicone short tip  
X-Flow

With 30ml balloon • Box of 5

- Length 41cm.

|        |              |                    |                 | Ch/Fr  |        |        |        |
|--------|--------------|--------------------|-----------------|--------|--------|--------|--------|
| Ways   | Balloon size | Tip                | Material        | 18     | 20     | 22     | 24     |
| 3 Ways | 30ml         | Dufour             | X-Flow Silicone | AB6A18 | AB6A20 | AB6A22 | AB6A24 |
|        |              | Delinotte          | X-Flow Silicone | AB6E18 | AB6E20 | AB6E22 | AB6E24 |
|        |              | Couvelaire         | X-Flow Silicone | AB6B18 | AB6B20 | AB6B22 | AB6B24 |
|        |              | Straight           | X-Flow Silicone | AB6C18 | AB6C20 | AB6C22 | AB6C24 |
|        |              | Straight short tip | X-Flow Silicone | AB6J18 | AB6J20 | AB6J22 | AB6J24 |

# ABVac

## Bladder evacuator\_

Intended to be used for evacuating irrigation solutions, blood clots, tissue collection and/or straining functions during prostate surgery or bladder surgery.

### Main Features

- Single Nozzle adaptor for Storz, Wolf, Olympus and ACMI Resectoscopes.
- Single hand use for an excellent eye - hand coordination.
- Excellent throw and suction.

| Order Ref. | Description       | Quantity |
|------------|-------------------|----------|
| 000PABVAC  | Bladder evacuator | 10 units |



4 • LOWER URINARY TRACT

4.3. Bladder injection

BoNee<sup>®</sup>

Needle for bladder injection



- For Injection inside the bladder muscle.
  - Short tip (4mm).
  - 22G needle
  - 5Ch/Fr tube.
  - Protection cap on the needle tip.
- Pointer stop: Reduce the bladder perforation risk.
  - Flexible body.
  - **Box of 1 or 5** • Sterile.

| Tube length | Use                                     | Box of 1 | Box of 5 |
|-------------|-----------------------------------------|----------|----------|
| 35cm        | With rigid cystoscope                   | NBI035   | BBI035   |
| 70cm        | With both flexible and rigid cystoscope | NBI070   | BBI070   |

22G needle, 5 Ch/Fr tube, 4mm tip.



## Bougies

- Dilation of the meatus before endoscopy.
- **Individually packaged** • Gradual dilation of strictures • Sterile.

Dilation of the meatus before endoscopy and gradual dilation of strictures are applicable for all bougies

### For routine urethral dilation

| Neoplex • Straight / tapered olive tip bougie - Length: 34cm |        |        |        |        |        |        |        |        |
|--------------------------------------------------------------|--------|--------|--------|--------|--------|--------|--------|--------|
| Ch/Fr                                                        | 6      | 7      | 8      | 9      | 10     | 12     | 14     | 16     |
| Order Ref.                                                   | AG5306 | AG5307 | AG5308 | AG5309 | AG5310 | AG5312 | AG5314 | AG5316 |
| Ch/Fr                                                        | 18     | 20     | 22     | 24     | 26     | 28     | 30     |        |
| Order Ref.                                                   | AG5318 | AG5320 | AG5322 | AG5324 | AG5326 | AG5328 | AG5330 |        |



## Bougies metric thread

For difficult-to-cross strictures (Order Ref. AG51)\_  
 Philips connectable bougies (Order Ref. AG54 only)\_

- **Individually packaged** • Sterile.

### Filiform bougies

|          |        | Male metric thread • Neoplex |        |        |        |
|----------|--------|------------------------------|--------|--------|--------|
|          |        | Ch/Fr                        |        |        |        |
| Tip      | Length | 4                            | 5      | 6      | 7      |
| Straight | 30cm   | AG5104                       | AG5105 | AG5106 | AG5107 |

### Philips Bougies

| Female metric thread • Neoplex - Length: 34cm |        |        |        |        |        |        |        |        |        |
|-----------------------------------------------|--------|--------|--------|--------|--------|--------|--------|--------|--------|
| Ch/Fr                                         | 8      | 10     | 12     | 14     | 16     | 18     | 20     | 22     | 24     |
| Order Ref.                                    | AG5408 | AG5410 | AG5412 | AG5414 | AG5416 | AG5418 | AG5420 | AG5422 | AG5424 |





# Foley catheters

## Silicone catheters



All lengths are indicative.



Tiemann



Tiemann+



Over guidewire



Folysil X-Tra



Neobladder 3 ways

• Box of 5 • Sterile.

| Type               | Balloon | Total length | Number of eyes | Ch/Fr  |        |        |        |        |        |        |
|--------------------|---------|--------------|----------------|--------|--------|--------|--------|--------|--------|--------|
|                    |         |              |                | 12     | 14     | 16     | 18     | 20     | 22     | 24     |
| Straight           | 10ml    | 40cm         | 2              | AW6112 | AW6114 | AW6116 | AW6118 | AW6120 | AW6122 | AW6124 |
|                    | 10ml    | 41cm         | 2              | AA6112 | AA6114 | AA6116 | AA6118 | AA6120 | AA6122 | AA6124 |
|                    | 30ml    | 41cm         | 2              |        |        |        | AA6C18 | AA6C20 | AA6C22 | AA6C24 |
| Female Straight    | 10ml    | 25cm         | 2              | AA6512 | AA6514 |        |        |        |        |        |
|                    | 15ml    | 25cm         | 2              |        |        | AA6516 | AA6518 | AA6520 | AA6522 | AA6524 |
| Over the guidewire | 10ml    | 40cm         | 2 + tip        | AA6412 | AA6414 |        |        |        |        |        |
|                    | 15ml    | 40cm         | 2 + tip        |        |        | AA6416 | AA6418 | AA6420 | AA6422 | AA6424 |
| Straight Grooved   | 10ml    | 40cm         | 2              | AA6612 | AA6614 |        |        |        |        |        |
|                    | 15ml    | 40cm         | 2              |        |        | AA6616 | AA6618 | AA6620 | AA6622 | AA6624 |
| Tiemann            | 10ml    | 40cm         | 1              | AA6312 | AA6314 |        |        |        |        |        |
|                    | 15ml    | 40cm         | 1              |        |        | AA6316 | AA6318 | AA6320 | AA6322 | AA6324 |
|                    | 10ml    | 42cm         | 1              | AW6312 | AW6314 | AW6316 | AW6318 | AW6320 | AW6322 | AW6324 |
| Tiemann+           | 10ml    | 40cm         | 1              | AA9312 | AA9314 |        |        |        |        |        |
|                    | 15ml    | 40cm         | 1              |        |        | AA9316 | AA9318 | AA9320 | AA9322 |        |
| Neobladder 3 ways  | 30ml    | 40cm         | 6 + tip        |        |        |        |        | AA6820 | AA6822 |        |
|                    | 30ml    | 40cm         | 4 + 2** + tip  |        |        |        | AA6918 | AA6920 | AA6922 |        |

\*\* 2 eyes under the balloon.

10ml balloon

\* Delivered without introducer.

## 4 • LOWER URINARY TRACT

### 4.5. Foley catheters

All lengths are indicative.

#### Paediatric foley catheters, up to 30 days, except LT\* use

• Box of 5 • Sterile.

| Type                             | Balloon | Length | Ch/Fr  |        |                 |
|----------------------------------|---------|--------|--------|--------|-----------------|
|                                  |         |        | 6      | 8      | 10              |
| Straight<br>(with<br>introducer) | 1.5ml   | 26cm   | AA6106 |        |                 |
|                                  | 1.5ml   | 31cm   | AW6106 |        |                 |
|                                  | 3ml     | 26cm   |        | AA6108 |                 |
|                                  | 3ml     | 31cm   |        | AW6108 | AA6110 / AW6110 |
| Over the<br>guidewire            | 1.5ml   | 26cm   | AA6406 |        |                 |
|                                  | 3ml     | 26cm   |        | AA6408 |                 |
|                                  | 3ml     | 31cm   |        |        | AA6410          |
| Tiemann                          | 3ml     | 27cm   |        | AA6308 |                 |
|                                  | 3ml     | 32cm   |        |        | AA6310          |
| *Long Term                       | 1.5ml   | 26cm   | AA7106 |        |                 |
|                                  | 3ml     | 26cm   |        | AA7108 |                 |
|                                  | 3ml     | 31cm   |        |        | AA7110          |

#### Folysil LT, for long term use (up to 12 weeks)

• Box of 5 • Sterile.

| Type              | Balloon | Length | Number<br>of eyes | Ch/Fr  |        |        |        |        |        |
|-------------------|---------|--------|-------------------|--------|--------|--------|--------|--------|--------|
|                   |         |        |                   | 12     | 14     | 16     | 18     | 20     | 22     |
| Straight          | 10ml    | 41cm   | 2                 | AA7112 | AA7114 | AA7116 | AA7118 | AA7120 | AA7122 |
| Female            | 10ml    | 25cm   | 2                 | AA7512 | AA7514 | AA7516 | AA7518 |        |        |
| Over<br>guidewire | 10ml    | 40cm   | 2 + tip           | AA7412 | AA7414 | AA7416 | AA7418 | AA7420 |        |

#### Folysil X-tra, with prefilled syringe for long term use (up to 12 weeks)

• Box of 5 • Sterile.

| Type           | Balloon | Length | Number<br>of eyes | Ch/Fr  |        |        |        |        |        |
|----------------|---------|--------|-------------------|--------|--------|--------|--------|--------|--------|
|                |         |        |                   | 12     | 14     | 16     | 18     | 20     | 22     |
| X-Tra Straight | 10ml    | 41cm   | 2                 | AA8112 | AA8114 | AA8116 | AA8118 | AA8120 | AA8122 |
| X-Tra Female   | 10ml    | 25cm   | 2                 | AA8512 | AA8514 | AA8516 | AA8518 |        |        |

#### Folysil X-Tra single pack, with prefilled syringe and empty syringe, for long term use (up to 12 weeks)

• Box of 1 • Sterile.

| Type                    | Balloon | Length | Number<br>of eyes | Ch/Fr  |        |        |        |        |        |
|-------------------------|---------|--------|-------------------|--------|--------|--------|--------|--------|--------|
|                         |         |        |                   | 12     | 14     | 16     | 18     | 20     | 22     |
| X-Tra Straight          | 10ml    | 41cm   | 2                 | AA8A12 | AA8A14 | AA8A16 | AA8A18 | AA8A20 | AA8A22 |
| X-Tra Female            | 10ml    | 25cm   | 2                 | AA8B12 | AA8B14 | AA8B16 | AA8B18 |        |        |
| X-Tra over<br>guidewire | 10ml    | 40cm   | 2 + tip           | AA8C12 | AA8C14 | AA8C16 | AA8C18 |        |        |

## Catheter Valve

Intended to be connected to an indwelling catheter, to drain the urine\_

• Box of 10 • Sterile.

Order Ref.

CV3808



## Urethral catheters without balloon

• Box of 5 • Sterile.



Couvelaire (Neoplex)

### Straight whistle tip catheters

|            |          |        |      | Ch/Fr  |        |        |        |        |        |        |        |        |        |        |
|------------|----------|--------|------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|
| Type       | Material | Length | Eyes | 6      | 8      | 10     | 12     | 14     | 16     | 18     | 20     | 22     | 24     | 26     |
| Couvelaire | Neoplex  | 40cm   | 2    | AD5D06 | AD5D08 | AD5D10 | AD5D12 | AD5D14 | AD5D16 | AD5D18 | AD5D20 | AD5D22 | AD5D24 | AD5D26 |



Mercier

### Coude tip catheters

|         |          |        |      | Ch/Fr  |        |        |        |        |        |        |        |        |        |        |
|---------|----------|--------|------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|
| Type    | Material | Length | Eyes | 6      | 8      | 10     | 12     | 14     | 16     | 18     | 20     | 22     | 24     | 26     |
| Mercier | Neoplex  | 40cm   | 2    | ADN306 | ADN308 | ADN310 | ADN312 | ADN314 | ADN316 | ADN318 | ADN320 | ADN322 | ADN324 | ADN326 |



# Prostatic Catheters tip terminology



## Dufour tip

Coudé whistle tip



## Delinotte tip

Coudé cylindrical tip



## Couvelaire tip

Straight whistle tip.



## Straight tip



# Surgery



## SURGERY

### Surgery

- |    |                      |    |
|----|----------------------|----|
| 1. | Laparoscopic surgery | 75 |
| 2. | Surgical drainage    | 76 |

# Elefant®

## Disposable suction-irrigation device for laparoscopy\_



- Single use cannula:
- Ergonomic and lightweight.
- Yankauer-type tip (6 distal holes).
- “E” fiberglass (non-conductive).

| Elefant, without tubing • Box of 10 |        |        |         |
|-------------------------------------|--------|--------|---------|
| Cannula diameter                    | 5mm    | 10mm   | 5mm     |
| Cannula length                      | 35cm   | 35cm   | 45cm    |
| Order Ref.                          | ASP145 | ASP150 | ASP180* |

\* For robotic surgery.

| Elefant, with double irrigation line tubing • Box of 5 |                         |                         |
|--------------------------------------------------------|-------------------------|-------------------------|
| Cannula diameter                                       | 5mm                     | 5mm                     |
| Cannula length                                         | 35cm                    | 45cm                    |
| Tube diameter                                          | ext. 9.5mm<br>int 6.4mm | ext. 9.5mm<br>int 6.4mm |
| Order Ref.                                             | ASP190                  | ASP195*                 |

\* For robotic surgery.



Please contact your Coloplast Territory Manager regarding this product’s availability in your market.

1 • LAPAROSCOPIC SURGERY

Non-traumatic and  
“E” fiberglass  
cannula

Yankauer tip



Elefant cannula made  
of non-conductive fiberglass,  
unlike conventional irrigation devices  
made of stainless steel

| Elefant, with PVC tubing • Box of 5 |                         |                         |                         |
|-------------------------------------|-------------------------|-------------------------|-------------------------|
| Cannula diameter                    | 5mm                     | 10mm                    | 5mm                     |
| Cannula length                      | 35cm                    | 35cm                    | 45cm                    |
| Tube diameter                       | ext. 9.5mm<br>int 6.4mm | ext. 9.5mm<br>int 6.4mm | ext. 9.5mm<br>int 6.4mm |
| Tube length                         | 3m                      | 3m                      | 3m                      |
| Order Ref.                          | ASP165                  | ASP170                  | ASP185*                 |

\* For robotic surgery.



| Elefant, with marked cannula<br>Box of 5 |                         |
|------------------------------------------|-------------------------|
| 5mm                                      | 5mm                     |
| 35cm                                     | 45cm                    |
| ext. 9.5mm<br>int 6.4mm                  | ext. 9.5mm<br>int 6.4mm |
| 3m                                       | 3m                      |
| ASP265                                   | ASP285*                 |

\* For robotic surgery.

# Suction drainage

- Ch/Fr 9-24: **box of 10** • Ch/Fr 27-36: **box of 5** • Sterile.

| Multiperforated tubular drains • Material: Silicone • Length: 50cm |        |        |        |        |        |        |        |        |        |        |
|--------------------------------------------------------------------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|
| Ch/Fr                                                              | 9      | 12     | 15     | 18     | 21     | 24     | 27     | 30     | 33     | 36     |
| Order Ref.                                                         | GB6009 | GB6012 | GB6015 | GB6018 | GB6021 | GB6024 | GB6027 | GB6030 | GB6033 | GB6036 |

Radiopaque line.

Ch/Fr 9-24: box of 10 • Ch/Fr 27-36: box of 5.



## 2 • SURGICAL DRAINAGE


# Simple drainage

- GA1035, GA1245, GA62xx, GA69xx Ch/Fr 8 to 25: **box of 10**
- GA6106, GA69xx Ch/Fr 30 to 35: **box of 5**
- Sterile.

| Tubular drains • Material: Silicone • Length: 30cm |        |        |        |        |        |        |        |        |        |        |
|----------------------------------------------------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|
| Ch/Fr                                              | 8      | 10     | 12     | 14     | 16     | 18     | 20     | 22     | 25     | 30     |
| Order Ref.                                         | GA6208 | GA6210 | GA6212 | GA6214 | GA6216 | GA6218 | GA6220 | GA6222 | GA6225 | GA6230 |

| Ch/Fr      | 35     | 40     |
|------------|--------|--------|
| Order Ref. | GA6235 | GA6240 |

Radiopaque line.





| Tubular drains • Material: Silicone • Length: 50cm |        |        |        |        |        |        |        |        |        |        |
|----------------------------------------------------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|
| Ch/Fr                                              | 8      | 10     | 12     | 14     | 16     | 18     | 20     | 25     | 30     | 35     |
| Order Ref.                                         | GA6908 | GA6910 | GA6912 | GA6914 | GA6916 | GA6918 | GA6920 | GA6925 | GA6930 | GA6935 |



| Multitubular drain • Material: Radiopaque Silicone • Length: 40cm |        |
|-------------------------------------------------------------------|--------|
| Width (mm)                                                        | 35     |
| Order Ref.                                                        | GA1035 |

| Multitubular drain • Material: Silicone • Length: 40cm |        |
|--------------------------------------------------------|--------|
| Width (mm)                                             | 45     |
| Order Ref.                                             | GA1245 |




| Delbet drains • Material: Silicone • Length: 36cm |        |
|---------------------------------------------------|--------|
| Width (cm)                                        | 6      |
| Order Ref.                                        | GA6106 |

Box of 5. Radiopaque

# Biliary drainage

- Box of 5 • Sterile.




| Transcystic drains • Material: PVC • Length: 41.5cm |        |        |        |
|-----------------------------------------------------|--------|--------|--------|
| Bulb Ø Ch/Fr                                        | 7      | 8      | 10     |
| Tube Ø Ch/Fr                                        | 5      | 6      | 8      |
| Order Ref.                                          | GD4105 | GD4106 | GD4108 |

Pedinielli drains with malleable stylet, whistle tip, 2 eyes.  
Radiopaque line.



| Transcystic drains • Material: PVC • Length: 154cm |        |        |        |        |
|----------------------------------------------------|--------|--------|--------|--------|
| Bulb Ø Ch/Fr                                       | 7      | 8      | 10     | 12     |
| Tube Ø Ch/Fr                                       | 5      | 6      | 8      | 10     |
| Order Ref.                                         | GD4205 | GD4206 | GD4208 | GD4210 |

Escat drains, whistle tip, 2 eyes, with 3-way stopcock.  
Radiopaque line.



| Cœliodrain® • Material: Radiopaque Neoplex • Length: approximately 64cm |        |
|-------------------------------------------------------------------------|--------|
| Ch/Fr                                                                   | 5      |
| Order Ref.                                                              | GD4605 |

For cholangiography with malleable stylet and whistle tip.





# Men's Health



MEN'S HEALTH

Men's health

|    |                           |    |
|----|---------------------------|----|
| 1. | Erectile Dysfunction (ED) | 82 |
| 2. | Male Urinary Incontinence | 86 |

## Titan® Touch

### Inflatable penile implant

3-pieces inflatable penile implant with hydrophilic coating and deflate button pump design.

- Touch pump: one touch of the deflate pads allows for easy deflation of the device.
- The Bioflex® material of the cylinders allows for unmatched cylinder strength and rigidity.
- Reduced risk of auto-inflation with Lock-out™ valve on the reservoir.
- Titan HydroVANTAGE™ hydrophilic coating is on all device components, including rear tip extenders and connectors. Physicians have the ability to choose the aqueous solution dip, tailored to the clinical needs of their patients.
- Zero Degree angle proximal base input tubing for ease of surgical placement.
- Soft molded tip designed for suitable anatomic positioning.
- Individually packaged • Sterile.



## Titan® Touch

### Penile prosthesis

| Length (cm)                 | 14     | 16     | 18     | 20     | 22     |
|-----------------------------|--------|--------|--------|--------|--------|
| Order Ref. /<br>Infrapubic  | ES2814 | ES2816 | ES2818 | ES2820 | ES2822 |
| Order Ref. /<br>Penoscrotal | ES2914 | ES2916 | ES2918 | ES2920 | ES2922 |

## Titan® Touch Narrow Base

### For narrow corpora cavernosa

### Penile prosthesis

| Length (cm)                 | 11     | 14     | 16     | 18     |
|-----------------------------|--------|--------|--------|--------|
| Order Ref. /<br>Infrapubic  | EN2811 | EN2814 | EN2816 | EN2818 |
| Order Ref. /<br>Penoscrotal | EN2911 | EN2914 | EN2916 | EN2918 |




## 1 • ERECTILE DYSFUNCTION

# Titan® Cloverleaf Reservoir

With Lock-out anti-reflux valve.

| Reservoir contents | 75 cc  | 125 cc |
|--------------------|--------|--------|
| Order Ref.         | ER8075 | ER8125 |



# True-Lock® Plug

The True-Lock polysulfone plug and sleeve are utilized to close unused tubing on a portion of the TitanR IPP, where deemed appropriate.

| Order Ref. | 9471S |
|------------|-------|
|------------|-------|



# Assembly kits

Standard assembly kit.

| Order Ref. / Standard | 519150 |
|-----------------------|--------|
|-----------------------|--------|

### Contains

- Connector tool (1) to compress connector sleeves
- True-Lock® Connectors (connector sleeve (8), connector handle (2), connector cage (4))
- Keith needle (2)
- Pusher protector (1), for corporotomy closure and proximal insertion
- Barbed syringe tip (3), to fill the device
- Flexible tip needle (2), pink syringe tip to remove air from silicone tubing
- Silicone tubing (6)



Rear tip extenders kit.

| Order Ref. / Standard    | 519310 |
|--------------------------|--------|
| Order Ref. / Narrow Base | 519300 |

To be used with the standard assembly kit.



# XL Cylinder Set

| Cylinders  | 24cm   | 26cm   | 28cm   |
|------------|--------|--------|--------|
| Order Ref. | E17924 | E17926 | E17928 |

Pump to be connected.



# Pumps

| Pumps      | Titan Pump | Titan Touch Pump |
|------------|------------|------------------|
| Order Ref. | 517770     | 517750           |

Titan Touch Pump



# Genesis®

## Malleable penile implant

Malleable penile implant with hydrophilic coating.

- **Kit includes:**
  - One pair of cylinders.
  - 3 pairs of caps (standard, +0.5cm and +1cm).
- Hydrophilic coating.
- **Individually packaged (pairs)** • Sterile.



| Diameter | Length cm |                           | Order Ref. |
|----------|-----------|---------------------------|------------|
|          | mm        | <div>Min.      Max.</div> |            |
| 9.5      |           | <div>14      23</div>     | 519250     |
| 11       |           | <div>16      25</div>     | 519260     |
| 13       |           | <div>18      27</div>     | 519270     |



1 • ERECTILE DYSFUNCTION

# Ancillaries for penile prostheses

## Hooks & dilators\_

**Furlow**  
Non sterile.

Order Ref. QB5555

- Furlow insertion tool



**Brooks**  
Non sterile.

- Brooks dilator



|                   | Ø   |     |      |      |      |      |      |
|-------------------|-----|-----|------|------|------|------|------|
| Order Ref. QB9350 | 8mm | 9mm | 10mm | 11mm | 12mm | 13mm | 14mm |

**Rosselló Dilators**  
Non sterile.

Order Ref. QB5558



**Retractor System without Frame**  
Sterile - Single use.

Order Ref. QB3321

- 6 - Retractor Stays, 12mm Blunt
- 1 - Retractor Stay, 5mm Sharp
- 1 - Penile Elevation Strap
- 2 - Disposable deavers
- 2 - Retractor Stays, 20mm 3-finger Rake



**Galaxy II Surgical Retractor**  
Sterile - Single use.

- Male kit

Order Ref. 000PJUNEGX-2020

- 6 - 12 mm blunt stays
- 1 - 5mm sharp stay
- 2 - 3 finger claw
- 2 - 7mm double hook sharp
- 1 - penile hammock
- 1 - snowman retractor frame
- 1 - catheter clip
- 1 - light/lux clip

## Virtue® male sling

### 4-arm male sling for male stress urinary incontinence\_

- Knitted monofilament polypropylene.
- Bidirectional compression and elevation of the bulbous urethra with 2 transobturator arms and 2 prepubic arms.
- **Individually packaged** • Sterile.



| Order Ref. | Description                              | Sling dimensions |
|------------|------------------------------------------|------------------|
| 50020      | Virtue: 1 sling, 1 single use introducer | 46 x 31cm        |





# Women's Health



WOMEN'S HEALTH

Women's Health

|    |                                   |    |
|----|-----------------------------------|----|
| 1. | Stress Urinary Incontinence (SUI) | 90 |
| 2. | Pelvic Organ Prolapse (POP)       | 92 |


**Altis®**

## Single Incision Sling System

The Altis Single Incision Sling System is indicated for the treatment of female stress urinary incontinence (SUI) resulting from urethral hypermobility and/or intrinsic sphincter deficiency (ISD). It is an implantable, synthetic, knitted, macroporous, low-elasticity, monofilament polypropylene, single incision midurethral sling.



| Order Ref. | Description                         | Mesh size   |            |
|------------|-------------------------------------|-------------|------------|
| 519650     | Altis® single incision sling system | 7.75cm long | 1.1cm wide |




1 • STRESS URINARY INCONTINENCE

# Aris®

## Transobturator Kit

A midurethral sling made from synthetic, knitted, low-elasticity monofilament polypropylene intended to be placed via the transobturator approach. The Aris sling is designed to support the urethra in the surgical treatment of female stress urinary incontinence. Aris Introducers are surgical instruments designed to facilitate the correct placement of the Aris Transobturator Sling.



| Order Ref. | Description                                                     | Mesh size  |
|------------|-----------------------------------------------------------------|------------|
| UR3105     | Aris® Transobturator kit<br>Contents: 1 sling and 3 introducers | 60 x 1.1cm |
| UR3101     | 1 Aris® Transobturator sling                                    | 60 x 1.1cm |



# Restorelle®

## Mesh


Restorelle mesh is a synthetic, macroporous mesh constructed of medical grade knitted, monofilament polypropylene. Restorelle L, XL and Y mesh are implanted transabdominally through an open abdominal, laparoscopic, or robot- assisted surgical approach.

Restorelle L (24cm x 8cm) and Restorelle XL (30cm x 30cm) are intended to provide anatomical support and suspension of the vaginal apex or uterine cervix or uterine isthmus, with or without uterus preservation (sacrocolpopexy/sacrocolpocervicopexy or sacrohysteropexy).



### Restorelle Y

| Order Ref. | Description   | Size     | Tail size |
|------------|---------------|----------|-----------|
| 501420     | Restorelle® Y | 24 x 4cm | 12cm      |
| 501430     | Restorelle® Y | 27 x 4cm | 15cm      |



### Flat Mesh

| Order Ref. | Description    | Dimensions |
|------------|----------------|------------|
| 501440     | Restorelle® L  | 24 x 8cm   |
| 501330     | Restorelle® XL | 30 x 30cm  |



# Guarantee and conditions of use

## Coloplast A/S - Guarantee and conditions of use

Without prejudice to your legal guarantee against all the consequences of defects or hidden defects, our commercial guarantee shall apply under the following conditions:

The duration of the normal guarantee of the products against defects in manufacturing workmanship and against defects in materials shall be a maximum of three months from the date of delivery (the "Warranty Period"), with the onus being on the Purchaser to prove the said defects.

Coloplast's guarantee shall be strictly limited to the products supplied and shall have as its sole effect the repair or replacement at the Company's expense and in its workshops of all those products rendered unusable as a result of defects. This is subject to Coloplast receiving written notice of the defect within 10 days from the date the defect was or should have been discovered by the Purchaser, within the Warranty Period. The return of the products must be done within the deadline communicated by Coloplast. Any return not respecting these conditions, cannot be the subject of an exchange. Products returned and replaced free of charge shall remain the property of the Company.

The guarantee shall not apply to replacements or repairs resulting from normal wear and tear of the products, to deterioration or accidents due to negligence, to any failure to respect the handling, transportation, storage or sterilization conditions, or to any failure to monitor, maintain or use the products correctly. EXCEPT AS STATED IN THIS GUARANTEE, THERE ARE NO OTHER EXPRESS, IMPLIED OR STATUTORY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

This guarantee shall not extend to any consequential or any other loss or indirect damage that may result from using a Coloplast product that has a manufacturing or material defect, whether for breach of any guarantee or warranty, for breach of any other term, or for ordinary negligence, unless it is shown from the circumstances of the defect that Coloplast has been guilty of gross negligence. The provisions of the United Nations Convention on the International Sale of Goods are expressly excluded from any purchase or sales contract of Coloplast products. Except as set forth in this paragraph, no compensation shall be payable by way of damages for any reason whatsoever. The storage conditions of the products shall be the sole responsibility of the Purchaser. We recommend that the products be stored in their original packaging with this packaging intact, away from light in a dry place, and at an ambient temperature.

The fact that a product is mentioned in this catalogue does not necessarily mean that it can be supplied. Committed to continually improving product quality, Coloplast reserves the right to make changes at any time, without notice, to its products, their specifications, packaging and presence in its range.

Contact your local representative for availability of the products in your country.


The purchase or use of our products implies acceptance of the general conditions of sale set out in the price list.




# Terminology




Proximal end




Whistle tip




Distal tip




Straight cylindrical tip



Solid tip




Coudé tip



Hollow tip



Straight olive tip



Open tip



Coudé olive tip

Instrument gauges • The catheter gauges are:  
 CHARRIERE: abbreviation «Ch/Fr» = 1/3mm  
 MILLIMETER: abbreviation «mm» = 1mm  
 GAUGE: abbreviation G. Gauge 14 = 1,63mm  
 INCH: abbreviation " = 25,4mm

# Material used

Coloplast uses various types of raw materials, depending on the various diagnostic or therapeutic indications for its devices. The biocompatibility of the medical devices is in accordance to EN ISO 10993-1: 2020 Standard.

Before being used in manufacturing, raw materials are controlled in our laboratories batch by batch, using the control methods defined by our internal procedure.

All materials and components used by Coloplast are controlled and used in accordance with procedures and the requirements of standard EN ISO 13485: 2016.

The manufacturing processes depend on the material and include extrusion, injection or compression-molding, assembly and soaking.

We use two main types of material:

- Silicones
- Thermoplastics: Vortek®, Biosoft®, Neoplex, Polyamide, PVC and Polyurethane.

**SILICONE:** this is a flexible material with optimum biocompatibility. Silicone comes in different degrees of rigidity; it may be transparent, may contain radiopaque lines and may be opaque to X-rays.

After use, it degrades only to silica (sand) and water, with no chlorine or carbon derivatives. It is very stable over time.

**VORTEK:** An original formulation designed by Coloplast. The material is biocompatible; it comes in different degrees of rigidity and is opaque to X-rays. Mastery of the co-extrusion of this material has made it possible to produce dual-hardness JJs which are rigid in contact with a guide wire to facilitate glide and assembly, while having a flexible exterior for patient comfort.

**BIOSOFT:** a flexible biocompatible polymer which is particularly well tolerated by patients.

**NEOPLEX:** a specific formulation developed by Coloplast based on PVC and plasticizer.

**POLYAMIDE (PA):** particularly rigid material suitable for accessories such as funnels, tips, sheaths, etc.

**PROLYPROPYLENE (PP):** is a synthetic material made of macromolecules. It is low density and exhibits excellent mechanical properties over time, with a high level of chemical inertness and biocompatibility.

**POLYVINYLCHLORIDE (PVC):** a transparent polymer which may contain radiopaque lines and/or pigments.

It has been used for many years for the manufacture of various medical devices.

**POLYURETHANE (PU):** a material which does not contain plasticizer and which is very resistant, and with different range of elasticity available.



# Numerical index

## AA

|      |        |
|------|--------|
| AA8A | 69     |
| AA8B | 69     |
| AA8C | 69     |
| AA61 | 68, 69 |
| AA63 | 68, 69 |
| AA64 | 68, 69 |
| AA65 | 68     |
| AA66 | 68     |
| AA71 | 69     |
| AA74 | 69     |
| AA75 | 69     |
| AA81 | 69     |
| AA85 | 69     |
| AA93 | 68     |

## AB

|      |    |
|------|----|
| AB6A | 63 |
| AB6B | 63 |
| AB6C | 63 |
| AB6E | 63 |
| AB6G | 61 |
| AB6H | 61 |
| AB6J | 63 |
| AB6R | 58 |
| AB6S | 58 |
| AB60 | 61 |
| AB61 | 60 |
| AB62 | 59 |
| AB63 | 58 |
| AB64 | 58 |
| AB65 | 60 |

## AC

|      |    |
|------|----|
| AC4B | 31 |
| AC4D | 31 |
| AC4E | 31 |
| AC5C | 45 |
| AC41 | 55 |
| AC42 | 31 |
| AC43 | 31 |
| AC44 | 37 |
| AC45 | 31 |

|      |    |
|------|----|
| AC46 | 31 |
| AC47 | 55 |
| AC50 | 43 |
| AC51 | 43 |
| AC52 | 43 |
| AC53 | 43 |
| AC54 | 43 |
| AC55 | 43 |
| AC59 | 45 |
| AC67 | 55 |
| AC68 | 55 |
| AC70 | 55 |
| ACA1 | 37 |
| ACA2 | 37 |
| ACB1 | 25 |
| ACB2 | 25 |
| ACB3 | 25 |
| ACB4 | 25 |
| ACB5 | 24 |
| ACB6 | 24 |
| ACB7 | 35 |
| ACB8 | 35 |
| ACB9 | 26 |
| ACBC | 24 |
| ACBG | 24 |
| ACBH | 31 |
| ACBM | 25 |
| ACBS | 25 |
| ACN6 | 43 |
| ACP2 | 44 |
| ACP3 | 44 |
| ACP4 | 44 |
| ACP6 | 44 |
| ACR2 | 45 |
| ACXL | 20 |

## AD

|      |    |
|------|----|
| AD5D | 70 |
| ADN3 | 70 |

## AE

|      |    |
|------|----|
| AE00 | 41 |
| AE0A | 41 |

|      |        |
|------|--------|
| AE0C | 41, 49 |
| AE0D | 41     |
| AE0E | 41     |
| AE0N | 40     |
| AE0P | 40     |
| AE0T | 40     |
| AE0X | 40     |
| AE04 | 42, 49 |
| AE05 | 42, 49 |
| AE11 | 41, 49 |
| AEAD | 40     |
| AEHA | 39     |
| AEJ4 | 42, 49 |
| AEJ5 | 42, 49 |

## AG

|      |    |
|------|----|
| AG51 | 66 |
| AG53 | 66 |
| AG54 | 66 |

## AJ

|      |    |
|------|----|
| AJ4A | 22 |
| AJ4C | 22 |
| AJ4U | 35 |
| AJ4W | 33 |
| AJ4Y | 34 |
| AJ41 | 22 |
| AJ42 | 22 |
| AJ45 | 22 |
| AJ47 | 22 |
| AJ48 | 22 |
| AJ50 | 22 |
| AJ66 | 51 |
| AJ67 | 51 |
| AJ92 | 54 |

## AK

|      |    |
|------|----|
| AK20 | 56 |
| AK30 | 57 |
| AK31 | 57 |
| AK32 | 57 |
| AK33 | 57 |
| AK34 | 57 |
| AK53 | 57 |

## AL

|      |    |
|------|----|
| AL00 | 62 |
| ALXL | 20 |

## AS

|      |        |
|------|--------|
| ASP1 | 74, 75 |
| ASP2 | 75     |
| ASXL | 20     |

## AW

|      |    |
|------|----|
| AW61 | 68 |
| AW63 | 68 |
| AXXL | 20 |

## BB

|      |    |
|------|----|
| BBIO | 65 |
|------|----|

## BC

|      |    |
|------|----|
| BCAA | 28 |
| BCAB | 28 |
| BCAD | 28 |
| BCAE | 28 |
| BCAG | 28 |
| BCAJ | 28 |
| BCAM | 28 |
| BCAN | 28 |
| BCAS | 28 |
| BCCG | 36 |
| BCCS | 36 |
| BCFA | 27 |
| BCFG | 27 |
| BCFJ | 27 |
| BCFM | 27 |
| BCFP | 27 |
| BCFS | 27 |
| BCGA | 29 |
| BCGG | 29 |
| BCGS | 29 |
| BCHA | 23 |
| BCHG | 23 |
| BCHJ | 23 |
| BCHS | 23 |

## BJ

BJ30 .....55

## BN

BNBA.....32

BNBB.....32

BNBJ.....32

BNBK.....32

BNFA.....27

BNFG.....27

BNGA.....29

BNNA.....28

## CD

CD1A.....53

CD1B.....53

CD1C.....53

CD10.....53

CD11.....53

CD12.....53

## CP

CP18.....56

## CV

CV38.....70

## D

DOG1.....19

DOG2.....19

## E

E179.....83

## EN

EN28.....82

EN29.....82

## ER

ER80.....83

ER81.....83

## ES

ES28.....82

ES29.....82

## EX

EXN4.....18

EXN7.....18

EXN9.....18

EXP4.....17

EXT2.....18

EXT4.....18

EXT6.....18

## F

FSI0.....11

FSI1.....11

FSI2.....11

FSI3.....11

FSI5.....11

FSI8.....11

FSV0.....11

FSV2.....11

FSV3.....11

FSV5.....11

FSV8.....11

FSX0.....11

FSX2.....11

FSX3.....11

FSX5.....11

FSX8.....11

## GA

GA10.....77

GA12.....77

GA61.....77

GA62.....77

GA69.....77

## GB

GB60.....76

## GD

GD41.....78

GD42.....78

GD46.....78

## N

NBIO.....65

## P

PB01.....56

PB02.....56

PB12.....56

## Q

QB33.....85

QB55.....85

QB93.....85

## RA

RAC01.....50

RAC02.....50

## RB

RBA0.....50

RBC0.....50

RBD0.....46

## RC

RCJ0.....51

RCM0.....51

## RE

REA0.....51

## RJ

RJA0.....47

RJA1.....47

RJB0.....47

RJB1.....47

RJE1.....52

## RM

RMA1.....48

RMA2.....48

## T

TFLD01.....11

## U

UR31.....91

## X

XB6L.....58

XB6M.....60

XB6N.....59

XB61.....60

XB62.....59

XB63.....58

## 000

000P1B.....14

000P1C.....13

000P1S.....14

000P1U.....12

000PAB.....64

000PCP.....11

000PFL.....11

000PFS.....11

000PJU.....85

000PKB.....21

000PKE.....11

## Others

5002.....86

5013.....92

5014.....92

5177.....83

5191.....83

5192.....84

5193.....83

5196.....90

9471.....83

# Alphabetical index

## A

|                                         |    |
|-----------------------------------------|----|
| ABVac .....                             | 64 |
| Altis.....                              | 90 |
| Ancillaries for penile prostheses ..... | 85 |
| Aris .....                              | 91 |
| Assembly kits .....                     | 83 |

## B

|                                                                              |    |
|------------------------------------------------------------------------------|----|
| Balloon catheters in silicone .....                                          | 51 |
| Biliary drainage .....                                                       | 78 |
| Biosoft duo - Double loop ureteral stents .....                              | 28 |
| Biosoft duo - Multi-length hydro-coated<br>double loop ureteral stents ..... | 29 |
| BoNee .....                                                                  | 65 |
| Bougies .....                                                                | 66 |
| Bougies metric thread.....                                                   | 66 |
| BubbleView - Single-Use .....                                                | 13 |

## C

|                               |       |
|-------------------------------|-------|
| Catheters.....                | 63    |
| Catheter Valve .....          | 70    |
| Chevassu .....                | 45    |
| Coloplast TFL Drive.....      | 10-11 |
| Connection .....              | 51    |
| Connectors and spigot .....   | 56-57 |
| CoralView - Single-Use.....   | 12    |
| Couvelaire tip .....          | 60    |
| Cystodrain integral Set ..... | 53    |

## D

|                                                                             |       |
|-----------------------------------------------------------------------------|-------|
| Delinotte tip.....                                                          | 59    |
| Dilators.....                                                               | 50    |
| Double loop ureteral stents - For pyeloplasty.....                          | 34    |
| Double loop ureteral stents - For ureteral stenosis .....                   | 33    |
| Double loop ureteral stents - Without eyes<br>on the straight section ..... | 35    |
| Dormia® DualPerc .....                                                      | 16-17 |
| Dormia® Front .....                                                         | 19    |
| Dormia® N.stone® .....                                                      | 18    |
| Dormia® No-Tip .....                                                        | 18    |
| Drainage catheters .....                                                    | 43    |
| Dufour tip .....                                                            | 58    |

## E

|               |       |
|---------------|-------|
| Elefant ..... | 74-75 |
|---------------|-------|

## F

|                          |       |
|--------------------------|-------|
| Floppy tip.....          | 45    |
| Foley catheters.....     | 68-69 |
| Freudenberg Stylet ..... | 62    |

## G

|                 |    |
|-----------------|----|
| Genesis.....    | 84 |
| Guidewires..... | 49 |

## H

|                                      |    |
|--------------------------------------|----|
| Hooks & dilators.....                | 85 |
| Hydrophilic Nitinol guidewires ..... | 40 |

## I

|                                                                                                                                              |    |
|----------------------------------------------------------------------------------------------------------------------------------------------|----|
| Imajin .....                                                                                                                                 | 22 |
| Imajin hydro-coated.....                                                                                                                     | 23 |
| Interventional catheters - Floppy tip hydro-coated<br>ureteric catheter, Chevassu catheters for Retrograde<br>UreteroPyelography (RUP) ..... | 45 |
| Interventional catheters - These catheters are suitable<br>for use with a guidewire during endoscopic procedures ..                          | 44 |

## J

|                               |    |
|-------------------------------|----|
| J and Malecot catheters ..... | 51 |
|-------------------------------|----|

## K

|                                         |    |
|-----------------------------------------|----|
| Kolibri set with J catheter .....       | 47 |
| Kolibri set with Malecot catheter ..... | 48 |
| K-Box .....                             | 21 |

## L

|                  |    |
|------------------|----|
| Lunderquist..... | 42 |
|------------------|----|

## N

|               |    |
|---------------|----|
| Needles ..... | 50 |
|---------------|----|

|                                               |    |                                                |       |
|-----------------------------------------------|----|------------------------------------------------|-------|
| <b>P</b>                                      |    | <b>V</b>                                       |       |
| Philips connectable bougies . . . . .         | 66 | ViewHub S13-T . . . . .                        | 14    |
| Polyamide & Polyurethane                      |    | Virtue male sling . . . . .                    | 86    |
| double loop ureteral stent sets . . . . .     | 32 | Vortek - Double loop ureteral stents . . . . . | 24-26 |
| Prostatic Catheters tip terminology . . . . . | 71 | Vortek hydro-coated . . . . .                  | 27    |
| PTFE-coated guidewires . . . . .              | 41 | Vortek Single loop ureteral stents . . . . .   | 37    |
| Pumps . . . . .                               | 83 | Vortek Tumor stents . . . . .                  | 36    |
| PU-R, PU-S . . . . .                          | 31 |                                                |       |
| <b>R</b>                                      |    | <b>X</b>                                       |       |
| Radifocus™ . . . . .                          | 40 | XL Cylinder Set . . . . .                      | 83    |
| Restorelle . . . . .                          | 92 |                                                |       |
| ReTrace . . . . .                             | 20 |                                                |       |
| <b>S</b>                                      |    |                                                |       |
| Seldinger guidewires . . . . .                | 41 |                                                |       |
| Set with J catheter . . . . .                 | 52 |                                                |       |
| Schüller . . . . .                            | 42 |                                                |       |
| Simple drainage . . . . .                     | 77 |                                                |       |
| Soprano® . . . . .                            | 39 |                                                |       |
| Straight tip . . . . .                        | 61 |                                                |       |
| Stainless steel guidewires . . . . .          | 42 |                                                |       |
| Suction drainage . . . . .                    | 76 |                                                |       |
| Supraflow . . . . .                           | 54 |                                                |       |
| <b>T</b>                                      |    |                                                |       |
| Terumo guidewires . . . . .                   | 40 |                                                |       |
| Titan Cloverleaf Reservoir . . . . .          | 83 |                                                |       |
| Titan Touch . . . . .                         | 82 |                                                |       |
| Titan Touch Narrow Base . . . . .             | 82 |                                                |       |
| True-Lock® Plug . . . . .                     | 83 |                                                |       |
| <b>U</b>                                      |    |                                                |       |
| Urethral catheters without balloon . . . . .  | 70 |                                                |       |
| Ureteral dilators . . . . .                   | 46 |                                                |       |
| Urinary diversion . . . . .                   | 55 |                                                |       |

---

## Our mission

Making life easier for people with intimate healthcare needs

Ostomy Care | Continence Care | Advanced Wound Care | Interventional Urology | Voice and Respiratory Care

Coloplast and the Coloplast logo are trademarks of Coloplast A/S.  
[www.coloplast.com](http://www.coloplast.com) © 2025 Coloplast A/S. All rights reserved. PM-37190.

