

# Men's Health



MEN'S HEALTH

Men's health

|    |                           |    |
|----|---------------------------|----|
| 1. | Erectile Dysfunction (ED) | 82 |
| 2. | Male Urinary Incontinence | 86 |

## Titan® Touch

### Inflatable penile implant

3-pieces inflatable penile implant with hydrophilic coating and deflate button pump design.

- Touch pump: one touch of the deflate pads allows for easy deflation of the device.
- The Bioflex® material of the cylinders allows for unmatched cylinder strength and rigidity.
- Reduced risk of auto-inflation with Lock-out™ valve on the reservoir.
- Titan HydroVANTAGE™ hydrophilic coating is on all device components, including rear tip extenders and connectors. Physicians have the ability to choose the aqueous solution dip, tailored to the clinical needs of their patients.
- Zero Degree angle proximal base input tubing for ease of surgical placement.
- Soft molded tip designed for suitable anatomic positioning.
- Individually packaged • Sterile.



## Titan® Touch

### Penile prosthesis

| Length (cm)              | 14     | 16     | 18     | 20     | 22     |
|--------------------------|--------|--------|--------|--------|--------|
| Order Ref. / Infrapubic  | ES2814 | ES2816 | ES2818 | ES2820 | ES2822 |
| Order Ref. / Penoscrotal | ES2914 | ES2916 | ES2918 | ES2920 | ES2922 |

## Titan® Touch Narrow Base

### For narrow corpora cavernosa

### Penile prosthesis

| Length (cm)              | 11     | 14     | 16     | 18     |
|--------------------------|--------|--------|--------|--------|
| Order Ref. / Infrapubic  | EN2811 | EN2814 | EN2816 | EN2818 |
| Order Ref. / Penoscrotal | EN2911 | EN2914 | EN2916 | EN2918 |




## 1 • ERECTILE DYSFUNCTION

# Titan<sup>®</sup> Cloverleaf Reservoir

With Lock-out anti-reflux valve.

|                    |        |        |
|--------------------|--------|--------|
| Reservoir contents | 75 cc  | 125 cc |
| Order Ref.         | ER8075 | ER8125 |



# True-Lock<sup>®</sup> Plug

The True-Lock polysulfone plug and sleeve are utilized to close unused tubing on a portion of the TitanR IPP, where deemed appropriate.

|            |       |
|------------|-------|
| Order Ref. | 9471S |
|------------|-------|



# Assembly kits

Standard assembly kit.

|                       |        |
|-----------------------|--------|
| Order Ref. / Standard | 519150 |
|-----------------------|--------|

### Contains

- Connector tool (1) to compress connector sleeves
- True-Lock<sup>®</sup> Connectors (connector sleeve (8), connector handle (2), connector cage (4))
- Keith needle (2)
- Pusher protector (1), for corporotomy closure and proximal insertion
- Barbed syringe tip (3), to fill the device
- Flexible tip needle (2), pink syringe tip to remove air from silicone tubing
- Silicone tubing (6)



Rear tip extenders kit.

|                          |        |
|--------------------------|--------|
| Order Ref. / Standard    | 519310 |
| Order Ref. / Narrow Base | 519300 |

To be used with the standard assembly kit.



# XL Cylinder Set

|            |        |        |        |
|------------|--------|--------|--------|
| Cylinders  | 24cm   | 26cm   | 28cm   |
| Order Ref. | E17924 | E17926 | E17928 |

Pump to be connected.



# Pumps

|            |            |                  |
|------------|------------|------------------|
| Pumps      | Titan Pump | Titan Touch Pump |
| Order Ref. | 517770     | 517750           |

Titan Touch Pump



# Genesis®

## Malleable penile implant

Malleable penile implant with hydrophilic coating.

- **Kit includes:**
  - One pair of cylinders.
  - 3 pairs of caps (standard, +0.5cm and +1cm).
- Hydrophilic coating.
- **Individually packaged (pairs)** • Sterile.



| Diameter | Length cm |           | Order Ref. |
|----------|-----------|-----------|------------|
|          | mm        | Min. Max. |            |
| 9.5      | 14        | 23        | 519250     |
| 11       | 16        | 25        | 519260     |
| 13       | 18        | 27        | 519270     |



## 1 • ERECTILE DYSFUNCTION

# Ancillaries for penile prostheses

## Hooks & dilators\_

**Furlow**  
Non sterile.

Order Ref. QB5555

- Furlow insertion tool



**Brooks**  
Non sterile.

- Brooks dilator



|                   | Ø   |     |      |      |      |      |      |
|-------------------|-----|-----|------|------|------|------|------|
| Order Ref. QB9350 | 8mm | 9mm | 10mm | 11mm | 12mm | 13mm | 14mm |

**Rosselló Dilators**  
Non sterile.


Order Ref. QB5558



**Retractor System without Frame**  
Sterile - Single use.

Order Ref. QB3321

- 6 - Retractor Stays, 12mm Blunt
- 1 - Retractor Stay, 5mm Sharp
- 1 - Penile Elevation Strap
- 2 - Disposable deavers
- 2 - Retractor Stays, 20mm 3-finger Rake



**Galaxy II Surgical Retractor**  
Sterile - Single use.

- Male kit

Order Ref. 000PJUNEGX-2020

- 6 - 12 mm blunt stays
- 1 - 5mm sharp stay
- 2 - 3 finger claw
- 2 - 7mm double hook sharp
- 1 - penile hammock
- 1 - snowman retractor frame
- 1 - catheter clip
- 1 - light/lux clip

## Virtue<sup>®</sup> male sling

### 4-arm male sling for male stress urinary incontinence\_

- Knitted monofilament polypropylene.
- Bidirectional compression and elevation of the bulbous urethra with 2 transobturator arms and 2 prepubic arms.
- **Individually packaged** • Sterile.



| Order Ref. | Description                              | Sling dimensions |
|------------|------------------------------------------|------------------|
| 50020      | Virtue: 1 sling, 1 single use introducer | 46 x 31cm        |

---

## Our mission

Making life easier for people with intimate healthcare needs

CARJ0CA 2024229

[Ostomy Care](#) | [Continence Care](#) | [Advanced Wound Care](#) | [Interventional Urology](#) | [Voice and Respiratory Care](#)

Coloplast and the Coloplast logo are trademarks of Coloplast A/S.  
[www.coloplast.com](http://www.coloplast.com) © 2025 Coloplast A/S. All rights reserved. PM-37190.

