

Surgery



SURGERY

Surgery

1.	Laparoscopic surgery	75
2.	Surgical drainage	76

Elefant®

Disposable suction-irrigation device for laparoscopy_



- Single use cannula:
- Ergonomic and lightweight.
- Yankauer-type tip (6 distal holes).
- “E” fiberglass (non-conductive).

Elefant, without tubing • Box of 10			
Cannula diameter	5mm	10mm	5mm
Cannula length	35cm	35cm	45cm
Order Ref.	ASP145	ASP150	ASP180*

* For robotic surgery.

Elefant, with double irrigation line tubing • Box of 5		
Cannula diameter	5mm	5mm
Cannula length	35cm	45cm
Tube diameter	ext. 9.5mm int 6.4mm	ext. 9.5mm int 6.4mm
Order Ref.	ASP190	ASP195*

* For robotic surgery.



Please contact your Coloplast Territory Manager regarding this product’s availability in your market.

1 • LAPAROSCOPIC SURGERY

Non-traumatic and
“E” fiberglass
cannula

Yankauer tip



Elefant cannula made
of non-conductive fiberglass,
unlike conventional irrigation devices
made of stainless steel

Elefant, with PVC tubing • Box of 5			
Cannula diameter	5mm	10mm	5mm
Cannula length	35cm	35cm	45cm
Tube diameter	ext. 9.5mm int 6.4mm	ext. 9.5mm int 6.4mm	ext. 9.5mm int 6.4mm
Tube length	3m	3m	3m
Order Ref.	ASP165	ASP170	ASP185*

* For robotic surgery.



Elefant, with marked cannula Box of 5	
5mm	5mm
35cm	45cm
ext. 9.5mm int 6.4mm	ext. 9.5mm int 6.4mm
3m	3m
ASP265	ASP285*

* For robotic surgery.

Suction drainage

- Ch/Fr 9-24: **box of 10** • Ch/Fr 27-36: **box of 5** • Sterile.

Multiperforated tubular drains • Material: Silicone • Length: 50cm										
Ch/Fr	9	12	15	18	21	24	27	30	33	36
Order Ref.	GB6009	GB6012	GB6015	GB6018	GB6021	GB6024	GB6027	GB6030	GB6033	GB6036

Radiopaque line.

Ch/Fr 9-24: box of 10 • Ch/Fr 27-36: box of 5.



2 • SURGICAL DRAINAGE


Simple drainage

- GA1035, GA1245, GA62xx, GA69xx Ch/Fr 8 to 25: **box of 10** • GA6106, GA69xx Ch/Fr 30 to 35: **box of 5** • Sterile.

Tubular drains • Material: Silicone • Length: 30cm										
Ch/Fr	8	10	12	14	16	18	20	22	25	30
Order Ref.	GA6208	GA6210	GA6212	GA6214	GA6216	GA6218	GA6220	GA6222	GA6225	GA6230

Ch/Fr	35	40
Order Ref.	GA6235	GA6240

Radiopaque line.





Tubular drains • Material: Silicone • Length: 50cm										
Ch/Fr	8	10	12	14	16	18	20	25	30	35
Order Ref.	GA6908	GA6910	GA6912	GA6914	GA6916	GA6918	GA6920	GA6925	GA6930	GA6935



Multitubular drain • Material: Radiopaque Silicone • Length: 40cm	
Width (mm)	35
Order Ref.	GA1035

Multitubular drain • Material: Silicone • Length: 40cm	
Width (mm)	45
Order Ref.	GA1245




Delbet drains • Material: Silicone • Length: 36cm	
Width (cm)	6
Order Ref.	GA6106

Box of 5. Radiopaque

Biliary drainage

- Box of 5 • Sterile.




Transcystic drains • Material: PVC • Length: 41.5cm			
Bulb Ø Ch/Fr	7	8	10
Tube Ø Ch/Fr	5	6	8
Order Ref.	GD4105	GD4106	GD4108

Pedinielli drains with malleable stylet, whistle tip, 2 eyes.
Radiopaque line.



Transcystic drains • Material: PVC • Length: 154cm				
Bulb Ø Ch/Fr	7	8	10	12
Tube Ø Ch/Fr	5	6	8	10
Order Ref.	GD4205	GD4206	GD4208	GD4210

Escat drains, whistle tip, 2 eyes, with 3-way stopcock.
Radiopaque line.



Cœliodrain® • Material: Radiopaque Neoplex • Length: approximately 64cm	
Ch/Fr	5
Order Ref.	GD4605

For cholangiography with malleable stylet and whistle tip.

Our mission

Making life easier for people with intimate healthcare needs

Ostomy Care | Continence Care | Advanced Wound Care | Interventional Urology | Voice and Respiratory Care

Coloplast and the Coloplast logo are trademarks of Coloplast A/S.
www.coloplast.com © 2025 Coloplast A/S. All rights reserved. PM-37190.

