

Urology - Kidney & Bladder Health



Urology

1 • Endourology	10
1.1. Thulium Fiber Laser	10
1.2. Visualization	12
1.3. Stone extractors	16
1.4. Access sheaths	20
1.5. K-Box	21
1.6. Double loop ureteral stents	22
1.7. Single loop ureteral stents	37
1.8. Guidewires	39
1.9. Ureteral catheters	43
1.10. Ureteral dilators	46
2 • Percutaneous access	47
2.1. Percutaneous nephrostomy set	47
2.2. Percutaneous nephrostomy accessories	49
2.3. Direct puncture	52
2.4. Suprapubic drainage	53
2.5. Urinary Diversion	55
3 • Connectors and spigot	56
3.1. Connectors and Spigot	56
4 • Lower urinary tract	58
4.1. Prostatic catheters	58
4.2. Bladder evacuator	64
4.3. Bladder injection	65
4.4. Urethral dilation	66
4.5. Foley catheters	68
4.6. Catheters without balloon	70

Coloplast TFL Drive

Power made Safer*

Laser System

Coloplast TFL Drive is a 60W laser system intended for multiple indications: Lithotripsy, BPH and Soft Tissue. It's the only TFL laser to have **safety monitoring thanks** to preset safety ranges and defined power limits. Intuitive Interface are designed by surgeon for surgeons and aiming to adapt power to the size of the fiber for specific clinical condition.

The Coloplast TFL Drive interface allows you to adapt energy, frequency and power values to the size of fiber used.

The **TFL technology** is the latest cutting-edge laser technology, bringing highly effective intra-operative performance. Finer and Faster dusting capabilities combined with less stone retropulsion of TFL allow decrease treatment time.

Coloplast TFL Drive in Key Points

All-in-one Solution

- Stone
- BPH
- Soft Tissue

Intuitive Interface

- Defined aiming at optimizing the effectiveness

Ergonomic design

- Easy transportability
- Limited footprint

Efficiency

- Finer and faster dusting capabilities
- Less retropulsion
- High ablation speed

More Focus

- Quiet System
- Noise level <70dBA

Profesional education

- Support you to learn this technology with Europe's leading experts

Technical Data

Laser information

- Type TFL (Thulium Fiber Laser)
- Pulse energy 0.02 – 6 J
- Repetition rate 1 – 2500Hz
- Max. Power 60W (CW and pulse mode)
- Pulse Duration 50µs – CW

Electrical Data

- Voltage 100-240 vac
- Frequency 50/60Hz
- Power Consumption 1000 VA

- Wavelength 1940nm ± 20nm
- Noise Level max. <70dBA
- Aiming Beam Green 532nm (adjustable) <5mW, class 3R

Dimensions

- Size 470 (W) x 810 (D) x 1156 (H) mm
470 (W) x 810 (D) x 940 (H) mm (display closed)
- Weight 100kg
- Cooling system Air cooling system



* User interface designed with warning limits for safer use.

1 • ENDOUROLOGY


1.1. Thulium Fiber Laser

Ordering Information

Order Ref.	Description
TFLD01	Coloplast TFL Drive Laser including: Coloplast TFL Drive Laser System - Pedal Double footswitch - Optical blast shield - Safety goggles (x2) - Power cord - User manual - Accessory case containing: Fiber stripper for optical fibers 300 to 1000µm, Fiber stripper for optical fibers 100 to 400µm, interlock connector, optical blast shield, and key for use with Coloplast TFL Drive.

Fibers


TFL is suitable with use with the smallest fiber 150 µm, allowing to free up space inside the ureterscope working channel to increase irrigation and thus visibility.



Fiber size	Order Ref.			Ordre Unit
	Single use	Reusable 5x	Reusable 10x	
150µm	FSI150			5 per box
200µm	FSI200	FSV200	FSX200	5 per box
272µm	FSI272	FSV272	FSX272	5 per box
272µm	FBI272 Ball tip			5 per box
365µm	FSI365	FSV365	FSX365	5 per box
550µm	FSI550	FSV550	FSX550	5 per box
800µm	FSI800	FSV800	FSX800	5 per box
1000µm	FSI000	FSV000	FSX000	5 per box


Spare Parts

Order Ref.	Description
000PCPGL01	Laser protection goggles
000PCPCL01	Ceramic fiber cleaver with slipcase
000PCPIC01	Interlock connector
000PFS0104	Fiber Stripper 100 to 400µm
000PFS0310	Fiber Stripper 300 to 1000µm
000PCPDFS1	Pedal Double footswitch
000PCPOBS1	Optical Blast Shield
000PCPCASE	Accessories case
000PKEY001	Key
000PFLGCASE	TFL Flight case



CoralView - Single-Use Digital Flexible Ureteroscope*

- The CoralView is used by physicians to access, visualize and perform procedures in the urinary tract.
- Each unit is supplied with a Leak resistant and rotatable T-valve.
- 120° field of view, 3-50mm depth of view.
- Accurate color reproduction.
- Autohold technology self-locking lever.
- **Box of 5 • Sterile.**




Order Ref.	Description	Distal end outer diameter	Nominal insertion tube outer diameter	Maximum insertion portion width	Deflection mode
000P1U10B	CoralView Ureteroscope – Standard	7.2 Ch/Fr	8.25 Ch/Fr	2.9 mm	Standard
000P1U10BR	CoralView Ureteroscope – Reverse	7.2 Ch/Fr	8.25 Ch/Fr	2.9 mm	Reverse
000P1U20B	CoralView Ureteroscope – Standard	6.7 Ch/Fr	7.5 Ch/Fr	2.7 mm	Standard
000P1U20BR	CoralView Ureteroscope – Reverse	6.7 Ch/Fr	7.5 Ch/Fr	2.7 mm	Reverse

* For France and Spain only.

Deflection

Reverse deflection Moving the lever downward will make the tip bend down




1 • ENDOUROLOGY

1.2 Visualization

BubbleView - Single-Use

Digital Flexible Cystoscope*


- The BubbleView is used by physicians to access, visualize and perform procedure in the bladder.
- Depth of view: 3mm - 100mm.
- Field of view: 140°.
- Digital scope, HD 920k pixels.
- **Box of 5 • Sterile.**



Order Ref.	Description	Nominal Insertion tube outer diameter	Maximum insertion portion width	Working channel width	Deflection mode
000P1C50C	BubbleView Cystoscope – Standard	14.4 Ch/Fr	17.4 Ch/Fr	6.6 Ch/Fr	Standard
000P1C50CR	BubbleView Cystoscope – Reverse	14.4 Ch/Fr	17.4 Ch/Fr	6.6 Ch/Fr	Reverse

* For France and Spain only.

Standard deflection *Moving the lever downward will make the tip bend up*




ViewHub S13-T

Video Processor*

Compatible with both CoralView ureteroscope and BubbleView cystoscope, adapting to various scenarios in outpatient and operating rooms.

Features

- 1920 x 1080 pixels.
- Snapshot.
- Image zooming.
- Video recording.
- Image freezing.
- Brightness control.
- White balance.



Order Ref.	Description
000P1S13T	ViewHub S13 - Video processor

* For France and Spain only.

Accessories



Order Ref.	Description
000P1BRACKET	Bracket to put screen on IV pole
000P1TROLLEY	Trolley without holder
000P1TROLLEYPLUS	Trolley with holder (removable part)

Dormia® DualPerc

A unique concept_

Dormia® DualPerc offers dual functionality thanks to its quick and easy switch between flexible and rigid configurations.

When used with both configurations, the device is designed to remove stones through a unique percutaneous access, reducing multiple renal punctures and/or the need of additional devices.



Rigid configuration

Working length: 35 cm




Nitinol basket

- Super elasticity and maximum shape memory
- 4 wires cage for capture and retrieval of stones of different sizes and volumes
- Tipless design for direct contact against tissues with minimal trauma
- Specific design for adapted opening in small cavities

Stainless-steel cannula

- Rounded distal tip for smoother contact with the tissue and instruments, designed to minimize trauma
- Radiopaque tube for accurate visualization under x-ray



Flexible mandrel

- Diameter of the sheath and mandrel permits full deflection of the flexible cystoscope during stone extraction
- Working length designed for optimal compatibility with flexible cystoscope

Flexible configuration


Full length: 64 cm

3 Ch/Fr flexible sheath

1 • ENDOUROLOGY

1.3. Stone extractors

Order Ref.	Description	Rigid Shaft Length	Rigid Shaft Diameter	Flexible Shaft Length	Flexible Shaft Diameter	Basket Diameter	Basket Length	Ordre Unit
EXP448	Dormia® DualPerc Stone Extractor 4.8 Ch/Fr	35cm	4.8 Ch/Fr	64cm	3 Ch/Fr	12.5mm	23mm	5 per box
EXP470	Dormia® DualPerc Stone Extractor 7 Ch/Fr	35cm	7 Ch/Fr	64cm	3 Ch/Fr	12.5mm	23mm	5 per box




Dormia® No-Tip Nitinol • 4 twisted wires

- Individually packaged • Sterile • Detachable handle.

Order Ref.	Ch/Fr	Length	Basket size
EXN434	3	90cm	14mm

Designed for semi-rigid and rigid ureteroscopy.





Nitinol • 4 round wires.

- Individually packaged • Sterile • Detachable handle.

Order Ref.	Ch/Fr	Length	Basket size
EXN734	2.2	120cm	11mm
EXN934	1.5	120cm	9mm


Designed for flexible, semi-rigid and rigid ureteroscopy.
Supplied with a waterproof connector (Tuohy-Borst).



Dormia® N.stone® Nitinol • 4 helical wires

- Individually packaged • Sterile • Detachable handle.

Order Ref.	Ch/Fr	Length	Basket size
EXT224	4	90cm	15mm
EXT424	3	90cm	15mm
EXT624	2.5	90cm	12.5mm



1 • ENDOUROLOGY

1.3. Stone extractors

Dormia® Front **Nitinol • Frontal design**

• Individually packaged • Sterile.

Order Ref.	Ch/Fr	Length	Basket size
DOG158	1,5	120cm	8mm
DOG151	1,5	120cm	11mm
DOG208	2	120cm	8mm
DOG201	2	120cm	11mm


2 basket sizes opening, 8mm and 11mm to treat multiple stone sizes and volumes

Reinforced tip with shrink tubing made of polyester to assure durability.

Large opening to the front as in the grasper



Secure wrapping like a basket



ReTrace®

Ureteral access sheaths


- Radiopaque ring.
- Hydro-coated reinforced sheath.
- Compatible with a 0.038" stiff guidewire.
- Individually packaged • Sterile.



Single wire concept: the working guidewire turns into a safety wire.

Flexible tip concept: designed to respect the anatomy and prevent ureteral injuries.

Order Ref.	Sheath inner diameter (Fr)	Sheath outer diameter (Fr)	Length (cm)
ASXL10	10	12	28
ACXL10	10	12	35
AXXL10	10	12	45
ALXL10	10	12	55
ASXL12	12	14	28
ACXL12	12	14	35
AXXL12	12	14	45
ALXL12	12	14	55



- Connector: ergonomic shape.

1 • ENDOUROLOGY

1.5. K-Box

K-Box®

The training model

Unique 3D model for flexible ureteroscopy simulation and training.


It allows to learn:

- To manipulate the scope.
- To manipulate the scope with single use devices as baskets to catch and relocate stones.
- To perform step by step techniques.
- To use laser in the different settings.



Order Ref.

000PKBOX01





Silicone double loop ureteral stents

The remarkable biocompatibility and flexibility of silicone makes it the material of choice for long-term implantation: improved patient comfort* and less encrustation material compared to Polyurethane.

Duration of implantation: up to 12 months.



Kits with
simple pusher

- **Kit composition:** each stent is supplied in a kit with a pusher (simple or steerable) and with or without a guidewire depending on the model.
Ø Guide: 0.035".
Steerable pusher: with or without radiopaque ring.

- **Box of 1** (with a guidewire) • Sterile.
- **Box of 5** (without a guidewire) • Sterile.



Kits with
steerable pusher

Ch/Fr	Length	O/O	O/O					O/C		C/C
		Fixed core guidewire	Fixed core guidewire	Movable core guidewire	Terumo guidewire	Poseidon® guidewire	Without guidewire	Fixed core guidewire	Movable core guidewire	Fixed core guidewire
6	16cm	AJ4A61	AJ4261			AJ5061		AJ4861	AJ4761	AJ4161
	20cm		AJ4262	AJ4562	AJ42D2	AJ5062	AJ4C62			
	22cm	AJ4A67	AJ4267		AJ42D7		AJ4C67			
	24cm	AJ4A63	AJ4263	AJ4563	AJ42D3	AJ5063	AJ4C63	AJ4863	AJ4763	AJ4163
	26cm	AJ4A64	AJ4264	AJ4564	AJ42D4	AJ5064	AJ4C64	AJ4864	AJ4764	AJ4164
	28cm	AJ4A65	AJ4265	AJ4565	AJ42D5	AJ5065	AJ4C65	AJ4865	AJ4765	AJ4165
	30cm		AJ4266		AJ42D6	AJ5066	AJ4C66			
7	16cm	AJ4A71	AJ4271			AJ5071		AJ4871		AJ4171
	20cm	AJ4A72	AJ4272		AJ42E2	AJ5072	AJ4C72			
	22cm	AJ4A77	AJ4277		AJ42E7		AJ4C77			
	24cm	AJ4A73	AJ4273	AJ4573	AJ42E3	AJ5073	AJ4C73	AJ4873	AJ4773	AJ4173
	26cm	AJ4A74	AJ4274	AJ4574	AJ42E4	AJ5074	AJ4C74	AJ4874	AJ4774	AJ4174
	28cm	AJ4A75	AJ4275	AJ4575	AJ42E5	AJ5075	AJ4C75	AJ4875	AJ4775	AJ4175
	30cm	AJ4A76	AJ4276	AJ4576	AJ42E6	AJ5076	AJ4C76	AJ4876	AJ4776	AJ4176
8	16cm		AJ4281			AJ5081				
	20cm		AJ4282			AJ5082				
	22cm	AJ4A87	AJ4287		AJ42F7		AJ4C87			
	24cm	AJ4A83	AJ4283		AJ42F3	AJ5083	AJ4C83			
	26cm	AJ4A84	AJ4284	AJ4584	AJ42F4	AJ5084	AJ4C84	AJ4884		
	28cm	AJ4A85	AJ4285	AJ4585	AJ42F5	AJ5085	AJ4C85	AJ4885		
	30cm		AJ4286	AJ4586	AJ42F6	AJ5086	AJ4C86			

* Compared to Percuflex™ plus (Boston Scientific) stent.

1 • ENDOUROLOGY

1.6. Double loop ureteral stents

Imajin[®] hydro-coated

Silicone double loop ureteral stents

The remarkable biocompatibility of silicone makes it the material of choice for long-term implantation: improved patient comfort* and less encrustation material compared to Polyurethane.

Hydro-coated stents to facilitate advancement.

Duration of implantation: up to 12 months.

- **Kit composition:** each stent is supplied in a kit with a steerable pusher and with or without a guidewire depending on the model.

Ø Guide: 0.035".

Steerable pusher: with or without radiopaque ring.

- **Box of 1** (with a guidewire) • Sterile.
- **Box of 5** (without a guidewire) • Sterile.



		o/o			
Ch/Fr	Length	Fixed core guidewire	Terumo guidewire	Poseidon [®] guidewire	Without guidewire
6	16cm	BCHA61	BCHG61	BCHS61	BCHJ61
	20cm	BCHA62	BCHG62	BCHS62	BCHJ62
	22cm	BCHA67	BCHG67	BCHS67	BCHJ67
	24cm	BCHA63	BCHG63	BCHS63	BCHJ63
	26cm	BCHA64	BCHG64	BCHS64	BCHJ64
	28cm	BCHA65	BCHG65	BCHS65	BCHJ65
	30cm	BCHA66	BCHG66	BCHS66	BCHJ66
7	16cm	BCHA71	BCHG71	BCHS71	BCHJ71
	20cm	BCHA72	BCHG72	BCHS72	BCHJ72
	22cm	BCHA77	BCHG77	BCHS77	BCHJ77
	24cm	BCHA73	BCHG73	BCHS73	BCHJ73
	26cm	BCHA74	BCHG74	BCHS74	BCHJ74
	28cm	BCHA75	BCHG75	BCHS75	BCHJ75
	30cm	BCHA76	BCHG76	BCHS76	BCHJ76
8	16cm	BCHA81	BCHG81	BCHS81	BCHJ81
	20cm	BCHA82	BCHG82	BCHS82	BCHJ82
	24cm	BCHA83	BCHG83	BCHS83	BCHJ83
	26cm	BCHA84	BCHG84	BCHS84	BCHJ84
	28cm	BCHA85	BCHG85	BCHS85	BCHJ85
	30cm	BCHA86	BCHG86	BCHS86	BCHJ86

* Compared to Percuflex[™] plus (Boston Scientific) stent.

Duration
of implantation:
up to 12 months.

Vortek®

Double loop ureteral stents

The double layer structure makes it easy to insert over the guidewire, the stiffness of each size has been adapted in order to keep the mechanical properties.




- **Kit composition:** each stent is supplied in a kit with a pusher (simple or steerable), with or without a guidewire and with or without a Flush catheter depending on the model.

Ø Guide: 0.035".

Steerable pusher: with or without radiopaque ring.

- **Box of 1**, except for Order Ref. ACB6xx (box of 5) • Sterile.



Kits with
simple pusher

Ch/Fr	Length	O/O			O/C
		Fixed core guidewire	Without guidewire	Terumo guidewire + Flush ACP6	Fixed core guidewire
Stent 4.8 Ch/Fr + Pusher 4.8 Ch/Fr 75cm	12cm	ACB550			ACBC50
	16cm	ACB551			ACBC51
	20cm	ACB552			ACBC52
	22cm	ACB557			ACBC57
	24cm	ACB553			ACBC53
	26cm	ACB554			ACBC54
	28cm	ACB555			ACBC55
	30cm	ACB556			ACBC56
Stent 6 Ch/Fr + Pusher 6 Ch/Fr 40cm	24cm	ACB563	ACB663		ACBC63
	26cm	ACB564	ACB664		ACBC64
	28cm	ACB565	ACB665		ACBC65
	30cm	ACB566	ACB666		ACBC66
Stent 7 Ch/Fr + Pusher 7 Ch/Fr 40cm	24cm	ACB573			ACBC73
	26cm	ACB574	ACB674	ACBG74	ACBC74
	28cm	ACB575	ACB675	ACBG75	ACBC75
	30cm	ACB576	ACB676		ACBC76
Stent 8 Ch/Fr + Pusher 7 Ch/Fr 40 cm	26cm	ACB584			
	28cm	ACB585			
	30cm	ACB586			

1 • ENDOUROLOGY

1.6. Double loop ureteral stents

		 Kits with simple pusher				 Kits with steerable pusher			
		O/O							
Ch/Fr	Length	Fixed core guidewire	Fixed core guidewire + Flush ACP2	Movable core guidewire	Movable core guidewire + Flush ACP6	Terumo guidewire	Poseidon URO guidewire	Without guidewire	Without guidewire + Flush ACP2
Stent 4.8 Ch/Fr + Pusher 4.8 Ch/Fr 75cm	12cm					ACB1B0	ACBS40		
	16cm					ACB1B1	ACBS41		
	20cm					ACB1B2	ACBS42		
	22cm					ACB1B7	ACBS47		
	24cm		ACB253			ACB1B3	ACBS43		
	26cm		ACB254			ACB1B4	ACBS44		ACB454
	28cm		ACB255			ACB1B5	ACBS45		ACB455
	30cm		ACB256			ACB1B6	ACBS46		
Stent 4.8 Ch/Fr + Pusher 6 Ch/Fr 40cm	12cm	ACB1C0				ACB1G0	ACBS50		
	16cm	ACB1C1				ACB1G1	ACBS51		
	20cm	ACB1C2				ACB1G2	ACBS52		
	22cm	ACB1C7				ACB1G7	ACBS57		
	24cm	ACB1C3	ACB2C3			ACB1G3	ACBS53		ACB4C3
	26cm	ACB1C4	ACB2C4			ACB1G4	ACBS54		ACB4C4
	28cm	ACB1C5	ACB2C5			ACB1G5	ACBS55		ACB4C5
	30cm	ACB1C6				ACB1G6	ACBS56		ACB4C6
Stent 6 Ch/Fr + Pusher 6 Ch/Fr 40cm	20cm	ACB162				ACB1D2	ACBS62		
	22cm	ACB167				ACB1D7	ACBS67	ACBM67	
	24cm	ACB163	ACB263	ACB363		ACB1D3	ACBS63	ACBM63	ACB463
	26cm	ACB164	ACB264	ACB364		ACB1D4	ACBS64	ACBM64	ACB464
	28cm	ACB165	ACB265	ACB365		ACB1D5	ACBS65	ACBM65	ACB465
	30cm	ACB166	ACB266	ACB366		ACB1D6	ACBS66	ACBM66	
Stent 7 Ch/Fr + Pusher 7 Ch/Fr 40cm	20cm	ACB172				ACB1E2	ACBS72		
	22cm	ACB177				ACB1E7	ACBS77		
	24cm	ACB173	ACB273	ACB373		ACB1E3	ACBS73	ACBM73	ACB473
	26cm	ACB174	ACB274	ACB374	ACB3C4	ACB1E4	ACBS74	ACBM74	ACB474
	28cm	ACB175	ACB275	ACB375	ACB3C5	ACB1E5	ACBS75	ACBM75	ACB475
	30cm	ACB176	ACB276	ACB376		ACB1E6	ACBS76	ACBM76	ACB476
Stent 8 Ch/Fr + Pusher 7 Ch/Fr 40cm	24cm	ACB183	ACB283			ACB1F3	ACBS83		
	26cm	ACB184	ACB284			ACB1F4	ACBS84		
	28cm	ACB185	ACB285			ACB1F5	ACBS85		
	30cm	ACB186				ACB1F6	ACBS86		

Duration
of implantation:
up to 12 months.

Vortek®

Double loop ureteral stents

Soft outer layer for the patient comfort.

The double layer structure makes it easy to insert over the guidewire, the stiffness of each size has been adapted in order to keep the mechanical properties.



Ch/Fr	Length	O/C Fixed core guidewire
Stent 4.8 Ch/Fr + Pusher 4.8 Ch/Fr 75cm	12cm	ACB950
	16cm	ACB951
	20cm	ACB952
	22cm	ACB957
	24cm	ACB953
	26cm	ACB954
	28cm	ACB955
	30cm	ACB956
Stent 4.8 Ch/Fr + Pusher 6 Ch/Fr 40cm	24cm	ACB9C3
	26cm	ACB9C4
	28cm	ACB9C5
Stent 6 Ch/Fr + Pusher 6 Ch/Fr 40cm	24cm	ACB963
	26cm	ACB964
	28cm	ACB965
	30cm	ACB966
Stent 7 Ch/Fr + Pusher 7 Ch/Fr 40cm	24cm	ACB973
	26cm	ACB974
	28cm	ACB975
	30cm	ACB976
Stent 8 Ch/Fr + Pusher 7 Ch/Fr 40cm	26cm	ACB984
	28cm	ACB985

Kits with
simple pusher

Kits with
steerable pusher

1 • ENDOUROLOGY

1.6. Double loop ureteral stents

Vortek® hydro-coated

Double loop ureteral stents

Soft outer layer for the patient comfort. The double layer structure makes it easy to insert over the guidewire.

The hydrophilic coating facilitates the advancement of the stent.


- **Kit composition:** each stent is supplied in a kit with a pusher (simple or steerable), with or without a guidewire and with or without a Flush catheter depending on the model.

Ø Guide: 0.035".

Steerable pusher: with or without radiopaque ring.

- **Box of 1** (with a guidewire) • Sterile.
- **Box of 5** (without a guidewire) • Sterile.

Duration
of implantation:
up to 12 months.



Kits with
simple pusher

Kits with
steerable pusher

		O/O			O/O					
Ch/Fr	Length	Fixed core guidewire	Terumo guidewire	Terumo guidewire + Flush ACP6	Fixed core guidewire	Fixed core guidewire + Flush ACP2	Terumo guidewire	Terumo guidewire + Flush ACP2	Poseidon URO guidewire	Without guidewire
Stent 4.8 Ch/Fr + Pusher 4.8 Ch/Fr 75cm	12cm				BCFA50		BCFG50			BCFJ50
	16cm				BCFA51		BCFG51			BCFJ51
	20cm				BCFA52		BCFG52			BCFJ52
	22cm				BCFA57		BCFG57		BCFS57	BCFJ57
	24cm				BCFA53		BCFG53		BCFS53	BCFJ53
	26cm				BCFA54		BCFG54		BCFS54	BCFJ54
	28cm				BCFA55		BCFG55		BCFS55	BCFJ55
	30cm				BCFA55		BCFG56			BCFJ56
Stent 4.8 Ch/Fr + Pusher 6 Ch/Fr 40cm	12cm				BCFA40		BCFG40			BCFJ40
	16cm				BCFA41		BCFG41			BCFJ41
	20cm				BCFA42		BCFG42			BCFJ42
	22cm				BCFA47		BCFG47		BCFS47	BCFJ47
	24cm				BCFA43		BCFG43		BCFS43	BCFJ43
	26cm				BCFA44		BCFG44		BCFS44	BCFJ44
	28cm				BCFA45		BCFG45		BCFS45	BCFJ45
	30cm				BCFA46		BCFG46			BCFJ46
Stent 6 Ch/Fr + Pusher 6 Ch/Fr 40cm	20cm				BCFA62		BCFG62			BCFJ62
	22cm				BCFA67	BCFM67	BCFG67		BCFS67	BCFJ67
	24cm				BCFA63	BCFM63	BCFG63		BCFS63	BCFJ63
	26cm	BNFA64	BNFG64	BNFQ64	BCFA64	BCFM64	BCFG64	BCFP64	BCFS64	BCFJ64
	28cm	BNFA65	BNFG65	BNFQ65	BCFA65	BCFM65	BCFG65	BCFP65	BCFS65	BCFJ65
	30cm		BNFG66	BNFQ66	BCFA66		BCFG66		BCFS66	BCFJ66
Stent 7 Ch/Fr + Pusher 7 Ch/Fr 40cm	20cm				BCFA72		BCFG72			BCFJ72
	22cm				BCFA77	BCFM77	BCFG77		BCFS77	BCFJ77
	24cm				BCFA73	BCFM73	BCFG73		BCFS73	BCFJ73
	26cm	BNFA74	BNFG74	BNFQ74	BCFA74	BCFM74	BCFG74	BCFP74	BCFS74	BCFJ74
	28cm	BNFA75	BNFG75	BNFQ75	BCFA75	BCFM75	BCFG75	BCFP75	BCFS75	BCFJ75
	30cm		BNFG76	BNFQ76	BCFA76		BCFG76		BCFS76	BCFJ76

Duration
of implantation:
up to 12 months.

Biosoft® duo

Double loop ureteral stents

Produced by coextrusion of a flexible material and a rigid material, Biosoft duo was designed for easy insertion, while retaining great flexibility for the patient comfort.



- **Kit composition:** each stent is supplied in a kit with a pusher (simple or steerable), with or without a guidewire and with or without a Flush catheter depending on the model.
Ø Guide: 0.035".
Steerable pusher: with or without radiopaque ring.
- **With a guidewire:** box of 1 • **Without a guidewire:** Box of 5.
- Sterile.

Kits with
simple pusher

Kits with
steerable pusher

		O/O	O/O								O/C	
Ch/Fr	Length (cm)	Fixed core guidewire	Fixed core guidewire	Fixed core guidewire + Flush ACP2xx	Movable core guidewire	Terumo guidewire	Terumo guidewire + Flush ACP2xx	Terumo guidewire + Flush ACP6xx	Poseidon URO guidewire*	Without guidewire	Fixed core guidewire	Movable core guidewire
Stent 4.8 Ch/Fr + Pusher 4.8 Ch/Fr 75cm	12cm	BNAA50		BCAM50	BCAD50	BCAGB0	BCAG50	BCAN50	BCAS50	BCAJ50	BCAB50	BCAE50
	16cm	BNAA51		BCAM51	BCAD51	BCAGB1	BCAG51	BCAN51	BCAS51	BCAJ51	BCAB51	BCAE51
	20cm	BNAA52		BCAM52	BCAD52	BCAGB2	BCAG52	BCAN52	BCAS52	BCAJ52	BCAB52	BCAE52
	22cm	BNAA57		BCAM57	BCAD57	BCAGB7	BCAG57	BCAN57	BCAS57	BCAJ57	BCAB57	BCAE57
	24cm	BNAA53		BCAM53	BCAD53	BCAGB3	BCAG53	BCAN53	BCAS53	BCAJ53	BCAB53	BCAE53
	26cm	BNAA54		BCAM54	BCAD54	BCAGB4	BCAG54	BCAN54	BCAS54	BCAJ54	BCAB54	BCAE54
	28cm	BNAA55		BCAM55	BCAD55	BCAGB5	BCAG55	BCAN55	BCAS55	BCAJ55	BCAB55	BCAE55
	30cm	BNAA56		BCAM56	BCAD56	BCAGB6	BCAG56	BCAN56	BCAS56	BCAJ56	BCAB56	BCAE56
Stent 4.8 Ch/Fr + Pusher 6 Ch/Fr 40cm	12cm		BCAA40	BCAM40	BCAD40	BCAGC0	BCAG40	BCAN40	BCAS40	BCAJ40	BCAB40	BCAE40
	16cm		BCAA41	BCAM41	BCAD41	BCAGC1	BCAG41	BCAN41	BCAS41	BCAJ41	BCAB41	BCAE41
	20cm		BCAA42	BCAM42	BCAD42	BCAGC2	BCAG42	BCAN42	BCAS42	BCAJ42	BCAB42	BCAE42
	22cm		BCAA47	BCAM47	BCAD47	BCAGC7	BCAG47	BCAN47	BCAS47	BCAJ47	BCAB47	BCAE47
	24cm		BCAA43	BCAM43	BCAD43	BCAGC3	BCAG43	BCAN43	BCAS43	BCAJ43	BCAB43	BCAE43
	26cm		BCAA44	BCAM44	BCAD44	BCAGC4	BCAG44	BCAN44	BCAS44	BCAJ44	BCAB44	BCAE44
	28cm		BCAA45	BCAM45	BCAD45	BCAGC5	BCAG45	BCAN45	BCAS45	BCAJ45	BCAB45	BCAE45
	30cm		BCAA46	BCAM46	BCAD46	BCAGC6	BCAG46	BCAN46	BCAS46	BCAJ46	BCAB46	BCAE46
Stent 6 Ch/Fr + Pusher 6 Ch/Fr 40cm	20cm		BCAA62			BCAGD2	BCAG62		BCAS62	BCAJ62	BCAB62	
	22cm		BCAA67			BCAGD7	BCAG67		BCAS67	BCAJ67		
	24cm	BNAA63	BCAA63	BCAM63		BCAGD3	BCAG63		BCAS63	BCAJ63	BCAB63	BCAE63
	26cm	BNAA64	BCAA64	BCAM64	BCAD64	BCAGD4	BCAG64	BCAN64	BCAS64	BCAJ64	BCAB64	BCAE64
	28cm	BNAA65	BCAA65	BCAM65	BCAD65	BCAGD5	BCAG65	BCAN65	BCAS65	BCAJ65	BCAB65	BCAE65
	30cm		BCAA66		BCAD66	BCAGD6	BCAG66		BCAS66	BCAJ66	BCAB66	
Stent 7 Ch/Fr + Pusher 7 Ch/Fr 40cm	20cm		BCAA72			BCAGE2	BCAG72		BCAS72	BCAJ72		
	22cm		BCAA77			BCAGE7	BCAG77		BCAS77	BCAJ77	BCAB77	
	24cm	BNAA73	BCAA73	BCAM73		BCAGE3	BCAG73		BCAS73	BCAJ73	BCAB73	
	26cm	BNAA74	BCAA74	BCAM74	BCAD74	BCAGE4	BCAG74	BCAN74	BCAS74	BCAJ74	BCAB74	BCAE74
	28cm	BNAA75	BCAA75	BCAM75	BCAD75	BCAGE5	BCAG75	BCAN75	BCAS75	BCAJ75	BCAB75	BCAE75
	30cm		BCAA76		BCAD76	BCAGE6	BCAG76		BCAS76	BCAJ76	BCAB76	
Stent 8 Ch/Fr + Pusher 8 Ch/Fr 40 cm	26cm	BNAA84	BCAA84			BCAGF4	BCAG84		BCAS84	BCAJ84	BCAB84	
	28cm	BNAA85	BCAA85		BCAD85	BCAGF5	BCAG85		BCAS85	BCAJ85	BCAB85	
	30cm		BCAA86		BCAD86	BCAGF6	BCAG86		BCAS86	BCAJ86	BCAB86	
Stent 9 Ch/Fr + Pusher 8 Ch/Fr 40 cm	26cm		BCAA94									
	28cm		BCAA95		BCAD95							
	30cm		BCAA96		BCAD96							

* Only Italy Poseidon.

1 • ENDOUROLOGY

1.6. Double loop ureteral stents

Duration
of implantation:
up to 12 months.

Biosoft® duo

Multi-length hydro-coated double loop ureteral stents



Produced by coextrusion of a flexible material and a rigid material, Biosoft duo was designed for easy insertion, while retaining great flexibility for the patient comfort.

- **Kit composition:** each stent is supplied in a kit with a pusher (simple or steerable) and a guidewire.

Ø Guide: 0.035".

Steerable pusher: with radiopaque ring.

- **Box of 1 • Sterile.**



Kits with
simple pusher


			O/O
Ch/Fr Stent	Ch/Fr Pusher	Length Pusher	Fixed core guidewire
4.8	4.8	75cm	BNGA58
6	6	40cm	BNGA68
7	7	40cm	BNGA78




Kits with
steerable pusher

			O/O		
Ch/Fr Stent	Ch/Fr Pusher	Length Pusher	Fixed core guidewire	Terumo guidewire	Poseidon URO guidewire
4.8	6	40cm	BCGA58	BCGG58	BCGS58
6	6	40cm	BCGA68	BCGG68	BCGS68
7	7	40cm	BCGA78	BCGG78	BCGS78

The global stent is between 22 to 28 cm when the multi-length is straightened.





1 • ENDOUROLOGY

1.6. Double loop ureteral stents

Duration
of implantation:
up to 30 days.

PU-R Double loop ureteral stents

Double loop ureteral stents in rigid polyurethane

- **Kit composition:** each stent is supplied in a kit with a pusher (simple or steerable) and with or without a guidewire depending on the model.

Ø Guide: 0.035".

Steerable pusher: with or without radiopaque ring.

- **Box of 1** (with a guidewire) • Sterile.
- **Box of 5** (without a guidewire) • Sterile.

		Kits with simple pusher		Kits with steerable pusher		
		O/O	O/O			O/C
Ch/Fr	Length	Fixed core guidewire	Fixed core guidewire	Movable core guidewire	Without guidewire	Fixed core guidewire
6	24cm	AC4B63	AC4663	AC4D63		AC4563
	26cm	AC4B64	AC4664	AC4D64	AC4E64	AC4564
	28cm	AC4B65	AC4665	AC4D65	AC4E65	AC4565
	30cm	AC4B66	AC4666	AC4D66		AC4566
7	24cm	AC4B73	AC4673	AC4D73	AC4E73	
	26cm	AC4B74	AC4674	AC4D74	AC4E74	AC4574
	28cm	AC4B75	AC4675	AC4D75	AC4E75	AC4575
	30cm	AC4B76	AC4676		AC4E76	AC4576
8	28cm	AC4B85				
	30cm	AC4B86				

PU-S Double loop ureteral stents

Double loop ureteral stents in soft polyurethane

Duration
of implantation:
up to 90 days.

- **Kit composition:** each stent is supplied in a kit with a pusher (simple or steerable) and a guidewire.

Ø Guide: 0.035".

Steerable pusher: with or without radiopaque ring.

- **Box of 5** • Sterile.

		Kits with simple pusher		Kits with steerable pusher
		O/O		O/C
Ch/Fr	Length	Fixed core guidewire	Fixed core guidewire	Fixed core guidewire
6	24cm		AC4263	
	26cm	ACBH64	AC4264	AC4364
	28cm	ACBH65	AC4265	
7	24cm		AC4273	
	26cm	ACBH74	AC4274	AC4374
	28cm	ACBH75	AC4275	



Duration
of implantation:
up to 30 days.

Polyamide & Polyurethane double loop ureteral stent sets



- **Set composition:** each stent is supplied in the same peel pouch with a simple pusher and with or without a guidewire depending on the model.

Ø Guide: 0.035".

- **Box of 5** • Sterile.

Sets with
simple pusher

Ch/Fr	Length	O/O		O/C	
		Fixed core guidewire	Without guidewire	Fixed core guidewire	Without guidewire
4.8	24cm	BNBA53	BNBJ53	BNBB53	BNBK53
	26cm	BNBA54	BNBJ54	BNBB54	BNBK54
	28cm	BNBA55	BNBJ55	BNBB55	BNBK55
	30cm	BNBA56	BNBJ56		
6	24cm	BNBA63	BNBJ63	BNBB63	BNBK63
	26cm	BNBA64	BNBJ64	BNBB64	BNBK64
	28cm	BNBA65	BNBJ65	BNBB65	BNBK65
	30cm	BNBA66	BNBJ66		
7	24cm	BNBA73	BNBJ73	BNBB73	BNBK73
	26cm	BNBA74	BNBJ74	BNBB74	BNBK74
	28cm	BNBA75	BNBJ75	BNBB75	BNBK75
	30cm	BNBA76	BNBJ76		
8	26cm	BNBA84	BNBJ84	BNBB84	BNBK84
	28cm	BNBA85	BNBJ85	BNBB85	BNBK85
	30cm	BNBA86	BNBJ86	BNBB86	BNBK86

Duration
of implantation:
up to 12 months.

Double loop ureteral stents

For ureteral stenosis

STENOSTENT • Silicone for long-term implantation.
12 Ch/Fr reinforced body for maximum resistance to stenosis, 8 Ch/Fr loops for patient comfort.

- **Kit composition:** each stent is supplied in a kit with a steerable pusher (with radiopaque ring) and a guidewire.
- **Ø Guide:** 0.035" (fixed core guidewire).
- **Box of 1** • Sterile.




Kits with
steerable pusher

	Q/O				
Length	16cm	24cm	26cm	28cm	30cm
Order Ref.	AJ4W81	AJ4W83	AJ4W84	AJ4W85	AJ4W86

Duration
of implantation:
up to 12 months.

Double loop ureteral stents

For pyeloplasty_

Pyelostent silicone for long-term implantation.

12 Ch/Fr reinforced part in renal pelvis for a better healing, 8 Ch/Fr loops for patient comfort.

- **Kit composition:** each stent is supplied in a kit with a steerable pusher (with radiopaque ring) and a guidewire.
Ø Guide: 0.035" (fixed core guidewire).
- **Box of 1** • Sterile.



Kits with
steerable pusher

	O/O		
Length	26cm	28cm	30cm
Order Ref.	AJ4Y84	AJ4Y85	AJ4Y86

1 • ENDOUROLOGY

1.6. Double loop ureteral stents

Duration
of implantation:
up to 12 months.

Double loop ureteral stents

Without eyes on the straight section

Vortek: its two-layer structure makes it easy to insert.

Silicone: the biocompatibility of silicone makes it the material of choice for long-term implantation.

- **Kit composition:** each stent is supplied in a kit with a pusher (simple or steerable) and a guidewire.

Ø Guide: 0.035" (fixed core guidewire).

Steerable pusher: with or without radiopaque ring.

- **Box of 1** • Sterile.



Sets with
simple pusher

		O/O Vortek	
Ch/Fr	Length	Steerable pusher	Simple pusher
6	16cm		
	20cm		
	22cm		
	24cm	ACB763	
	26cm	ACB764	ACB864
	28cm	ACB765	ACB865
	30cm	ACB766	ACB866
7	16cm		
	20cm		
	22cm		
	24cm	ACB773	
	26cm	ACB774	ACB874
	28cm	ACB775	ACB875
	30cm	ACB776	ACB876
8	16cm		
	20cm		
	22cm		
	24cm	ACB783	
	26cm	ACB784	
	28cm	ACB785	
	30cm	ACB786	

		O/O Silicone	
		Steerable pusher	Simple pusher
		AJ4U61	AJ4B61
		AJ4U62	AJ4B62
		AJ4U67	AJ4B67
		AJ4U63	AJ4B63
		AJ4U64	AJ4B64
		AJ4U65	AJ4B65
		AJ4U66	AJ4B66
		AJ4U71	AJ4B71
		AJ4U72	AJ4B72
		AJ4U77	AJ4B77
		AJ4U73	AJ4B73
		AJ4U74	AJ4B74
		AJ4U75	AJ4B75
		AJ4U76	AJ4B76
		AJ4U81	AJ4B81
		AJ4U82	AJ4B82
		AJ4U87	AJ4B87
		AJ4U83	AJ4B83
		AJ4U84	AJ4B84
		AJ4U85	AJ4B85
		AJ4U86	AJ4B86

Duration
of implantation:
up to 12 months.

Vortek® Tumor stents

Without eyes on the straight section

The inner layer is reinforced compared to the standard Vortek ureteral stent. It allows passing through stricture and the stiffness of material offers an excellent resistance to compression.

- **Kit composition:** each stent is supplied in a kit with a steerable pusher (with radiopaque ring) and a guidewire.
Ø Guide: 0.035" (hydrophilic nitinol guidewire).
- **Box of 1** • Sterile.




Kits with
steerable pusher

		O/O	
Ch/Fr	Length	Terumo guidewire	Poseidon® guidewire
7	26cm	BCCG74	BCCS74
	28cm	BCCG75	BCCS75
	30cm	BCCG76	BCCS76
8	26cm	BCCG84	BCCS84
	28cm	BCCG85	BCCS85
	30cm	BCCG86	BCCS86

Vortek® Single loop ureteral stents

Device for drainage and urine collection

- **Set composition:**
 - Single loop ureteral stent 90cm.
 - Seldinger fixed core 0.035" guidewire.
 - Clamp.
 - A connector with a Luer end for connection to a Luer-bag.
- **Box of 1** • Sterile.



Ch/Fr	Q/C		
	6	7	8
Order Ref.	AC4406	AC4407	AC4408

Eyes on loop and body.

Ch/Fr	Q/O		
	6	7	8
Order Ref.	ACA206	ACA207	ACA208

Eyes on loop and body.

Ch/Fr	Q/O		
	6	7	8
Order Ref.	ACA106	ACA107	ACA108

Eyes on loop only.

1 • ENDOUROLOGY

1.8. Guidewires

Soprano®

Exclusive Hybrid Design

Soprano® is a guidewire with an exclusive hybrid design that combines the access of a nitinol wire with the handling of a PTFE-coated wire

Soprano®

The hydrophilic flexible tip reduces surface friction, and its round distal tip creates a lubricious surface for ease of insertion and navigation.



Soprano® Dual Flex

The Soprano Dual Flex features a flexible proximal end designed to ease passage and to minimize the risk of scope damage.



Order Ref.	Description	Diameter	Length	Tip shape	Tip style	Quantity
AEHA35	Soprano®	0.035" (0.89mm)	150cm	Straight	5cm hydrophilic straight distal tip with the first 3cm flexible	Box of 5
AEHB35	Soprano® Dual Flex	0.035" (0.89mm)	150 cm	Straight	5cm hydrophilic straight distal tip with the first 3cm flexible 7.5cm flexible (proximal)	Box of 5

Hydrophilic Nitinol guidewires

Radifocus™ Terumo guidewires*

Its configuration and the combination of Nitinol and hydrophilic properties allow the Terumo guidewires to reach difficult access sites.

- Flexible tip length: 3cm.
- Box of 5 • Sterile.


Diameter	Length	Rigidity	Tip	Order Ref.
0.035"	150cm	Standard	Straight	AE0N35
			Angled	AE0T35
		Stiff	Straight	AE0P35
			Angled	AE0X35



Poseidon guidewires**

A hydrophilic coating is applied to the guide wire to lower the friction and improve the navigation and trackability.

- Flexible tip length: 3cm.
- Box of 5 • Sterile.



Diameter	Length	Rigidity	Tip	Order Ref.
0.035"	150cm	Standard	Straight	000P1HYDS35
			Angled	000P1HYDA35
		Stiff	Straight	000P1HYSS35
			Angled	000P1HYSA35



* Not authorized in Italy.
** For Italy only.

PTFE-coated guidewires

Seldinger guidewires

- Box of 5 • Sterile.

Fixed core guidewires

Straight tip



	150cm	150cm	90cm
Diameter	Standard	Stiff	Standard
0.032"	AE0032		AE0E32
0.035"	AE0035	AE0C35	
0.038"	AE0038		AE0E38

5cm flexible tip except AE0C35 11cm.

Fixed core guidewires

"J" tip



	80cm
Diameter	Stiff
0.038"	AE1138

5cm flexible tip.

Movable core guidewires

Straight tip



	150cm	
Diameter	Standard	Stiff
0.032"	AE0A32	
0.035"	AE0A35	AE0D35
0.038"	AE0A38	

5cm flexible tip.

Stainless steel guidewires

Schüller and Lunderquist guidewires

• Box of 5 • Sterile.



Schüller

Diameter	100cm	
	Straight	"J"
	AE0435	AEJ435

7cm flexible tip.



Lunderquist

Diameter	80cm	
	Straight	"J"
	AE0535	AEJ535

6cm flexible tip for AE0535.
7cm flexible tip for AEJ535.

1 • ENDOUROLOGY

1.9. Ureteral catheters

Drainage catheters

Wide range of tips

- Each catheter is supplied with a stylet (metal or polyester) and a colour code connector.
- Material: Neoplex.
- Length: 70cm for AC55xx, AC52xx, ACN6xx and 69cm for AC54xx, AC51xx, AC53xx and AC50xx.
- Catheters are individually packed.
- **Box of 10** • Sterile.

Straight open whistle tip

Eyes: 2 - Stylet: metal							
Ch/Fr	3	4	5	6	7	8	9
Order Ref.	AC5503	AC5504	AC5505	AC5506	AC5507	AC5508	AC5509

Eyes: 6 - Stylet: metal							
Ch/Fr	3	4	5	6	7	8	9
Order Ref.	ACN603	ACN604	ACN605	ACN606	ACN607	ACN608	ACN609

Eyes: 2 - Stylet: polyester							
Ch/Fr	3	4	5	6	7	8	9
Order Ref.	AC5203	AC5204	AC5205	AC5206	AC5207	AC5208	AC5209



Straight closed olive tip

Eyes: 2 - Stylet: metal						
Ch/Fr	3	4	5	6	7	8
Order Ref.	AC5403	AC5404	AC5405	AC5406	AC5407	AC5408

Eyes: 2 - Stylet: polyester						
Ch/Fr	3	4	5	6	7	8
Order Ref.	AC5103	AC5104	AC5105	AC5106	AC5107	AC5108



Straight closed cylindrical tip

Eyes: 2 - Stylet: metal						
Ch/Fr	3	4	5	6	7	8
Order Ref.	AC5303	AC5304	AC5305	AC5306	AC5307	AC5308

Eyes: 2 - Stylet: polyester						
Ch/Fr	3	4	5	6	7	8
Order Ref.	AC5003	AC5004	AC5005	AC5006	AC5007	AC5008



Interventional catheters

These catheters are suitable for use with a guidewire during endoscopic procedures

Made of PEBA.

- Each catheter is supplied with a metal stylet and a Luer syringe connector.
- Length: 74cm for ACP20x, APC30x and 73cm for ACP40x, ACP60x.
- Catheters are individually packed.
- **Box of 5** • Sterile.

Flush straight open tip ureteric catheters without side eye

Ch/Fr	4.8	6	7	8
Guidewire recommended	0.038"	0.038"	0.038"	0.038"
Order Ref.	ACP205	ACP206	ACP207	ACP208



Straight open tip ureteric catheters with 2 side eyes

Ch/Fr	3	4	4.8	6	7
Guidewire recommended	0.021"	0.035"	0.038"	0.038"	0.038"
Order Ref.	ACP303	ACP304	ACP305	ACP306	ACP307




Open coude tip ureteric catheters with 2 eyes

Ch/Fr	3	4	4.8	6	7
Guidewire recommended	0.021"	0.035"	0.038"	0.038"	0.038"
Order Ref.	ACP403	ACP404	ACP405	ACP406	ACP407



Open coude tip ureteric catheters without side eye

Ch/Fr	3	4	4.8	6	7
Guidewire recommended	0.021"	0.035"	0.038"	0.038"	0.038"
Order Ref.	ACP603	ACP604	ACP605	ACP606	ACP607



1 • ENDOUROLOGY


1.9. Ureteral catheters

Interventional catheters

Floppy tip hydro-coated ureteric catheter

The flexible open tip is designed for smooth passage through the ureter.
The hydrophilic coating over 15cm makes catheter insertion and advancement easier.
The grey part at the distal end improves the visualisation of the catheter tip.

- Each catheter is supplied with a polyester stylet and a Luer syringe connector.
- Material: Neoplex.
- Length: 70cm.
- Catheters are individually packed.
- **Box of 5** • Sterile.



Straight open ureteric catheters without side eye

Ch/Fr	5	6	7
Guidewire recommended	0.035"	0.035"	0.035"
Order Ref.	ACR205	ACR206	ACR207

Chevassu catheters for Retrograde UreteroPyelography (RUP)

The shaped bulb of Chevassu catheters allows effective injection of contrast medium.

- Each catheter is supplied with a metal stylet and a Luer syringe connector.
- Material: Neoplex.
- 73cm for AC5C07 and 70cm for AC59xx.
- Catheters are individually packed.
- **Box of 5** • Sterile.

Chevassu ureteric catheters compatible with a 0.035" guidewire

Bulb Ch/Fr	7
Body Ch/Fr	5
Order Ref.	AC5C07



Chevassu ureteric catheters


Bulb Ch/Fr	4	5	6	7	8	9	10	11	12
Body Ch/Fr	3	3	4	5	5	6	7	7	7
Order Ref.	AC5904	AC5905	AC5906	AC5907	AC5908	AC5909	AC5910	AC5911	AC5912



Ureteral dilators

Used for dilatation of the ureter during ureteroscopic procedures

- Compatible with a 0.038" guidewire.
- Diameter (Ch/Fr) indicated at each tip.
- **Box of 10** • Sterile.



Order Ref.	Ch/Fr	Length
RBD010	8 and 10	48cm
RBD014	12 and 14	48cm


2 • PERCUTANEOUS ACCESS

2.1. Percutaneous nephrostomy set

Kolibri® set with J catheter

Complete range of drainage set with J catheter

- **Kit composition:**
 - 2 Chiba needles: 18 and 22G.
 - 1 Seldinger or Schüller J-tip guidewire.
 - 1 set of dilators according to the catheter size with or without working sheath.
 - 1 Vortek J catheter with sleeve (29cm) except for 14Ch/Fr (32cm).
 - 1 Luer-lock stopcock and PVC connector tubing.
- **Box of 1 • Sterile.**




		Ch/Fr						
J catheter		6	7	8	9	10	12	14
Seldinger guidewire 0.038"	Simple dilator	RJA006*		RJA008* RJA108		RJA110	RJA112	RJA114
	Simple dilators + dilator with working sheath			RJA208		RJA210	RJA212	
Schüller guidewire 0.035"	Simple dilator	RJB006*	RJB007*	RJB108	RJB109	RJB110		
	Simple dilators + dilator with working sheath			RJB208		RJB210		

* Catheter without stylet.

Kolibri® set with Malecot catheter

Complete range of drainage set with Malecot catheter

- **Kit composition:**
 - 2 Chiba needles: 18 and 22G.
 - 1 Seldinger J-tip guidewire 0.038".
 - 1 set of dilators according to the catheter size with or without working sheath.
 - 1 Vortek Malecot catheter with sleeve (36cm with the connector, except for 8Ch/Fr 34cm).
 - 1 Luer-lock stopcock and PVC connector tubing.
- **Individually packaged** • Sterile.



		Ch/Fr			
Malecot catheter		8	10	12	14
Seldinger guidewire	Simple dilator	RMA108	RMA110	RMA112	RMA114
	Simple dilators + dilator with working sheath	RMA208	RMA210	RMA212	RMA214


2 • PERCUTANEOUS ACCESS

2.2. Percutaneous nephrostomy accessories


Guidewires

Catheter insertion and changing_


- Box of 5 • Sterile.



Type	Tip	Diameter	Length	Order Ref.
Seldinger (stiff)	"J"	0.038"	80cm	AE1138
	Straight	0.035"	150cm	AE0C35



Type	Tip	Diameter	Length	Order Ref.
Schüller	"J"	0.035"	100cm	AEJ435
	Straight	0.035"	100cm	AE0435



Type	Tip	Diameter	Length	Order Ref.
Lunderquist	"J"	0.035"	80cm	AEJ535
	Straight	0.035"	80cm	AE0535


Needles

Puncture of the renal cavities

- Box of 5 • Sterile.

Order Ref.	Tip	Diameter
RAC018	Chiba	18G
RAC022	Chiba	22G

Length: 20cm.




Dilators

Tract preparation simple or splittable working sheath to protect the parenchyma and make withdrawal easier

- Box of 5 • Sterile.



	Simple dilator				
Ch/Fr	8	10	12	6, 8, 10	6, 8, 10, 12, 14
Order Ref.	RBA008	RBA010	RBA012	RBB010	RBB014



	Simple dilator and dilator with working splittable sheath			
	Simple	With working sheath	Simple	With working sheath
Ch/Fr	6	8	6 and 8	10
Order Ref.	RBC008		RBC010	

2 • PERCUTANEOUS ACCESS

2.2. Percutaneous nephrostomy accessories

J and Malecot catheters

For percutaneous drainage of the renal cavities_

- Box of 5 • Sterile.

J catheters

J catheters are supplied with sleeve + 1 Luer-lock stopcock and PVC connector tubing.

	Without stylet			With stylet				
Ch/Fr	6	7	8	8	9	10	12	14
Order Ref.	RCJ006	RCJ007	RCJ008	RCJ108	RCJ109	RCJ110	RCJ112	RCJ114

Length: without connector: 27cm except for RCJ114: 30cm.

Malecot catheters

Malecot catheter is supplied with sleeve + 1 Luer-lock stopcock and PVC connector tubing. Except RCM016/18/20

	With stylet						
Ch/Fr	8	10	12	14	16	18	20
Order Ref.	RCM008	RCM010	RCM012	RCM014	RCM016	RCM018	RCM020

Length: without connector 36cm except 34cm for RCM008 and RCM016/18/20.

Balloon catheters in silicone

For increased comfort during drainage_

- All silicone for biocompatibility and patient comfort.
- Radiopaque line helps to see product on X-Ray.
- Box of 10 • Sterile.

Ch/Fr	10	12	14	16	18	20
Catheter with 2 eyes above the balloon	AJ6610	AJ6612	AJ6614	AJ6616	AJ6618	AJ6620
Catheter with 2 eyes above and under the balloon	AJ6710	AJ6712	AJ6714	AJ6716	AJ6718	AJ6720

Length: 40cm.

Connection

- Box of 10 • Sterile.

Removable Luer-lock and universal urine bag adapter
REA000



Set with J catheter

Percutaneous set. Each component of the kit is supplied in the same peel pouch_

- **Set composition:**
 - Chiba needle 18G.
 - Seldinger guidewire 0.038".
 - 1 set of dilators depending on the catheter size.
 - Vortek J catheter with stylet (except RJE106).
 - Length:29cm (except RJE114: 32.5cm).
- **Box of 5 • Sterile.**

Ch/Fr	6	8	10	12	14
Order Ref.	RJE106	RJE108	RJE110	RJE112	RJE114

2 • PERCUTANEOUS ACCESS

2.4. Suprapubic drainage

Cystodrain® integral Set

Suprapubic drainage sets with Silicone Integrated balloon catheter

Punction set	Light punction set	Dilation set	Replacement set
Silicone Integral balloon catheter	Silicone Integral balloon catheter	Silicone Integral balloon catheter	Silicone Integral balloon catheter
Splittable trocar (12cm)	Splittable trocar (12cm)	Dilator	Adhesive tape
Locking sleeve		PTFE Guidewire (90 cm)	Spigot
Adhesive tape		Adhesive tape	Optional PTFE guidewire (90 cm)
Spigot		Spigot	
Scalpel			
Box of 5 • Sterile	Box of 5 • Sterile	Box of 5 • Sterile	Box of 1 or 5 • Sterile

Silicone Integral balloon catheter: Channels: 2 – Eyes: 2 – Balloon size: 3 (Ch.10) / 5 ml – PTFE Guidewire: 0.038" (0.97mm).



Punction sets references

Trocar Ø int/ext Ch,Fr	12/14	13/16	16/17
Catheter Ø Ch/Fr	10	12	14
Punction set Order Ref.	CD1010	CD1012	CD1014
Light punction set Order Ref.	CD1210	CD1212	CD1214

Dilatation set references

Catheter Ch/Fr	12	14
Set Order Ref.	CD1112	CD1114

Simple replacement sets

Catheter Ch/Fr	10	12	14	16	18	20
Replacement set Order Ref. (with guidewire)	CD1A10	CD1A12	CD1A14	CD1A16		
Replacement set Order Ref. (without guidewire)	CD1B10	CD1B12	CD1B14	CD1B16	CD1B18	CD1B20
	CD1C10	CD1C12	CD1C14	CD1C16	CD1C18	CD1C20

Box of 1

Box of 5

Supraflow


The benefits of suprapubic puncture
with the simplicity of a balloon catheter_

- Complete set.
- Silicone Foley catheter.
- A peelable sheath trocar.
- Scalpel.
- **Box of 5** • Sterile.

Complete set references

Trocar Ø int/ext mm	5.5/6.8	6.75/8.0
Catheter Ø Ch/Fr	12	16
Ballon size (ml)	10	15
Set Order Ref.	AJ9212	AJ9216

Silicone Foley catheter additional info: Length: 40cm – Channels: 2 – Eyes: 2.



2 • PERCUTANEOUS ACCESS

2.5. Urinary Diversion

Urinary diversion

- Individually packaged, except for AC41xx and AC65xx (box of 5) • Sterile.



Length: 155cm • Material: PVC								
Ch/Fr	5	6	7	8	9	10	11	12
Order Ref.	AC4105	AC4106	AC4107	AC4108	AC4109	AC4110	AC4111	AC4112

Catheter for uretero-sigmoidostomy, 6 eyes, straight proximal end, radiopaque line.



Length: 120cm • Material: PVC				
Ch/Fr	8	10	12	14
Order Ref.	AC4708	AC4710	AC4712	AC4714

Catheter for entero-cystoplasty, with eyes, funnel proximal end, radiopaque line.



Length: 45cm • Material: Silicone								
Ch/Fr	8	10	12	14	16	18	20	22
Order Ref.	AC6708	AC6710	AC6712	AC6714	AC6716	AC6718	AC6720	AC6722

Ureterostomy catheter, 5 eyes, bevel tip, funnel proximal end, radiopaque line. Round flange.



Length: 45cm • Material: Silicone							
Ch/Fr	8	10	12	14	16	20	22
Order Ref.	AC6808	AC6810	AC6812	AC6814	AC6816	AC6820	AC6822

Ureterostomy catheter, 5 eyes, bevel tip, funnel proximal end, radiopaque line. Oval flange.



Length: 45cm • Material: Silicone			
Ch/Fr	8	10	12
Order Ref.	AC7008	AC7010	AC7012

Ureterostomy catheter, 5 eyes, bevel tip, straight proximal end, radiopaque line. Round flange.



Total Length: 52cm • Material: Silicone			
Ch/Fr	8	10	12
Order Ref.	BJ3008	BJ3010	BJ3012

Ureterostomy and entero-cystoplasty catheter, single loop, open and tapered distal tip, radiopaque. Oval flange.

Connectors and spigot

- Box of 10 • Sterile.
- CP18200000: Box of 50 • Sterile.

Connectors for tube



Material: Polypropylene			
External diameter (mm)	4-6-4	7-11-7	10-14-10
Internal diameter (mm)	3-3	5-5	8-8
Order Ref.	PB0101	PB0102	PB0103



Material: Polypropylene			
External diameter (mm)	4-6 / 7-11	4-6 / 14-10	7-11 / 14-10
Internal diameter (mm)	3-5	3-8	5-8
Order Ref.	PB0201	PB0202	PB0203



Material: Polypropylene				
External diameter (mm)	3 x 4-5	3 x 6-8	3 x 9-11	3 x 11-13
Internal diameter (mm)	2-2-2	4-4-4	7-7-7	8-8-8
Order Ref.	PB1201	PB1203	PB1204	PB1205



Tuohy Borst Adapter

Material: Polycarbonate • Acrylic • Silicone	
Order Ref.	AK2000

Spigot



Material: ABS	
Order Ref.	CP18200000

3 • CONNECTORS AND SPIGOT

3.1. Connectors and Spigot

Connectors for catheters



Material: Thermoplastic Elastomer Styrene (TPE)		
Ch/Fr	3-10	9-12
Order Ref.	AK4010	AK4012



Material: PVC	
Order Ref.	AK5300



Material : Thermoplastic Elastomer Styrene (Sheath) Acrylic copolymer (Luer male) • Thermoplastic (Luer female)	
Order Ref.	AK4100



Material: Thermoplastic Elastomer Styrene (Sheath) Acrylic copolymer (Luer male) • Thermoplastic and Silicone (Luer connector)		
Ch/Fr	3-8	9-12
Order Ref.	AK4200	AK4201



Material: Thermoplastic Elastomer Styrene (Sheath) • Acrylic copolymer (Luer male)	
Order Ref.	AK4300



Material: Polyamide • Silicone		
Ch/Fr	3-8	9-12
Order Ref.	AK3400	AK3401

Dufour tip



- Length: 41cm.
- Box of 5 • Sterile.



Silicone X-Flow
Hydro X-Flow



Silicone X-Flow
Short Tip

Ways	Balloon size	Material	Ch/Fr			
			18	20	22	24
3 Ways	30ml	X-Flow Silicone Dufour short tip	AB6S18	AB6S20	AB6S22	AB6S24
	50ml	X-Flow Silicone Dufour short tip	AB6R18	AB6R20	AB6R22	AB6R24
		X-Flow Silicone Dufour tip	AB6318	AB6320	AB6322	AB6324
		Hydro X-Flow Silicone Dufour tip	XB6318	XB6320	XB6322	XB6324
	80ml	Hydro X-Flow Silicone Dufour tip	XB6L18	XB6L20	XB6L22	XB6L24
2 Ways	50ml	X-Flow Silicone Dufour tip	AB6418	AB6420	AB6422	AB6424

4 • LOWER URINARY TRACT
 4.1. Prostatic catheters

Delinotte tip
 Coudé tip_

- Length: 41cm.
- **Box of 5** • Sterile.



Silicone X-Flow
 Hydro X-Flow

			Ch/Fr			
Ways	Balloon size	Material	18	20	22	24
3 Ways	50ml	X-Flow Silicone Delinotte tip	AB6218	AB6220	AB6222	AB6224
		Hydro X-Flow Silicone Delinotte tip	XB6218	XB6220	XB6222	XB6224
	80ml	Hydro X-Flow Silicone Delinotte tip	XB6N18	XB6N20	XB6N22	XB6N24

Couvelaire tip

Whistle tip_

- Length: 41cm.
- Box of 5 • Sterile.



Silicone X-Flow
Hydro X-Flow

			Ch/Fr			
Ways	Balloon size	Material	18	20	22	24
3 Ways	50ml	X-Flow Silicone Couvelaire tip	AB6118	AB6120	AB6122	AB6124
		Hydro X-Flow Silicone Couvelaire tip	XB6118	XB6120	XB6122	XB6124
	80ml	Hydro X-Flow Silicone Couvelaire tip	XB6M18	XB6M20	XB6M22	XB6M24
2 Ways	50ml	X-Flow Silicone Couvelaire tip	AB6518	AB6520	AB6522	AB6524

4 • LOWER URINARY TRACT

4.1. Prostatic catheters

Straight tip



- Length: 41cm.
- **Box of 5** • Sterile.



Silicone

X-Flow



Silicone short tip

X-Flow



Over guidewire

			Ch/Fr			
Ways	Balloon size	Material	18	20	22	24
3 Ways	50ml	X-Flow Silicone straight tip	AB6018	AB6020	AB6022	AB6024
		X-Flow Silicone straight short tip	AB6H18	AB6H20	AB6H22	AB6H24
		X-Flow Silicone over the guide short tip*	AB6G18	AB6G20	AB6G22	AB6G24

* Compatible with guidewire.

Freudenberg Stylet

- Intended to stiffen or elbow a prostatic catheter in case of urethral insertion difficulty.
- **Individually packaged** • Sterile • Reusable.



Freudenberg stylet	
Order Ref.	AL0000

4 • LOWER URINARY TRACT

4.1. Prostatic catheters

Catheters

With 30ml balloon



- Sterile.



Silicone
X-Flow



Silicone short tip
X-Flow

With 30ml balloon • Box of 5

- Length 41cm.

				Ch/Fr			
Ways	Balloon size	Tip	Material	18	20	22	24
3 Ways	30ml	Dufour	X-Flow Silicone	AB6A18	AB6A20	AB6A22	AB6A24
		Delinotte	X-Flow Silicone	AB6E18	AB6E20	AB6E22	AB6E24
		Couvelaire	X-Flow Silicone	AB6B18	AB6B20	AB6B22	AB6B24
		Straight	X-Flow Silicone	AB6C18	AB6C20	AB6C22	AB6C24
		Straight short tip	X-Flow Silicone	AB6J18	AB6J20	AB6J22	AB6J24

ABVac

Bladder evacuator

Intended to be used for evacuating irrigation solutions, blood clots, tissue collection and/or straining functions during prostate surgery or bladder surgery.

Main Features

- Single Nozzle adaptor for Storz, Wolf, Olympus and ACMI Resectoscopes.
- Single hand use for an excellent eye - hand coordination.
- Excellent throw and suction.

Order Ref.	Description	Quantity
000PABVAC	Bladder evacuator	10 units



4 • LOWER URINARY TRACT
4.3. Bladder injection

BoNee®
Needle for bladder injection



- For Injection inside the bladder muscle.
- Short tip (4mm).
- 22G needle
- 5Ch/Fr tube.
- Protection cap on the needle tip.
- Pointer stop: Reduce the bladder perforation risk.
- Flexible body.
- **Box of 1 or 5** • Sterile.

Tube length	Use	Box of 1	Box of 5
35cm	With rigid cystoscope	NBI035	BBI035
70cm	With both flexible and rigid cystoscope	NBI070	BBI070

22G needle, 5 Ch/Fr tube, 4mm tip.



Bougies

- Dilation of the meatus before endoscopy.
- **Individually packaged** • Gradual dilation of strictures • Sterile.

Dilation of the meatus before endoscopy and gradual dilation of strictures are applicable for all bougies

For routine urethral dilation

Neoplex • Straight / tapered olive tip bougie - Length: 34cm								
Ch/Fr	6	7	8	9	10	12	14	16
Order Ref.	AG5306	AG5307	AG5308	AG5309	AG5310	AG5312	AG5314	AG5316
Ch/Fr	18	20	22	24	26	28	30	
Order Ref.	AG5318	AG5320	AG5322	AG5324	AG5326	AG5328	AG5330	



Bougies metric thread

For difficult-to-cross strictures (Order Ref. AG51)_
 Philips connectable bougies (Order Ref. AG54 only)_

- **Individually packaged** • Sterile.

Filiform bougies

		Male metric thread • Neoplex			
		Ch/Fr			
Tip	Length	4	5	6	7
Straight	30cm	AG5104	AG5105	AG5106	AG5107

Philips Bougies

Female metric thread • Neoplex - Length: 34cm									
Ch/Fr	8	10	12	14	16	18	20	22	24
Order Ref.	AG5408	AG5410	AG5412	AG5414	AG5416	AG5418	AG5420	AG5422	AG5424



Foley catheters

Silicone catheters



All lengths are indicative.



Tiemann



Tiemann+



Over guidewire



Folysil X-Tra



Neobladder 3 ways

• Box of 5 • Sterile.

Type	Balloon	Total length	Number of eyes	Ch/Fr						
				12	14	16	18	20	22	24
Straight	10ml	40cm	2	AW6112	AW6114	AW6116	AW6118	AW6120	AW6122	AW6124
	10ml	41cm	2	AA6112	AA6114	AA6116	AA6118	AA6120	AA6122	AA6124
	30ml	41cm	2				AA6C18	AA6C20	AA6C22	AA6C24
Female Straight	10ml	25cm	2	AA6512	AA6514					
	15ml	25cm	2			AA6516	AA6518	AA6520	AA6522	AA6524
Over the guidewire	10ml	40cm	2 + tip	AA6412	AA6414					
	15ml	40cm	2 + tip			AA6416	AA6418	AA6420	AA6422	AA6424
Straight Grooved	10ml	40cm	2	AA6612	AA6614					
	15ml	40cm	2			AA6616	AA6618	AA6620	AA6622	AA6624
Tiemann	10ml	40cm	1	AA6312	AA6314					
	15ml	40cm	1			AA6316	AA6318	AA6320	AA6322	AA6324
	10ml	42cm	1	AW6312	AW6314	AW6316	AW6318	AW6320	AW6322	AW6324
Tiemann+	10ml	40cm	1	AA9312	AA9314					
	15ml	40cm	1			AA9316	AA9318	AA9320	AA9322	
Neobladder 3 ways	30ml	40cm	6 + tip					AA6820	AA6822	
	30ml	40cm	4 + 2** + tip				AA6918	AA6920	AA6922	

** 2 eyes under the balloon.

10ml balloon

* Delivered without introducer.

4 • LOWER URINARY TRACT

4.5. Foley catheters

All lengths are indicative.

Paediatric foley catheters, up to 30 days, except LT* use

• Box of 5 • Sterile.

Type	Balloon	Length	Ch/Fr		
			6	8	10
Straight (with introducer)	1.5ml	26cm	AA6106		
	1.5ml	31cm	AW6106		
	3ml	26cm		AA6108	
	3ml	31cm		AW6108	AA6110 / AW6110
Over the guidewire	1.5ml	26cm	AA6406		
	3ml	26cm		AA6408	
	3ml	31cm			AA6410
Tiemann	3ml	27cm		AA6308	
	3ml	32cm			AA6310
*Long Term	1.5ml	26cm	AA7106		
	3ml	26cm		AA7108	
	3ml	31cm			AA7110

Folysil LT, for long term use (up to 12 weeks)

• Box of 5 • Sterile.

Type	Balloon	Length	Number of eyes	Ch/Fr					
				12	14	16	18	20	22
Straight	10ml	41cm	2	AA7112	AA7114	AA7116	AA7118	AA7120	AA7122
Female	10ml	25cm	2	AA7512	AA7514	AA7516	AA7518		
Over guidewire	10ml	40cm	2 + tip	AA7412	AA7414	AA7416	AA7418	AA7420	

Folysil X-tra, with prefilled syringe for long term use (up to 12 weeks)

• Box of 5 • Sterile.

Type	Balloon	Length	Number of eyes	Ch/Fr					
				12	14	16	18	20	22
X-Tra Straight	10ml	41cm	2	AA8112	AA8114	AA8116	AA8118	AA8120	AA8122
X-Tra Female	10ml	25cm	2	AA8512	AA8514	AA8516	AA8518		

Folysil X-Tra single pack, with prefilled syringe and empty syringe, for long term use (up to 12 weeks)

• Box of 1 • Sterile.

Type	Balloon	Length	Number of eyes	Ch/Fr					
				12	14	16	18	20	22
X-Tra Straight	10ml	41cm	2	AA8A12	AA8A14	AA8A16	AA8A18	AA8A20	AA8A22
X-Tra Female	10ml	25cm	2	AA8B12	AA8B14	AA8B16	AA8B18		
X-Tra over guidewire	10ml	40cm	2 + tip	AA8C12	AA8C14	AA8C16	AA8C18		

Catheter Valve

Intended to be connected to an indwelling catheter, to drain the urine_

• Box of 10 • Sterile.

Order Ref.

CV3808



4 • LOWER URINARY TRACT

4.6. Catheters without balloon

Urethral catheters without balloon

• Box of 5 • Sterile.



Couvelaire (Neoplex)

Straight whistle tip catheters

				Ch/Fr										
Type	Material	Length	Eyes	6	8	10	12	14	16	18	20	22	24	26
Couvelaire	Neoplex	40cm	2	AD5D06	AD5D08	AD5D10	AD5D12	AD5D14	AD5D16	AD5D18	AD5D20	AD5D22	AD5D24	AD5D26



Mercier

Coude tip catheters

				Ch/Fr										
Type	Material	Length	Eyes	6	8	10	12	14	16	18	20	22	24	26
Mercier	Neoplex	40cm	2	ADN306	ADN308	ADN310	ADN312	ADN314	ADN316	ADN318	ADN320	ADN322	ADN324	ADN326

Prostatic Catheters tip terminology



Dufour tip

Coudé whistle tip



Delinotte tip

Coudé cylindrical tip



Couvelaire tip

Straight whistle tip.



Straight tip



Our mission

Making life easier for people with intimate healthcare needs

CARJ0CA 2024229

Ostomy Care | Continence Care | Advanced Wound Care | Interventional Urology | Voice and Respiratory Care

Coloplast and the Coloplast logo are trademarks of Coloplast A/S.
www.coloplast.com © 2025 Coloplast A/S. All rights reserved. PM-37190.

