

RenalGuard

by RenalGuard Solutions, Inc.

RENALGUARD THERAPY®

A Therapeutic Approach to Reducing Contrast-Induced Nephropathy

RenalGuard Therapy® is designed to reduce the toxic effects that contrast media can have on the kidneys, which may lead to a reduction in the incidence of Contrast-Induced Nephropathy (CIN) in at-risk patients. This therapy is based on existing published literature, including the PRINCE study¹, which support the theory that creating and maintaining a high urine output through the kidneys allows the body to rapidly eliminate contrast, reducing its toxic effects.

RenalGuard Therapy entails the use of a physician-prescribed loop diuretic which may induce the required high urine output. RenalGuard System is designed to measure urine output and replaces it in real-time with an equal volume of sterile saline. This matched fluid replacement aims to minimize the risk of over- or under-hydration which can lead to increased patient risks.

Please contact us at EPS Vascular for more information about RenalGuard.



your niche. our core.

Please contact us or visit www.EPSvascular.com for more information.

© EPS Vascular AB. All rights reserved. | P +46 (0) 42 33 65 11 | info@epsvascular.com

